



T.C
İSTANBUL MEDİPOL UNIVERSITY
INSTITUTE OF HEALTH SCIENCES

MASTER THESIS

**THE COMPARISON OF PERIPHERAL DOSE IN
STEREOTACTIC BRAIN IRRADIATION WITH THE USE OF
DIFFERENT TREATMENT TECHNIQUES**

HANİEH MABHOUTI

DEPARTMENT OF MEDICAL PHYSICS

SUPERVISOR
Doç. Dr. HİLAL ACAR DEMİR

İSTANBUL – 2017

THESIS APPROVAL FORM

Institute : İstanbul Medipol University Institute of Health Sciences
Level of Programme : Master (X) Doctorate ()
Department : Health Physics
Student : Hanieh MABHOUTI
Title : The Comparison of Peripheral Dose in Stereotactic Brain Irradiation
with The Use of Different Treatment Techniques
Place : Medipol Mega University Hospital
Examination date : 27.03.2017

We have reviewed the present thesis in regard to content and quality and approved as a Master Thesis.

Supervisor

Doç.Dr. Hilal ACAR DEMİR

Institution

İstanbul Medipol University

Signature



Members of Examination Jury

Prof.Dr. Hale Başak ÖZKÖK

İstanbul Medipol University


Prof.Dr. Hatice Bilge BECERİR

İstanbul University



Based on above jury's decision this MSc thesis the board of the Health Science Institute dated on 30.03.2017 and cited 2017/11-01... has approved to this thesis according the thesis Writing Guidance Regulations.


Professor Nesrin EMEKLİ
Director of Institute of Health Sciences



DECLARATION

I, Hanieh Mabhouti, declare that this thesis has been obtained and presented in accordance with academic rules and ethical conduct. Also, it has been submitted in fulfillment of the requirements for the award of the master of science – research, in the Department of Medical Physics, University of Istanbul Medipol, is wholly my own work unless otherwise referenced or acknowledged.

Hanieh MABHOUTI

A handwritten signature in black ink, consisting of several overlapping loops and a long horizontal stroke at the bottom, positioned below the printed name.

ACKNOWLEDGEMENTS

It would not have been possible to complete this work without the help of many people. This integrated position provided me with insight into the power and usefulness of research and development. I would like to thank Dr. Hilal Acar Demir for supporting this initial attempt, and for her support as my supervisor throughout my Master's degree. Thank you for seeing such great potential in me.

I would especially like to appreciate our amazing head of department Prof.Dr Hale Başak Çağlar for providing us all the facilities and her extraordinary support. She genuinely cares about her students and believes that they can succeed.

This thesis would not have been possible without the support, and encouragement I have received from my educators Esra Küçükorkoç, Mine Doyuran, Hande Sertkaya, Reshad Rızazadeh, Doğu Canoğlu, Yasin Ceylan and Sercan Güneş.

I have always been strongly aware of the sacrifices my parents have made so that I would have the opportunity to pursue an advanced degree. I do not have words to adequately describe my deep thanks for all that my dad has provided me with both financially and morally over the course of my life, though I hope to show him in the years to come.

I have benefited greatly from the many great friendships I have made over these two years actively working in Medipol University. Thank you to Mehmet Sıddık Cebe, Pelin Paçacı, Esra Serin, Emre Şanlı, Sevtap Yıldırım, Dursun İştmez for all their help and consultations.

And finally, I would like to thank and dedicate this thesis to my dear educator Nadir Küçük for all his encouragement, support, and endless enthusiasm for life. His patience and motivation during my Masters have been incredible. I am so grateful to know you. You are so enthusiastic about what you do and so knowledgeable too.

CONTENTS

Page number

THESIS APPROVAL FORM.....	i
DECLARATION.....	ii
ACKNOWLEDGEMENTS.....	iii
ABBREVIATIONS AND SYMBOLS.....	vii
LIST OF TABLES.....	ix
LIST OF FIGURES.....	x
LIST OF PHOTOGRAPHS.....	xii
1. ABSTRACT.....	1
2. ÖZET.....	2
3. INTRODUCTION.....	3
4. GENERAL DEFINITIONS.....	5
4.1. Stereotactic Radiotherapy.....	5
4.1.1. Rigid Immobilization.....	6
4.1.2. Conformal Dose Distribution.....	7
4.1.3. Different types of treatment planning evaluation parameter.....	8
4.1.3.1. Conformity Index.....	8
4.1.3.1.1 Paddick Description.....	9
4.1.3.1.2 ICRU Description.....	9
4.1.3.2 Gradient Index.....	10
4.1.3.3 Homogeneity Index.....	11
4.1.4 Three-dimensional positioning accuracy.....	11
4.1.5 Real-time tumor tracking.....	12
4.1.6 Varian Eclipse 13.0.26 Treatment Planning System.....	13
4.1.6.1 Analytical Anisotropic Algorithm (AAA).....	16
4.1.6.2 Acuros XB Algorithm (AXB).....	17
4.1.7 MultiPlan Treatment Planning System.....	19
4.1.7.1 Beam Targeting.....	21
4.1.7.1.1 Isocentric Targeting.....	21
4.1.7.1.2 Conformal Targeting.....	21
4.1.7.1.3 Manual Beam Placement.....	21
4.1.7.2 Dose Optimization Algorithm.....	21

4.1.7.2.1 Simplex Optimization	22
4.1.7.2.2 Iterative Optimization	22
4.1.7.2.3 Sequential Optimization	23
4.1.7.3 Dose Calculation Algorithm.....	23
4.1.7.3.1 Ray-Tracing Dose Calculation Algorithm.....	23
4.1.7.3.2 Monte Carlo Dose Calculation Algorithm	24
4.1.7.3.3 Finite Size Pencil Beam Dose Calculation Algorithm	24
5. MATERIALS AND METHOD	26
5.1.MATERIALS	26
5.1.1 Philips Gemini TF PET-CT 16.....	26
5.1.2 Varian Truebeam 2.0 Stx Linear Accelerator.....	27
5.1.3 CyberKnife M6	29
5.1.4 Atom Dosimetry Verification Phantoms.....	31
5.1.5 Film Dosimetry	32
5.1.6 PTW Farmer Type Chamber TN30013	33
5.1.7 RW3 Slab Phantoms	33
5.2.METHOD.....	34
5.2.1 A CT Simulation of the Phantoms.....	34
5.2.2 Exporting Data into the Treatment Planning System	35
5.2.3 Contouring	35
5.2.4 Treatment Planning in Eclipse Treatment Planning System	36
5.2.5 Treatment Planning in Multiplan.....	37
5.2.6 Preparation of Gafchromic EBT 3 Films	39
5.2.6.1 Scanner Orientation	39
5.2.6.2 Preparation and Irradiation of Films	40
5.2.6.3 Scanning of EBT3 Films.....	41
5.2.7. Ion chamber measurement	42
5.2.8 Set up of CIRS Phantom	42
5.2.8.1. Irradiation of Truebeamstx plans	44
5.2.8.2. Irradiation of MultiPlan TPS plans.....	44
5.2.8.3. Evaluation of irradiated Gafchromic EBT3 films.....	44
6. RESULTS	46
6.1 Film Calibration Curve.....	46
6.2 Peripheral Dose Calculation Difference in TPS.....	47
6.3 Gamma Index Analysis for In field Volumes	49

6.4. MU comparison in TB and CK	51
6.5. Peripheral Dose Measurements in Truebeam STx	52
6.6. Peripheral Dose Measurements in CyberKnife M6	55
6.7. Comparison of Peripheral Doses in Different Modalities	59
7. DISCUSSION.....	65
8. CONCLUSION	68
9. REFERENCES	69
10. CURRICULUM VITAE	78



ABBREVIATIONS AND SYMBOLS

PD	: Peripheral Dose
TPS	: Treatment Planning System
SRS	: Stereotactic Radiosurgery
LQM	: Linear-Quadratic Model
CI	: Conformity Index
GI	: Gradient Index
HI	: Homogeneity Index
MV	: MegaVoltage
KV	: KiloVoltage
CBCT	: Cone-Beam Computed Tomography
CT	: Computed Tomography
MVCBCT	: MegaVoltage Cone-Beam Computed Tomography
KVCBCT	: KiloVoltage Cone-Beam Computed Tomography
IMRT	: Intensity Modulated Radiation Therapy
VMAT	: Volumetric Modulated Arc Therapy
DVO	: Dose Volume Optimization
PGO	: Plan Geometry Optimization
PRO	: Progressive Resolution Optimizer
MRDC	: Multi-Resolution Dose Optimization
MU	: Monitor Unit
DVH	: Dose Volume Histogram
AAA	: Analytical Anisotropic Algorithm
AXB	: Acuros XB
LBTE	: Linear Boltzmann Transport Equation
RT	: Ray-Tracing
MC	: Monte Carlo
FSPB	: Finite Size Pencil Beam
VIO	: Volume Of Interest
MLC	: Multi Leaf Collimator
PET	: Positron Emission Tomography

FFF	: Flattening Filter Free
TB	: TrueBeam
CK	: CyberKnife
PDD	: Percentage Depth Dose
HDMLC	: High Definition Multi Leaf Collimator
DCA	: Dynamic Conformal Arc
3DCRT	: Three-Dimensional Conformal Radiation Therapy
OAR	: Organ At Risk
PTV	: Planning Target Volume
GTV	: Gross Tumor Volume
CTV	: Clinical Target Volume
TLD	: ThermoLuminescent Detector
DICOM	: Digital Imaging and Communications in Medicine
V	: Vertex
NC	: NASal Cavity
PF	: Posterior Fossa

LIST OF TABLES

Page number

Table 6.2.1: Results of PD Differences Between AAA and AXB Algorithms in Eclipse TPS.	48
Table 6.2.2: Results of PD Differences Between Ray-Tracing and MC Algorithms in Multiplan TPS.	49
Table 6.3.1: Evaluation of Irradiated Volumes Using Gamma Index Score for Truebeam STx	50
Table 6.3.2: Evaluation of Irradiated Volumes Using Gamma Index Score for Cyberknife..	50
Table 6.4.1: Calculated MU values for different tumor size and various tumor location in the brain for AAA and AXB algorithms in Eclipse TPS	51
Table 6.4.2: Calculated MU values for different tumor size and various tumor location in the brain for RT, MC and FSPB algorithms in Multiplan TPS	52

LIST OF FIGURES

	Page number
Fig.4.1.2. Example of non-coplanar beams used in SRS treatments.....	8
Fig.4.1.3.1. Description of Conformity Index Value Parameters	9
Fig.4.1.5. Real-time tumor tracking for Cyberknife	13
Fig.4.1.6. Control points of VMAT optimization	15
Fig.4.1.6.1. Example of AAA dose calculation result.....	17
Fig.4.1.6.2. Example of AXB Dose Calculation Result	18
Fig.4.1.7. Multiplan Treatment Planning System.....	20
Fig.5.1.3. Skull Tracking Method.....	31
Fig.5.2.3.1. 4 Different PTV Volumes in Vertex Region.....	35
Fig.5.2.3.2. 4 Different PTV Volumes in NC Region.....	36
Fig.5.2.3.3. 4 Different PTV Volumes in PF Region.....	36
Fig.5.2.5.1. Calculated dose distribution with Ray-Tracing Algorithm.....	38
Fig.5.2.5.2. Calculated dose distribution with Monte Carlo Algorithm.....	38
Fig.5.2.5.3. Calculated dose distribution with FSPB Algorithm	39
Fig.5.2.6.1. Scanner Orientation	40
Fig.5.2.8.3: Evaluation of an irradiated plan in SNC patient software; Upper left section represents dose distribution of irradiated plan, upper right section represents dose distribution of the slice of the irradiated plan, lower left section represents superposition of the both dose distributions with the best alignment, lower right section represents dose agreement graphically	45
Fig 6.1.1: Film Calibration Curve for Truebeamstx with 6MV FFF Beam.....	46
Fig 6.1.2: Film Calibration Curve for Cyberknife with 6 MV FFF Beam	47
Fig 6.5.1: Tumors located in heterogeneous medium	52
Fig 6.5.2: Comparison of calculated PD with AAA in Eclipse and measured using TB for all tumor size in V region.....	53
Fig 6.5.3: Comparison of calculated PD with AXB in Eclipse and measured using TB for all tumor size in V region.....	53
Fig 6.5.4: Comparison of calculated PD with AAA in Eclipse and measured using TB for all tumor size in NC region.	54
Fig 6.5.5: Comparison of calculated PD with AXB in Eclipse and measured Using TB for all tumor size in NC region	55
Fig 6.6.1: Comparison of calculated PD with RT in Multiplan and measured using CK for all tumor size in V region.....	56

Fig 6.6.2: Comparison of calculated PD with MC in Multiplan and measured using CK for all tumor size in V region	56
Fig 6.6.3: Comparison of calculated PD with FSPB in Multiplan and measured using CK for all tumor size in V region	57
Fig 6.6.4: Comparison of calculated PD with RT in Multiplan and measured using CK for all tumor size in NC region	58
Fig 6.6.5: Comparison of calculated PD with MC in Multiplan and measured using CK for all tumor size in NC region.....	58
Fig 6.6.6: Comparison of calculated PD with FSPB in Multiplan and measured using CK for all tumor size in NC region.....	59
Fig 6.7.1 Calculated PD for small 0,6 cm tumor size in V region.	60
Fig 6.7.2 Measured PD for small 0,6 cm tumor size in V region.	60
Fig 6.7.3 Calculated PD for larger 2 cm tumor size in V region.	61
Fig 6.7.4 Measured PD for larger 2 cm tumor size in V region.....	61
Fig 6.7.5 Calculated PD for 1 cm tumor size in NC region.....	62
Fig 6.7.6 Measured PD for 1 cm tumor size in NC region.	63
Fig 6.7.7 Calculated PD for 2 cm tumor size in NC region.....	63
Fig 6.7.8 Measured PD for 2 cm tumor size in NC region.	64

LIST OF PHOTOGRAPHS

	Page number
Photo.4.1.1. a) Thermoplastic mask, b) Head frame	7
Photo.5.1.1. GEMINI TF 16 PET/CT system	27
Photo.5.1.4. Cirs Atom Phantom.....	32
Photo.5.2.1. Phantom Position During CT Simulation.....	34
Photo.5.2.6.2. EBT3 Gafchromic film for calibration.....	41
Photo.5.2.6.3. Scanning of Films	42
Photo.5.2.8.1. position of film in vertex part of Phantom	43
Photo 5.2.8.2. Position of film on nasal cavity part of Phantom	43



1. ABSTRACT

THE COMPARISON OF PERIPHERAL DOSE IN STEREOTACTIC BRAIN IRRADIATION WITH THE USE OF DIFFERENT TREATMENT TECHNIQUES

Brain stereotactic radiosurgery involves the use of precisely directed, single session radiation to create a desired radiobiologic response within the brain target with acceptable minimal effects on surrounding tissues. In this study, the comparison of peripheral dose (PD) was made for Truebeam STx and Cyberknife M6 treatment plans. For Truebeam STx, treatment planning was done using VMAT technique with 6 FFF beam in Vertex, Nasal Cavity, and Posterior fossa regions. The dose distribution was calculated using Eclipse Treatment Planning System (TPS) with AAA and AXB dose calculation algorithms. The treatment plans of the same targets were also done for Cyberknife M6 with Multiplan TPS using Ray-Tracing, Monte Carlo, and Finite Size Pencil Beam dose calculation algorithms. Using the same film batch, the net OD to dose calibration curves were obtained for plans which were done for both devices. Dose distributions of each plan for target site and its periphery were measured using EBT3 film and compared with TPS calculations. For cyberknife plans, the gamma analysis passing rates between measured and calculated dose distributions ranged between 92% and 99.9%. For Truebeam plans, the gamma analysis passing rates ranged between 93.8% and 100%. Although target dose distribution calculated more accurately by Acuros XB and Monte Carlo dose calculation algorithms, Finite Size Pencil Beam algorithm predicted dose distribution around the peripheral region of target lower than the other algorithms.

Keywords: Calculation Algorithm, Peripheral Dose, Stereotactik Radiosurgery

2. ÖZET

STEROTAKTİK BEYİN IŞINLAMALARINDA KRİTİK ORGANLARIN ALDIĞI DOZLARIN FARKLI TEDAVİ TEKNİKLERİ KULLANILARAK KARŞILAŞTIRILMASI

Beyin stereotaktik radyocerrahisi, tek fraksiyonda beyindeki bir hedefe kesin doğrulukta yönlendirilmiş radyasyonun hedefte istenilen radyobiyojik etkiyi oluşturması ve hedefi çevreleyen dokularda oluşacak etkiyi minimal seviyede tutmayı içerir. Bu çalışmada Truebeam STx ve CyberKnife M6 cihazlarında ışınlanan tedavi planlarında periferik dozlar karşılaştırıldı. Truebeam STx cihazı için, verteks, nazal kavite ve posterior fossa yerleşimli hedefler için 6FFF foton demetleri kullanılarak VMAT tekniği ile tedavi planları oluşturuldu. Oluşturulan planların doz dağılımları Eclipse TPS' de AAA ve AXB doz hesaplama algoritmaları kullanılarak hesaplandı. Aynı hedeflerin CyberKnife M6 cihazında ışınlanması için Multiplan TPS' de Ray Tracing, Monte Carlo ve Finite Size Pencil Beam doz hesaplama algoritmaları kullanılarak tedavi planları oluşturuldu. Aynı seriden alınmış filmler kullanılarak her bir cihaz için "Net OD – Doz" kalibrasyon eğrileri elde edildi. Oluşturulan planların doz dağılımları EBT3 film dozimetrisi ile ölçüldü. Her bir planın hedef ve hedefin periferindeki doz dağılımları EBT3 film kullanılarak ölçüldü ve TPS hesaplamaları ile karşılaştırıldı. CyberKnife planlarında gama indeks analizi geçme oranları %92 – 99.1 aralığında oldu. Truebeam planlarının gama indeks analizi geçme oranları ise %93.8 – 100 aralığında oldu. Acuros XB ve Monte Carlo algoritmaları kullanılan diğer algoritmalara göre hedefte daha doğru hesaplama gerçekleştirirken, Finite Size Pencil Beam algoritması periferik bölgede diğer algoritmalara göre daha düşük doz dağılım hesabı gerçekleştirdi.

Anahtar Kelimeler: Hesaplama algoritması, Periferik doz, Stereotaktik radyocerrahi

3. INTRODUCTION

Wilhelm Conrad Roentgen discovered X-rays in 1895. It was one of the major revolutions of that century. After the discovery of X-rays, it was shown that radiation could be used for diagnostic purposes. One year later, kidney stone and a penny in the throat of a child was shown by X-rays in one of the first radiology departments at the Glasgow Royal Infirmary. Hall-Edwards became one of the first people who used an X-ray to diagnose and discovered a needle embedded in a woman's hand. Scientists continued to find out how radiation works and how to measure the dose accurately between World War I and II(1).

It was eventually recognized that not only in diagnostic but also it can be used in the therapeutic field to treat malign lesions in the body. As the time progressed, especially computer and other new technologies were rapidly developed, and new treatment modalities were implemented to treat cancer patient(2). Nowadays it continues as a treatment for cancer in radiation therapy, where it is applied with strict safety precautions(3).

The primary goals of radiation therapy are to shrink tumors and kill cancer cells. The use of radiation therapy for the control of malignant diseases results in both normal and malignant tissues being irradiated and damaged. By using new radiation therapy techniques such as IMRT, Volumetric Arc Therapy or Stereotactic Radiation Therapy, we can concentrate precisely on the dose within target sparing the surrounding normal tissues(4).

Stereotactic Radiosurgery is a different aspect of radiotherapy that refers to the focused delivery of high dose in a single fraction. It works the same as another form of radiation treatment and directly affects the DNA of the tumor cells. It is essential to achieve a high conformity and steep dose gradient to reduce organ at risks exposure while using this technique. Several studies have been carried out with the aim of identifying and measuring doses to be deposited outside the treatment volume(5).

The peripheral dose (PD) is defined as the dose outside of the therapeutic radiation field edge. Three components affect PD. These are:

- Internally scattered radiation from within the treatment field
- Leakage from the linac head
- Scattered radiation from the collimators

The amount of PD is crucial since it increases the probability of secondary cancer risk in the patient. Some anatomic sites such as breast and thyroid could be effected from the PD more than the others and result in increasing secondary cancer risk probability. Some studies have shown that radiation-induced cancer risk increases linearly with increasing dose(6).

The main purpose of this study was to investigate the accuracy of dose calculation algorithms of Eclipse and Multiplan Treatment Planning System (TPS) for out of field doses. The calculated PD were compared with that of measured by the film.

4. GENERAL DEFINITIONS

4.1. Stereotactic Radiotherapy

Stereotactic Radiosurgery (SRS) is a particular type of radiation therapy that is most commonly used to treat tumors in the brain in which focused type of radiation are given to stereotactically determined tumor volume in one to five treatment sessions. SRS is not a form of unlike open surgery, and there is not any incision. It can deliver exactly-targeted radiation within narrow margins in fewer high-dose treatment fractions than conventional radiation therapy which can spare as much as possible healthy tissues according to the highly sophisticated computer-based radiation delivery system by using 3D imaging(7).

SRS works by damaging the DNA of the tumor. Therefore, the affected cells then lose the ability to reproduce which causes to slow or stop tumor growth. Radiobiological principles of SRS are currently not well understood. However, attempts have been made to rationalize the delivery of a single large fraction of dose to a small circumscribed lesion in the brain(8). It has been shown that irradiation with high doses, in a single fraction causes severe vascular damage in human tumor xenografts or animal tumors(9). High dose radiation such as $>10-12$ Gy in a single exposure is likely to cause significant vascular damage followed by indirect cell death. Some clinical studies have demonstrated that SRS of cranial tumors with high-dose lower fractions is very effective in achieving local tumor control. Linear Quadratic Model (LQM) is one of the best models for comparing different fractionated radiotherapy protocols in conventional radiation therapy, but there have been still some discussions if this model is applicable for SRS treatments. In 2008, Brenner et al showed that LQM could be used up to 15-18 Gy, but need to be aware that it has been progressively less accurate at doses above 10 Gy(10).

SRS treatment steps are different than conventional radiation therapy. It needs to work precisely and accurately in all steps. Critical steps in SRS treatment are listed below:


- Rigid immobilization of the patient
- Highly conformal dose distribution
- Different type of treatment planning evaluation parameters
- Three-dimensional positioning accuracy
- Real-time tumor tracking

4.1.1. Rigid Immobilization

Rigid immobilization of patients and accurate positioning of their targets have long been recognized as critically important aspects of quality of SRS treatment. Accurate positioning combined with strict immobilization permit reduced margins around the target, resulting in a decrease in dose to healthy tissue and a potential increase in dose to the target. The increasing use of computed tomography-based three-dimensional treatment planning programs has made highly conformal dose distributions possible, thus further emphasizing the need for accurate positioning. Special immobilization techniques for SRS treatment may allow reducing margin even zero for intracranial lesions(11). Thermoplastic mask and stereotactic head frames are commonly used for intracranial SRS treatments (Photo 4.1.1). Thermoplastic masks are the most simple and standard immobilization method in SRS treatment in all over the world. Immobilize patient with the stereotactic head frame is more rigid way than using the thermoplastic mask. Patient has local anesthesia during the procedure and frames are fixed to the cranial bone by using small screws(12).



Photo.4.1.1. a) Thermoplastic mask



b) Head frame

4.1.2. Conformal Dose Distribution

One of the essential functions of SRS is to achieve highly conformal dose distribution. In order to get it, many steps should be done. These steps are;

- Use a rigid/nonrigid immobilization device
- Determine accurate location of lesion
- Contouring lesion and adjacent anatomy from different imaging platforms
- Using multiple planar/non-coplanar beam arrangement
- Verify correct position of the lesion in a machine and accurately delivery of a conformal plan.

The beam geometry used in radiosurgery is very different from that used in conventional radiotherapy, in which a large number of beams spread with various angles as shown in (Fig 4.1.2). It could be obtained a very high maximum relative dose in the central part of the target and sharp dose gradients at field edges by using non-coplanar beam arrangement. A high maximum dose in the central parts of the target may be beneficial if the tumor contains anoxic, more radioresistant, cells in the central parts. A high dose in the anoxic cells will thus increase the chance of killing more the clonogenic tumor cells(13).




Fig.4.1.2. Example of non-coplanar beams used in SRS treatments

4.1.3. Different types of treatment planning evaluation parameter

The assessment of SRS treatment plans is entirely different from conventional ones. There are some additional criteria which must be checked to make sure if plans quality are acceptable for delivery. Also, dose constraints are entirely different from the conventional scheme because of the high dose and low fraction. The treatment planning indices which are used in SRS plan listed and discussed below:

- Conformity Index
- Gradient Index
- Homogeneity Index

4.1.3.1. Conformity Index

Radiosurgical targets are typically nonspherical, except for brain metastases. Some tumors exhibit more complex shapes than others(14). One goal of radiosurgery is to design a treatment plan in which the prescription isodose line covers the target with a minimal excess volume as shown in (Fig 4.1.3.1). To facilitate comparison and to evaluate compliance with clinical trial protocols, some of the plan conformity indexes has been proposed(15). Each index takes factors such as the overlap and size of the dose and target volumes into account. Several different conformity index

formulas have been reported to describe the conformity of the prescription isodose to the target volume(16). Two of them are commonly used for evaluating SRS plans, these are;

- Paddick Formula
- ICRU Formula




Fig.4.1.3.1. Description of Conformity Index Value Parameters

4.1.3.1.1 Paddick Description

$$\text{Paddick CI} = \frac{(\text{TVPIV})^2}{(\text{TV} \times \text{PIV})}$$

TVPIV=Target Volume covered by PIV

TV = Target Volume

PIV = Prescription Isodose Volume

4.1.3.1.2 ICRU Description

$$CI = \frac{TV}{PTV}$$

TV=Treated volume

PTV=Planning target volume

4.1.3.2 Gradient Index

A dose gradient index (Cozzi et al.) can be utilized to analyze treatment plans with equivalent conformities. The sharp dose gradient out of the radiosurgical target is one of the variables that makes radiosurgery achievable. Thus it is reasonable to quantify this variable and utilize it to make comparison between rival treatment plans explore optimum prescription isodoses, or make comparison of treatment modalities(14). The GI is explained as the ratio of the volume of half the prescription isodose to the volume of the prescription isodose. For a plan normalized to 100% isodose line, it is the ratio of 50% isodose volume to that of 100% isodose volume. Paddick determines GI formula like below:

$$GI = \frac{PIV_{half}}{PIV}$$

PIVhalf = Prescription isodose volume, at half the prescription isodose (e.g. at 25%)

PIV = Prescription isodose volume (e.g. at 50%)

4.1.3.3 Homogeneity Index

Homogeneity Index (HI) is an objective tool to analysis the uniformity of dose distribution in the target volume. In SRS treatment, we could create heterogeneity inside target volume. It is very useful for most cases to increase target dose and creates a strict dose drop off after the boundary of a target(17). RTOG prescribed HI index like below;

$$HI = \frac{MD}{PD}$$

MD = Maximum Dose within Target Volume

PD = Prescribed Dose to the Target

4.1.4 Three-dimensional positioning accuracy

Accurate and reproducible patient setups can potentially reduce harm to critical organs during radiation treatment. Patient set-up accuracy is especially critical for high -intensity treatments, such as SRS. Linear accelerator linac-based SRS treatments have been performed for over 25 years, however, treatment verification was indigent because of invisibility of megavoltage portal images(MV) and the unavailability of high-contrast kilovoltage (kV) imaging devices in the treatment room until last decade. Cone-beam computed tomography (CBCT) imaging is a relatively new capability in the treatment room for image-guided patient treatment that can be potentially used for high-precision SRS patient setup and treatment verification(18). The initial use of daily computed tomography (CT) has been for assessing internal organ position and defining the subsequent isocenter shifts to be performed at the treatment unit. CBCT is now one of the most common tools in the clinics to determine where the tumor is and how accurately can treat patients in the linear accelerators(19).

Megavoltage cone-beam CT (MVCBCT) uses the treatment beam of the accelerator and its portal imaging system to provide volumetric datasets with sufficient contrast for image-guidance. High energy photon beams are used to create MVCBCT in the machine. Clinical implementation of both kV-CBCT and MVCBCT systems make users to treat the patient safer than before(19). Increasing clinical evidence suggests that 6 degrees of freedom couch used in conjunction with cone-beam computed tomography(CBCT) may reduce patient set-up errors and improve target positioning for high-dose SRS treatments(20).

4.1.5 Real-time tumor tracking

High dose is delivered within a very short time in SRS treatment. Therefore intrafraction motion of tumor in SRS is more critical than conventional therapies. In recent years, some companies have produced a different type of tumor tracking methods which encompass tumor movement during the treatment. Immobilization of cranium is easier than another part of the body, but still, real-time tumor tracking is essential(21).

Both linac and robotic systems have different types of tumor tracking systems. Linacs use kV imaging or some software to track tumors. For example, imaging system could be utilized in fluoroscopic mode during the treatment and see real time if targets are inside irradiation volume. Cyberknife uses 6D Skull tracking algorithm which has frameless tracking of bony anatomy of the skull (Fig 4.1.5)(20). Users define alignment center during treatment planning, and the system generates DRRs to correlate with acquired images within treatment delivery. During the treatment, initial alignment performed by adjusting patient and couch then corrections applied to robot throughout treatment. The system assumes fixed relationship between target, align center and bony anatomy of the skull. 6D Skull tracking algorithm correlates acquired images of skull with DRRs. Correlation is based on pixel similarity criteria. The algorithm is used for intracranial, head and neck and some upper spine(C1-2) lesions (22).



Fig.4.1.5. Real-time tumor tracking for Cyberknife

4.1.6 Varian Eclipse 13.0.26 Treatment Planning System

Varian's Eclipse treatment planning system optimizes a radiotherapy treatment plan based on instructions, and information about the size, shape, and location of the tumor to be treated with radiation. The treatment plan tells the radiotherapy machine how to deliver the treatment, what gantry angles to use, how much radiation dose to deliver from each angle, and how the treatment beam should be shaped to create proper dose distribution. Such a sophisticated treatment technique like Intensity Modulated Radiation Therapy, Volumetric Modulated Arc Therapy, etc. Could be easily calculated by this software. It offers users a set of treatment plan templates that make the process efficient for creating a personalized plan for each patient.

Different types of algorithm are used for dose optimization and calculation. The planning system chooses dose optimization algorithm according to treatment technique which is decided at the beginning of the treatment plan(24). Dose Volume Optimization (DVO), Plan Geometry Optimization (PGO), Progressive Resolution Optimizer (PRO) and Multi-Resolution Dose Calculation (MRDC) are used in the system.

The algorithm used in Eclipse IMRT, DVO, determines the optimal field shape and intensity by iteratively conforming the dose distribution to the desired objectives until an optimum solution is reached. This algorithm based on convolution-superposition principle. The DVO algorithm can use a calculated plan dose as an intermediate dose when optimizing a plan. The DVO algorithm calculates the difference between the intermediate dose and the first round optimization result and uses this difference to compensate the optimization lead to the consequent iterations(25).

PGO algorithm enables the Eclipse Beam Angle Optimization, an integrated optimization option for the Eclipse treatment planning system. Beam Angle Optimization is an automated tool for selecting the suitable beam angles based on user-defined dose-volume objectives that speed up the planning process and is performed with the PGO algorithm, which relies on the Eclipse DVO algorithm. Global optimization creates the new field geometry, which can be either coplanar or non-coplanar, depending on user-defined optimization parameters. The PGO uses a fixed isocenter, which means that all fields in the initial field distribution share the same isocenter. The initial number of fields in both geometries is controlled with a parameter. You can also control the offset for the gantry angles in the coplanar field geometry with the coplanar offset angle parameter. This parameter does not affect the non-coplanar initial field distribution. In the non-coplanar field geometry, you can enter a limit for the elevation angle of the fields from the coplanar plan. The fields in the initial field distribution do not have elevation angle values higher than the specified limit.

PRO algorithm creates VMAT (RapidArc) plans based on dose-volume objectives. VMAT fields use Dynamic MLC, variable dose rate and variable gantry speeds(23). The PRO algorithm generates a sequence of control points which define MLC leaf positions and MU/deg as a function of gantry angle. MU/deg is encoded in dicom and the Varian system database with the cumulative meter set weight, which defines the increase in MU between control points about the total MU in the field. The initial conditions for the PRO algorithm are determined using control points to

represent each VMAT field (Fig 4.1.6) The PRO algorithm uses multi-resolution approach to optimize the plan. It means that the dose is modeled using first a lower number of dose calculation segments that are distributed evenly in each field(25). The number of dose calculation segments increases when moving from one multi-resolution level to another.




Fig.4.1.6. Control points of VMAT optimization

MRDC algorithm is used for fast dose estimation inside the DVO, PRO and PGO algorithms to improve the optimization accuracy, which can be seen as a good agreement between the optimization Dose Volume Histograms (DVH) and the last DVHs. The high speed of the MRDC algorithm allows the optimization algorithms to perform full dose computation during each iteration(26). The MRDC algorithm is based on the convolution superposition principle, and it uses 3D convolution scattering calculation. The scattering model is based on the 3D superposition of point spread functions in the patient model. The point spread functions are built from Monte Carlo calculations.

Eclipse treatment planning system uses different algorithms during the dose calculation. These are:

4.1.6.1 Analytical Anisotropic Algorithm (AAA)

Analytical Anisotropic Algorithm (AAA) is immensely utilized as a part of clinical practice. The AAA dose calculation model is a 3D pencil beam convolution superposition algorithm. The algorithm has separated models for initial photons, secondary central photons which are scattered, and electrons scattered from collimation systems. Tissue heterogeneities are considered as anisotropic in the whole 3D neighborhood using 13 lateral photons scatter kernels. The final dose distribution is obtained by superposition of the doses from the photon and electron convolution(27).

The AAA dose calculation model is comprised of two main components:

- Configuration Algorithm
- Actual Dose Calculation Algorithm.

The configuration algorithm is utilized to obtain physical parameters needed to characterise the fluence and the energy spectra of the photons and electrons that present in the clinical beam and their essential scattering characteristics in water equivalent media (28). Even though some parameters used for the dose calculation algorithm can be derived for acceptable accuracy from basic measurements of depth dose and lateral dose profile water equivalent phantoms, in practise experimental determination of whole parameters is impossible. This is settled in AAA model via all pre-computed parameters using MC simulations and modified parameters in order to agree with real measured clinical beam data in the course of beam data configurations (29). Once machine specific datas are determined according to proper procedures during the beam configuration phase, all parameters are stored. All these parameters are retrieved during the actual dose calculation process.

The broad clinical beam is divided into small, finite-sized beamlets to which the convolutions are applied. The final dose distribution in (Fig 4.1.6.1) is obtained by the superposition of the dose calculated with photon and electron convolutions for the individual beamlets. AAA accounts for tissue heterogeneity anisotropically

within interaction site in 3D volume. This is performed by the use of radiological scaling of the dose deposition functions and the electron density based scaling of the photon scatter kernels independently in four lateral directions(30).




Fig.4.1.6.1. Example of AAA dose calculation result

4.1.6.2 Acuros XB Algorithm (AXB)

Another important algorithm which is used in Eclipse planning system is AXB. AXB algorithm was developed to provide accurate and rapid dose calculations for treatments ranging from 4-25 MV. It utilizes advanced techniques in order to find the solution of Linear Boltzmann Transport Equation (LBTE). It also directly takes impact of heterogeneities into account of patient dose calculations (31). Different materials such as air, bone, lung and implants that are composing heterogeneities considerably effects the dose distribution within patient. This effect may be more important particularly when there are irregular or small fields. AXB can also obtain datas from source model of AAA and utilize them for its beam model in order to adopt it for the calculation (25). It supplies similar accuracy with MC for treatment planning in the full range of X-ray beams (Fig 4.1.6.2).




Fig.4.1.6.2. Example of AXB Dose Calculation Result

The LBTE is the essential equation describing macroscopic characteristic of particles when they move through a matter. For a specified portion of matter which is exposed to radiation, the solution of the equation can give the exact explanation of the dose inside the portion. There are two main way to obtain open for solution of LBTE: the first was is MC method. The second way is finding solution of LBTE with explicit numerical methods. MC and explicit LBTE are both convergent. The potential of accuracy for both ways is the same and converge on the same solution (32). But in practical terms neither of the two ways are precise and both ways generate errors. AXB of Eclipse Source Model have four constituents:

Primary Source: User defined annular or elliptic source placed on the target plane that models bremsstrahlung x-ray photons produced within the target which doesn't have interaction with the treatment head.

Extra Focal Source: Gaussian plane source placed below the flattening filter that models x-ray photons as a result of interaction within the treatment head out of the target (primary within FF, primary collimators, and secondary jaws).

Electron Contamination: Accounts for the dose stored in build up area that is not explained by primary and extra-focal source constituents.

Photons Scattered From Wedge: Accounts for scattering from hard wedges. Carried out with dual Gaussian model that gaussian kernel width get bigger as the distance to the wedge increase.

The Acuros XB patient transport have four individual courses:

- Transport of source model fluence into the patient
- Calculation of scattered photon fluence in the patient
- Calculation of scattered electron fluence in the patient
- Dose calculation

Application of courses from course 1 to course 3 is used to calculate the electron beam fluence of each voxel within the patient. When energy dependent electron fluence is solved the intended dose amount (dose to medium or dose to water) is computed in the course 4. Only course 1 is repeated for each beam and course 2 and course 4 applied only once despite number of beams. When it comes to VMAT, each beam has important number of orientations and Course 1 will be repeated for each orientation and from course 2 to course 4 will be applied only once. In course 1 the machine sources are modeled as external sources and ray tracing is implemented for calculation of collided x-ray photon fluence and electron fluence within the patient. In course 2 and course 3 AXB decomposes in space, angle, and energy, and iteratively find the solution of LBTE. In course 4 dose of any voxel in the problem is determined by the application of an energy dependent fluence to dose response function to the local energy dependent electron fluence in that voxel (32).

4.1.7 MultiPlan Treatment Planning System

The MultiPlan Treatment Planning System is a highly interactive, workflow-based software application designed for radiosurgery and high precision radiation therapy planning. The MultiPlan System provides tools necessary to generate quality treatment plans for delivery using the CyberKnife. It enables the creation of treatment plans using all tracking algorithm which is defined in CK system. It could

automatically generate precise contours for intracranial and pelvic anatomy using both model-based and atlas-based delineation methods. It has a Sequential Optimization Algorithm which provides an intuitive and intelligent plan optimization algorithm for rapidly developing custom-tailored treatment plans specific to the unique clinical objectives for each patient(33).

Multiplan uses different algorithms during the dose calculation. Ray tracing is one of the common algorithms which is used Multiplan calculations except for heterogeneity areas. It is a simple pencil beam algorithm, and it calculates dose in the ray of the beam. Monte Carlo (MC) algorithm which is often considered the gold standard for dose calculation produces results in minutes. The Finite Size Pencil Beam(FSPB) algorithm for photon consists of decomposing a radiation field into small beamlets. The FSPB algorithm is based on an analytical kernel describing the dose distribution of each beamlet. The parameters of this kernel are determined from the dose of broad beams, and total dose is calculated by summing the dose contributions for each beamlet. (Fig 4.1.7)(34)




Fig 4.1.7. Multiplan Treatment Planning System

4.1.7.1 Beam Targeting

Treatment planning process has various step in Multiplan. After the collimator size has been selected, beams must be directed at the treatment target. The system includes three methods for beam targeting(35).

4.1.7.1.1 Isocentric Targeting

Treatment planning system permits the user to target beams isocentrically. For each defined isocenter, a single beam from each treatment robot node is directed toward the specified isocenter. For multiple or single isocenter treatments, you can reduce the treatment time by only using beams from 1 or 2 paths to treat a defined isocenter(36).

4.1.7.1.2 Conformal Targeting

Treatment planning system includes algorithms to target beams conformally based on the shape of the tumor. To select the MC Dose Calculation as the input to a conformal treatment optimization procedure, an existing optimized treatment plan based on the Ray-Tracing (RT) dose calculation algorithm must first exist.

Only those beams with at least one Monitor Unit in the treatment plan that was generated initially using the RT dose calculation algorithm will be recalculated using the MC algorithm and be available for re-optimization(36).

4.1.7.1.3 Manual Beam Placement

Treatment planning system allows the user to place beam manually if needed.

4.1.7.2 Dose Optimization Algorithm

After beam targeting process, the optimization algorithms could be chosen, and calculation step started. This is performed in the two-step process. First, candidate beams are generated, and their dose distribution is calculated. Second, the

relative weight of each beam is optimized by minimizing a cost function defined by the similarity of the resulting dose distribution to dosimetric and MU goals specified by the use three algorithms are available for second part of the procedure(36).

4.1.7.2.1 Simplex Optimization

Simplex optimization is a standard linear algebra method for solving optimization problems. Considering the set of beams identified with subscript j and a set of volumes of interest (VOI) identified with subscript i . For each target VOI, and organ at risk VOI system creates a formula and solve it according to which was entered(36)(35).

4.1.7.2.2 Iterative Optimization

The difficulty of impracticability was coped with the introduction of the Iterative Optimization algorithm by annulling settled restraints. In this approach, the cost capacity is calculated as a weighted sum of the maximum and minimum dose deviations from user defined goals. The optimum set of bar weights is achieved through an iterative hunt technique similar to slope plummet. Initially, each candidate beam is assigned an MU setting based on the ratio of the total dose it delivers to all points within the target volume about to all points within the patient (the larger this ratio, the greater the initial MU setting). The iterative calculation then takes the effect of increasing and decreasing the MU setting into consideration by a fixed increment of each beam in turn and retains the setting that gives the largest reduction in the cost function. When no further reduction is possible by altering any beam weight, the increment is reduced, and the process repeats. The calculation also incorporates beam retargeting, where a proportion of the candidate beams with zero MU are retargeted at areas of low dose within the target volume to form new candidate beams. In common with the Simplex calculation, though, the reliance on simultaneous optimization of multiple conflicting dose objectives requires the use of non-intuitive manual weighting factors (36)(35).

4.1.7.2.3 Sequential Optimization

Optimization approaches where multiple objectives are grouped in a single cost function and optimized simultaneously. The Sequential Optimization algorithm is executed sequentially as a series of individual optimization steps. Each step performs a linear programming optimization applied to a single objective cost function, designed to correspond to a particular clinical purpose(35). The available objectives include maximization of target volume coverage by a defined dose level, target volume dose homogeneity, conformity of the dose distribution around the target volume, minimization of maximum or mean dose for critical structures, minimize the volume of critical structures that exceed a dose, maximize the volume of target structures that exceed a dose, and minimization of total monitor units. The sequence of steps is defined by the user to match the overall clinical goals(36).

4.1.7.3 Dose Calculation Algorithm

MultiPlan Treatment Planning System (Accuray) is based on inverse planning using linear optimization. Dose calculation algorithm of this treatment planning system are:

- Ray-Tracing Dose Calculation Algorithm
- Monte Carlo Dose Calculation Algorithm
- Finite Size Pencil Beam Dose Calculation Algorithm

4.1.7.3.1 Ray-Tracing Dose Calculation Algorithm

A Ray-tracing algorithm provides a fast dose calculation method based on measured beam data look up. Heterogeneity correction is performed using effective path length, and obliquity correction is carried out by casting multiple rays within each beam. The RT dose calculation algorithm uses 3 system-specific beam description comprised of data measured in water using a water phantom. Measured data are Tissue Phantom Ratio, Off-Center Ratio, and Output Factor. The effective depth is determined by summing the contribution of each voxel along the ray from

the source to plane containing the target voxel using the CT electron density about water(35).

4.1.7.3.2 Monte Carlo Dose Calculation Algorithm

The Monte Carlo dose calculation algorithm samples the interactions of photons entering the patient. To provide an accurate calculation of the dose deposited in the patient, sufficient photons are sampled. The algorithm includes 2 models: first, a source model that describes the distribution of the energies and second trajectories of photons exiting the linac. The patient geometry is modeled as 3D arrays of mass density and material type. Photons interact with the patient volume to form particles, such as electrons and positrons, and deposit energy. Each particle is tracked and, for each voxel, the deposition of its energy is calculated(35). The total dose recorded at each voxel is an estimate of the dose deposited by all photons in the original treatment beam. The Monte Carlo dose calculation algorithm carries out the necessary computations in a matter of minutes rather than hours without compromising dose calculation accuracy. The algorithm uses a beam commissioning procedure which requires several days of computation time to derive source model parameters automatically based on measured beam data; however, it is only performed during the commissioning process and can be broken up into a few hours at a time. The Monte Carlo dose calculation can be enabled for both fixed collimators and the Iris Collimator(34).

4.1.7.3.3 Finite Size Pencil Beam Dose Calculation Algorithm

Dose calculation for Multileaf Collimator plans is performed using a finite size pencil beam. Dose optimization for small field beamlets requires the computation of a large number of slight fields of typically a few mm to 1 cm in size. Some dose calculation methods based on analytical kernels such as FSPB were specifically designed for the purpose of beamlet-based. A beamlet corresponds to a complex fluence distribution: all particles that pass through a small rectangular opening in the field in a given plane. It becomes apparent that in the presence of collimator scatter, the fluence distribution of a beamlet depends on the position of all

collimators, usually multi-leaf collimators (MLCs) and jaws. Consequentially, the interpretation of a beamlet as a ‘unit fluence element’ is to some extent oversimplified, and in any case rather unwieldy for direct dose computation(34).

One of the marks of FSPB is that, can be a beam edge anywhere in the field, in particular, if the projections of organs at risk and the target overlap. Many MLC segments are used to irradiate the target volume. As a consequence, the dose, especially in the mostly shielded organs at risk, stems from out-of-field penumbra dose (leakage, scatter) to a large degree, rather than primary fluence. Here, it is essential to use the correct measure of the beam penumbra during the dose optimization to facilitate the accurate placement of dose gradients.

5. MATERIALS AND METHOD

5.1. Materials

This study has been done in Istanbul Medipol University Radiation Oncology Department. All the machines and other materials in this thesis belong to this institute which is listed below:

- Philips Gemini TF PET/CT16
- Varian truebeam 2.0 STx linear accelerator
- CyberKnife M6
- Atom dosimetry verification phantoms
- Film dosimetry
- PTW Farmer-type chamber TN30013
- RW3 slab phantoms

5.1.1 Philips Gemini TF PET-CT 16

CT scanner is widely used in radiation oncology clinics to simulate patient anatomy and preparation of treatment plan before delivery. The Gemini TF 16 is a diagnostic imaging system for fixed or mobile installations that combine Positron Emission Tomography (PET) and X-ray (CT) system (Photo 5.1.1). The CT subsystem produces cross-sectional images of the body by computer reconstruction of X-ray transmission data. Both subsystems (PET and CT) can also be operated independently as fully functional, diagnostic imaging systems including an application of CT scanner as a radiation therapy simulation scanner. CT scanner has connectivity to Aria (Varian) Oncology Information system and Eclipse treatment planning system.




Photo.5.1.1. GEMINI TF 16 PET/CT system

5.1.2 Varian Truebeam 2.0 Stx Linear Accelerator

Varian has recently presented a new class of linear accelerator referred to as TrueBeam linear accelerator (TB-LINAC). This platform delivers both traditionally flattened photon beams and flattening-filter-free (FFF) photon beams. The TB-LINAC is equipped with a newly designed waveguide, carousel assembly, monitoring control, and integrated imaging systems(37).

TB-LINAC is digital in design and in energy switch that is continuously adjustable to allow more beam energies, Ionization chamber with additional segments to allow better beam monitoring of flatness (or beam profile for FFF beams) and symmetry. The dose rates can be up to 1400 MU/min for 6MV FFF and 2400 MU/min for 10 MV FFF, respectively, which can significantly reduce the beam-on time, thereby limiting the dosimetric impact of intrafraction tumor motion. TB-LINAC is capable of delivering photon beams with 6 MV, 6 MV FFF, 10 MV, 10 MV FFF, 15 MV as well as electron beams(38).

Dosimetric characteristics of the TB-LINAC treatment units are systematically measured for commissioning terms of their percentage depth dose (PDD) curves, beam profiles, relative scatter factors, dosimetric leaf gaps, output

factors, and MLC leakage. High-resolution diode detectors and ion chambers are used to measure dosimetric data for a range of field sizes from 0.5x0.5 to 40x40 cm².

The FFF modality has been introduced to increase dose rate and reduce leaf transmission, head scatters, and leakage radiation with the removal of the flattening filters. There is a noticeable dose reduction outside of the field in FFF beams compared to FF beams, which can improve the target conformity and have sharper dose drop-off to limit radiation to distant organs(30). However, to achieve uniform dose to the tumor, in particular for the large tumors, it needs the distance-dependent modulation, since the intensity of the beam decreases sharply with the off-axis distance. This could increase MU and offset the advantage of using FFF beams(38).

Both Millennium and High Definition Multileaf Collimator (HDMLC) with 120 leaves can be installed on TB-LINAC. The tongue and groove and rounded leaf edges of Millennium MLC are still used in the updated design. In new design HDMLC, innermost 32 pairs of tungsten leaves are 2.5 mm wide, and the 28 outer pairs of leaves are 5.0 mm wide, as projected to isocenter. When the HDMLC is in use, irregularly shaped fields up to 40 cm wide and 22 cm long can be generated. Both dynamic and conformal arc treatments can be used with the TrueBeam STx, as well as static and step and shoot treatments(29).

Patient respiratory motion is monitored using an integrated 3D tracking system consisting of a stereo-view infrared video camera installed above the couch. A respiratory signal obtained from the RPM™ system (Varian Medical Systems, Inc) showing the vertical motion of an optical marker is placed on a patient's chest wall versus time(39).

Truebeam offers different types of treatment modality for the patients. Dynamic Conformal Arc(DCA), Intensity Modulated Radiation Therapy(IMRT) and Volumetric Arc Technique(VMAT) are commonly used in stereotactic treatments in Truebeam. Treatment is delivered continuously as the gantry rotates around the patient in DCA technique. However, during this rotation, the leaves simply conform to the shape of the target. IMRT is an advanced form of three-Dimensional

conformal radiotherapy (3DCRT). It uses sophisticated software and hardware to vary the shape and create the intensity of radiation delivered to different types of tumor in the body. VMAT is a different type of IMRT. The machine rotates around the patient and continuously reshapes and changes the intensity of the radiation beam. MLC is used to create a beam intensity according to tumor shapes(40).

5.1.3 CyberKnife M6

The CyberKnife (CK) Robotic Radiosurgery System is a medical device designed for SRS and SBRT treatments. It enables the delivery of radiation from multiple beams with stereotactic precision provided by image guidance. Throughout the treatment, target localization is achieved by automatic registration of two orthogonal live X-ray images with a library of digitally reconstructed radiographs generated from the patient's planning CT. There are different types of image registration methods in CK. Image registration has been used by gold markers which are implanted in soft tissue such as prostate cancer, however, skull or vertebral bodies based on bony anatomical landmarks and for lung tumors center of mass of lung lesions are used(41).

Cyberknife consist of a linac mounted on a robot. The robot has an arm which can move 6 degrees of freedom. The robot corrects translation errors. Rotational errors are displayed on the monitor and can be fixed by rotating the treatment couch. In a typical CK plan, hundreds of nonisocentric and noncoplanar radiation beams are pointed to the target, creating a highly conformal dose distribution with sharp dose drop-off at its periphery and low dose to adjacent organs at risk (OARs)(33). These characteristics make CK an ideal machine for treatments that require high spatial accuracy and high conformity, such as SRS and SBRT treatments. Targeting accuracy is within 1 mm for both static targets such as cranial or spinal tumors, as well as dynamic targets such as lung tumors. High targeting accuracy allows smaller planning target volume (PTV) margins(21)(42).

The treatment delivery system consists of a 6 MV linear accelerator mounted on a robotic arm that can deliver radiation from hundreds of robot positions spaced uniformly around the target. Three types of collimator system are used to control beam size. Fixed collimator system and a variable aperture Iris collimator, both with the choice of 12 circular beams of diameters ranging from 5-60 mm. Fixed collimators can be changed to vary the beam size as generated by the treatment plan. For each fixed collimator, the manipulator traverses a separate path. The Iris Collimator creates beams with characteristics virtually identical to those of fixed collimators. It consists of two banks of 6 tungsten segments each with each bank forming a hexagonal aperture. The two are offset by 30° about each other resulting in a dodecahedral (12-sided) aperture when viewed from one end of the collimator to the other. Incise MLC system is newly introduced in the latest model of CK. The MLC includes 52 flat-sided leaves, each of which is 90 mm tall and projects 3.85 mm width at the nominal treatment distance of 800 mm SAD. The design allows full overtravel and unrestricted interdigitation. It is a secondary collimator whose aperture is adjustable under computer control. Using tungsten leaves to adjust the aperture rapidly, the MLC can deliver variable shaped beams from each LINAC position. Leaf position is determined by primary motor encoders and is checked with a secondary optical camera system. It is exchangeable with the alternate fixed and variable circular aperture collimator systems(43).

The algorithm which is used to treat intracranial tumors called Skull Tracking Algorithm in CK. It is a frameless tracking of bony anatomy of the skull. It can be used intracranial, head and neck and some upper spine lesions. The system assumes fixing the relationship between target align center and bony anatomy of the skull. The user defines align center during treatment planning. The system generates DRRs to correlate with acquired images during treatment delivery. The correlation is based on pixel similarity criteria. The corrections are applied to robot throughout treatment. (Fig 5.1.3)

Lesions which are located in the body (except brain/head and neck), target tracking can be carried out different types of a tracking system. Fiducials could be used for soft tissue tumors. X-Sight Spine Tracking System can be used for targets in the spine. For moving targets, X-Sight Lung Algorithm utilized with Synchrony Respiratory Tracking System. The Synchrony Respiratory Tracking System continuously synchronizes treatment beam delivery to the motion of a target that is moving with respiration.




Fig.5.1.3. Skull Tracking Method

5.1.4 Atom Dosimetry Verification Phantoms

CIRS ATOM phantoms are a full line of an anthropomorphic phantom. It is a cross-sectional dosimetry phantom and possible to investigate organ dose, whole body effective dose or verification of delivery of therapeutic radiation doses. Phantom is sectional in design with traditional 25 mm thick sections. The sectional surfaces are extremely flat and smooth and do not require any special coatings or treatment. Tissue-equivalent epoxy resins are used in all aspects of the phantom.

CIRS technology offers superior tissue simulation for Dosimetry Verification Phantoms by covering a wider range of energy levels from diagnostic to therapeutic. Also, all bones are homogeneous, average bone composition. CIRS bone formulations offer distinct advantages over natural skeletons and other types of simulated bone. (Photo 5.1.4)




Photo.5.1.4. Cirs Atom Phantom

5.1.5 Film Dosimetry

The film we used for this study was GAFCHROMIC EBT3 (Lot number 11031501). Dimensions of each sheet were $20.3 \times 25.4 \text{ cm}^2$. The development in EBT3 is the symmetric structure which can prevent the possible errors in measurements of optical density due to scan side in EBT2. GAFCHROMIC EBT3 radiochromic dosimetry film is composed of a single active layer, $27 \mu\text{m}$ thickness, comprising the active constituent, marker dye, stabilizers, and other admixtures that makes the film's energy dependence low. The yellow marker dye reduce the ultraviolet or light sensitivity and used together with an Red Green Blue film scanner which allows advantage of multichannel dosimetry. The active layer is placed between two, 120

μm thick transparent polyester substratums; the symmetric structure excludes the requirement of side selection of the film according to light source of the scanner. The polyester substratum has a particular surface consisting microscopic silica particles which creates a gap between the film surface and the glass window in a flatbed scanner. When this gap is about 10 times of visible light's wavelength it will prevent the flatbed scanner to form Newton's Rings interference patterns in acquired images.(49)

5.1.6 PTW Farmer Type Chamber TN30013

The 30013 Farmer chamber is the standard ionization chamber for absolute dose measurements in radiation therapy. Correction factors needed to determine absorbed dose to water or air kerma are published in the various dosimetry protocols. Its waterproof design allows the chamber to be used in water or solid state phantoms. The product's type is vented cylindrical ionization chamber, and its Sensitive Volume is 0.6 cm^3 . The type of central electrode is Aluminum, and the acrylic walls are made of graphite(50)(51).

5.1.7 RW3 Slab Phantoms

The slab phantom is designed for the use with photon radiation. The phantoms are used for monitor unit calibration and quality assurance measurements. Depth dose measurements are made by varying the measuring depth. The size of the complete phantoms is $30 \text{ cm} \times 30 \text{ cm} \times 30 \text{ cm}$. Adapter plates for some of the detector types are available for the Phantom. Each plate is precisely machined to a thickness tolerance of only $\pm 0.1 \text{ mm}$ (52).

5.2. METHOD

5.2.1 A CT Simulation of the Phantoms

Head of CIRS phantom was used in this study. The thermoplastic mask was used to eliminate potential submillimetric movement of head phantom during the scan. Also, Portrait Type-S™ mounted to couch to prevent movement during the procedure. Initial alignment was conducted by manually aligning the built-in planar positioning lasers, to intersect at the phantom's brain center (Photo.5.2.1) Once positioning was satisfactory, a full helical CT scan of the head at 120 kVp and 500 mA at 1,25 mm increments was performed. A total of 250 images were collected in a time of 46 seconds through a 17.5 mm cone beam, and the thickness of the slices are 1mm. 3D image is created from the reconstruction of x-ray images acquired during the simulation(53).




Photo.5.2.1. Phantom Position During CT Simulation

5.2.2 Exporting Data into the Treatment Planning System

The imaging DICOM data from the CT scan was then sent to the Varian Medical Systems, Inc. (Palo Alto, CA) Eclipse treatment planning system 13.0. This system allowed for the automatic identification of the scanned object in all three planes: axial, coronal, and sagittal. Once the appropriate body contour had been identified, the center of the dosimeter was set to the origin of the planes

5.2.3 Contouring

In the planning process of SRS, Gross Tumor Volume (GTV)/Clinical Target Volume (CTV) are determined by the physician. GTV was contoured in all cases and to account for set-up errors the 0.1cm margin was added to create PTV for simulating real patient situation. Because of single fraction and rigid frame fixation, uncertainties are thought to be close to zero in all our cases.

Brain part of Phantom was divided into three section in this study called Vertex (V), Nasal Cavity(NC) and Posterior Fossa(PF) regions. Vertex region was chosen because of the thinner region of cranium and beams passes very short distances to create conformal dose distribution. NC was chosen due to the heterogeneity of that region and to see the accuracy of the TPS for calculating dose distribution in heterogeneous medium. And at last the reason for choosing PF is the proximity of critical organs in that region, and some beams passes longer distances (especially in CK) to produce high conformity in the tumor. In every region, 4 PTV was contoured in Eclipse Treatment Planning System with the diameter of 3cm, 2cm, 1cm, 0.6cm as shown in photo (Fig 5.2.1-3).




Fig.5.2.3.1. 4 Different PTV Volumes in Vertex Region




Fig.5.2.3.2. 4 Different PTV Volumes in NC Region




Fig.5.2.3.3. 4 Different PTV Volumes in PF Region

5.2.4 Treatment Planning in Eclipse Treatment Planning System

12 PTV was created for 3 different anatomical sites in the brain. 5 Gy dose was prescribed for all tumors. Prescription dose was slightly lower than real SRS dose in order not to exceed linear dose response of EBT Gafchromic EBT3 films. 6FFF photon beam with 1400 MU/min dose rate was used in all plans. VMAT treatment delivery technique was employed in all TB plans. After arcs had been placed by designating the couch angles and gantry angles without the risk of collision, the jaw positions and collimator angles of each arc were adjusted to cover the combined PTV. To minimize the impact of interleaf leakage and tongue and Groove effect, the collimator angles were set to non-zero values and adjusted to

match the PTV position during rotation. Non-coplanar beams were used in V and NC regions. 3 Arcs were used for V and NC tumors. Vertex PTVs consist of focusing 2 full arcs with the collimator of 10°/80° and one partial non-coplanar arc 181°-0° with collimator 30° and the couch rotation of 90°. NC PTVs consist of focusing 2 full arcs with the collimator of 10°/80° and one partial non-coplanar arc 25°-181° with the collimator of 30° and the couch rotation of 90°. PF plans involve 2 full arcs with the collimator of 5°/85°. PRO algorithm was used in Rapid Arc optimization. Actual fluence map was determined by calculation all plans with AAA and AXB algorithm. Plan normalization was made to ensure $\geq 95\%$ of target volume received prescribed dose. All plans were calculated with the grid size of 0,1 mm. All SRS plans were evaluated by the radiation oncologist and clinical physicist. Plans were approved and prepared for delivery.

5.2.5 Treatment Planning in Multiplan

All CT data and RT structures contoured in Eclipse TPS were transferred to Multiplan TPS. Prescribed dose was 5 Gy for all plans. Treatment plans were done for same 12 PTVs as described above. 6FFF photon beam with 1000 MU/min dose rate was used in all plans. Skull tracking method was selected for treatment modality. Imaging center was aligned according to protocol. Fixed collimator was used for all plans. 30mm fixed collimator was utilized for all 3cm PTVs in V and NC and PF. 20 mm were used for all 2cm PTVs, and 7.5 mm fixed collimator were utilized for both 1cm and 0.6cm PTVs. Sequential Optimisation Algorithm was used for all calculations. Multiple non-coplanar beams were used in the treatment plan. A maximum number of beams for all plans were below than 70. Plan normalization was made to ensure $\geq 95\%$ of target volume received prescribed dose All plans were calculated using ray tracing algorithm and high-resolution (Fig 5.2.5.1). Then the high resolution plans were recalculated using the Monte Carlo algorithm with same beam parameters. (Fig 5.2.5.2). The CK-MLC plans were generated by using FSPB (Fig 5.2.5.3). All CK SRS plans were evaluated by the radiation oncologist and clinical physicist. Plans were approved and prepared for delivery.




Photo.5.2.5.1. Calculated dose distribution with Ray-Tracing Algorithm




Fig.5.2.5.2. Calculated dose distribution with Monte Carlo Algorithm




Fig.5.2.5.3. Calculated dose distribution with FSPB Algorithm

5.2.6 Preparation of Gafchromic EBT 3 Films

5.2.6.1 Scanner Orientation

The scan response of radiochromic films is sensitive to the orientation of the film on the scanner. This behavior results from the anisotropic scattering of the photons emitted by the scanner when passing through the polymer network, and the polarization of the transmit light by the needle-like shape particles of film active component that are preferentially aligned parallel to the direction in which the film was coated that is parallel to the short edge of the film(54). (Fig 5.2.6.1)




Fig.5.2.6.1. Scanner Orientation

5.2.6.2 Preparation and Irradiation of Films

Different dimensions of EBT3 films was used for film calibration and treatment delivery. Films with dimensions of $5 \times 20 \text{ cm}^2$ and $10 \times 12 \text{ cm}^2$ were used for film calibration and treatment delivery respectively (Photo 5.2.6.2). The EBT3 calibration films were oriented perpendicular to the central axis of the beam and irradiated, using 6 FFF photon beams of Varian TruebeamStx at the center of $10 \times 10 \text{ cm}^2$ field at 10cm depth and 100 cm SSD in a RW3Solid Water phantom ($30 \times 30 \times 20 \text{ cm}^3$)(55).




Photo.5.2.6.2. EBT3 Gafchromic film for calibration

The same setup was prepared for Cyberknife machine except for field size and SSD parameters. 60 cm fix collimator and SAD=80 cm were used for CK film calibration. 12 different MU (from 2 -1200 MU) values were used for both TB and CK film calibrations.

5.2.6.3 Scanning of EBT3 Films

To minimize the effect of the lateral dependence artifacts (the nonuniform response of the readout due to the light scattering of the scanner lamp caused by particles in the film active layer), a $10 \times 12.5 \text{ cm}^2$ cardboard template was fitted to the scanner to position films at a reproducible central location of the scan surface that can be considered uniform as shown Photo 5.2.6.3. All experimental, calibration and background films were scanned at the same location and orientation of an Epson Expression 11000 XL (America Inc., Long Beach, CA) flatbed scanner with reflection mode. All the EBT films were scanned at the same location and orientation that they irradiated. The settings of 48-bit color and 75dpi were used, color correction was disabled, and files were saved in TIFF format for the scanner. The reading and scanning of all the films occurred at least one day (24h) after the irradiation.

The calibration and background EBT films were grouped in strips, and each strip was scanned at the scanner bed center, one strip at a time. The calibration and background films were grouped together and scanned at the scanner bed center. The optical densities of exposed films were converted to the doses using the SNC Patient software program (sun nuclear), and the calibration curve was obtained for measurement(55).




Photo.5.2.6.3. Scanning of Films

5.2.7. Ion chamber measurement

PTW TN30013 (Farmer-type chamber PTW, Freiburg, Germany) ion chamber was inserted in the solid water phantom at the center of $10 \times 10 \text{ cm}^2$ field and 10 cm depth with SSD 100 cm to check the linac output before the irradiation. IAEA-TRS398 protocol was used for all calculations.

5.2.8 Set up of CIRS Phantom

For each irradiation gafchromic film was cut and inserted into V, NC and parts of the phantom, where the PTVs contoured at the geometric center of the target (Photo 5.2.8). Coordinates of films were decided according to treatment planning data.




Photo.5.2.8.1.position of film in vertex part of Phantom



Photo 5.2.8.2.Position of film on nasal cavity part of Phantom

The phantom was then placed on the LINAC treatment couch and aligned on the table in the exact orientation as it was scanned on the CT couch utilizing the red scribe lines.

5.2.8.1. Irradiation of Truebeamstx plans

After set up of the phantom in the Linac couch, KV images were taken from AP and LR axis. KV images were matched with DRR images and couch was shifted according to matching values. After KV imaging, Cone Beam CT (CBCT) were taken, and soft tissues were matched in all treatments. Plans were ready after CBCT matching.

5.2.8.2. Irradiation of MultiPlan TPS plans

After preparation of the phantom in the couch, KV images were obtained. Images were matched with skull DRR images and couch was shifted according to matching values in 6D directions. Before treatment, another image was taken and checked if any difference in phantom final position or not. The best alignment was done by shifts which were applied and was verified by the clinical physicist to start the irradiation. During the treatment, 2 more images were taken. Alignment was done with fewer images because phantom is stationary during treatment and more images could affect the result of Gafchromic EBT3.

5.2.8.3. Evaluation of irradiated Gafchromic EBT3 films

The films which were irradiated in Truebeam STx and CyberKnife, immediately removed after exposure and placed in the dark room and allowed to rest 24 hours before analysis for consistency. Then all the films were analyzed in Epson scanner according to the procedure followed in Film calibration. After scanning process for each film, a 'tif' file was created. These tif files were converted to '.flm.' files using SNC Patient software. Then dose distribution of slices, where films were placed according to, were imported from Eclipse and Multiplan to SNC patient software. Planned and measured dose was evaluated according to Gamma Index

analysis. Plan evaluation criteria in Gamma Index analysis was chosen as DTA=3mm and DD=3% (Fig 5.2.8.3).




Fig.5.2.8.3: Evaluation of an irradiated plan in SNC patient software; Upper left section represents dose distribution of irradiated plan, upper right section represents dose distribution of the slice of the irradiated plan, lower left section represents superposition of the both dose distributions with the best alignment, lower right section represents dose agreement graphically

6. RESULTS

6.1 Film Calibration Curve

In this study, Gafchromic EBT3 films were used to measure experimental dose distributions. Film calibration curves were measured and drawn for both TB and CK as shown Fig 6.1.1 and Fig 6.1.2 respectively. Twelve different MU (from 2 - 1200 MU) values were used for both TB and CK film calibrations. A same calibration curve was used for all film measurements.




Fig 6.1.1: Film Calibration Curve for Truebeamstx with 6MV FFF Beam




Fig 6.1.2: Film Calibration Curve for Cyberknife with 6 MV FFF Beam

6.2 Peripheral Dose Calculation Difference in TPS

In this study PD calculated by two different TPS were compared with EBT3 measured PD. In each planning system, 0.01 cc dummy volumes were created both left and right side in axial slices. These dummy volumes were drawn in the same plane of the geometrical center of PTV. The distances between PTV lateral borders and dummy volume centers were 2 cm on each side. The mean doses in dummy volume were defined as PD. The differences of PD between AAA and AXB algorithms for various tumor volumes and treatment region in the brain were shown in Table 6.2.1. Calculation differences between AAA and AXB for both V and PF regions were nearly similar for all tumor sizes. As the tumor size grows up, the differences between calculation algorithms decreases. The differences were more significant in NC region than others. The difference of PD between Ray-Tracing and Monte Carlo algorithms for various tumor volumes and treatment region in the brain were shown in Table 6.2.2. Calculated differences in Multiplan were slightly higher

than Eclipse. Calculation differences between Ray-Tracing and Monte Carlo were nearly similar except for NC region. NC region differences were quite higher than the others.

Table 6.2.1: Results of PD Differences Between AAA and AXB Algorithms in Eclipse TPS.

Comparison of AAA and AXB		2 cm Left Side	2 cm Right Side
		% Difference	%Difference
Vertex			
	0.6 cm	3%	3%
	1 cm	2%	2%
	2 cm	2%	2%
	3 cm	1%	1%
Nazal Cavite			
	0.6 cm	6%	5%
	1 cm	5%	4%
	2 cm	2%	3%
	3 cm	4%	4%
Posterior Fossa			
	0.6 cm	3%	2%
	1 cm	3%	2%
	2 cm	3%	3%
	3 cm	1%	1%

Table 6.2.2: Results of PD Differences Between Ray-Tracing and MC Algorithms in Multiplan TPS.

Comparison of Ray-Tracing and Monte Carlo		2 cm Left Side	2 cm Right Side
		% Difference	%Difference
Vertex			
	0.6 cm	6%	6%
	1 cm	6%	5%
	2 cm	3%	3%
	3 cm	3%	3%
Nazal Cavite			
	0.6 cm	15%	17%
	1 cm	13%	14%
	2 cm	13%	13%
	3 cm	12%	12%
Posterior Fossa			
	0.6 cm	7%	6%
	1 cm	5%	5%
	2 cm	5%	4%
	3 cm	2%	2%

6.3 Gamma Index Analysis for In field Volumes

V and NC regions plans were irradiated using Cirs Phantom and Gafchromic film. All films were analyzed using film QA program which is called “SNC Patient.” Gamma Index values for both regions and different tumor sizes were listed Table 6.3.1 and Table 6.3.2 for TB and CK. All plans were passed over $\geq 90\%$.

Table 6.3.1: Evaluation of Irradiated Volumes Using Gamma Index Score for Truebeam STx

TRUEBEAMSTX		GAMMA INDEX VALUES (DTA=3mm/DD=3%)	
		AAA	AXB
Vertex			
	0.6 cm	99%	97.60%
	1 cm	98.70%	99.00%
	2 cm	96.10%	97.00%
	3 cm	97.60%	96.40%
NC			
	0.6 cm	100%	99.00%
	1 cm	99.90%	96.60%
	2 cm	96%	99.40%
	3 cm	93.80%	94.90%

Table 6.3.2: Evaluation of Irradiated Volumes Using Gamma Index Score for Cyberknife

CYBERKNIFE		GAMMA INDEX VALUES (DTA=3mm/DD=3%)		
		RT	MC	FSPB
Vertex				
	0.6cm	99%	97.50%	98.80%
	1cm	97.50%	98.60%	99.20%
	2cm	94.50%	96.20%	95.0%
	3cm	99.90%	95.60%	95.90%
NC				
	0.6cm	99%	97.60%	99.20%
	1cm	95.80%	95.90%	96.20%
	2cm	92%	94.20%	92.10%
	3cm	94.70%	99.40%	92.90%

6.4. MU comparison in TB and CK

Type of treatment technique may effect calculated MU values. Calculated MU values for different tumor size and various tumor location in the brain for different types of treatment techniques and algorithms were shown in Table 6.4.1 and Table 6.4.2.

Table 6.4.1: Calculated MU values for different tumor size and various tumor location in the brain for AAA and AXB algorithms in Eclipse TPS

		AAA	AXB
Vertex			
	0.6 cm	1198	1170
	1 cm	1284	1268
	2 cm	1373	1368
	3 cm	1435	1399
Nazal Cavite			
	0.6 cm	1402	1362
	1 cm	1446	1414
	2 cm	1460	1433
	3 cm	1567	1519
Posterior Fossa			
	0.6 cm	1249	1216
	1 cm	1313	1303
	2 cm	1371	1338
	3 cm	1583	1530

Table 6.4.2: Calculated MU values for different tumor size and various tumor location in the brain for RT, MC and FSPB algorithms in Multiplan TPS

MU		RT	MC	FSPB
Vertex				
	0.6 cm	2106.4	2070.8	1184.2
	1 cm	1260	902.7	789.2
	2 cm	929.5	911	860.6
	3 cm	916.3	831.6	810.9
Nazal Cavite				
	0.6 cm	2347.5	2281.1	1230.3
	1 cm	1132.2	1008	859.7
	2 cm	945.8	925.6	827.5
	3 cm	886.9	882.1	752.8
Posterior Fossa				
	0.6 cm	1960.8	1940.7	1349.7
	1 cm	1338.6	1332.1	1109.2
	2 cm	1155.3	1132.7	1017.8
	3 cm	1021.9	999.4	920.3

6.5. Peripheral Dose Measurements in Truebeam STx

The PD distributions in V and NC regions were measured for all tumor sizes in TB. Although V was located in a homogeneous part of the brain, NC was located in the heterogeneous medium which contains bone, air, and soft tissue as shown Fig 6.5.1 PD distributions calculated by AAA and AXB algorithms in Eclipse TPS and measured in TB for V region are shown in Fig 6.5.2, Fig 6.5.3, respectively.




Fig 6.5.1: Tumors located in heterogeneous medium




Fig 6.5.2: Comparison of calculated PD with AAA in Eclipse and measured using TB for all tumor size in V region




Fig 6.5.3: Comparison of calculated PD with AXB in Eclipse and measured using TB for all tumor size in V region

PD distributions calculated by AAA and AXB algorithms in Eclipse TPS and measured in TB for NC region are shown in Fig 6.5.4 and Fig 6.5.5, respectively.




Fig 6.5.4: Comparison of calculated PD with AAA in Eclipse and measured using TB for all tumor size in NC region.




Fig 6.5.5: Comparison of calculated PD with AXB in Eclipse and measured Using TB for all tumor size in NC region

6.6. Peripheral Dose Measurements in CyberKnife M6

The PD distributions of V and NC regions were measured for all tumor sizes in CK. All plans were calculated with MC and RT with using the same beam sets and prescribed to same isodose level. MLC plans were calculated with completely different plan geometry with using FSPB algorithm. PD which were calculated using RT, MC and FSPB algorithms of Multiplan and measured in CK for V region are shown in Fig 6.6.1, Fig 6.6.2 and Fig 6.6.3 respectively.




Fig 6.6.1: Comparison of calculated PD with RT in Multiplan and measured using CK for all tumor size in V region




Fig 6.6.2: Comparison of calculated PD with MC in Multiplan and measured using CK for all tumor size in V region




Fig 6.6.3: Comparison of calculated PD with FSPB in Multiplan and measured using CK for all tumor size in V region

PD distributions of RT, MC, and FSPB algorithms were calculated by Multiplan and measured in CK for NC region are shown in Fig 6.6.4, Fig 6.6.5 and Fig 6.6.6, respectively.




Fig 6.6.4: Comparison of calculated PD with RT in Multiplan and measured using CK for all tumor size in NC region




Fig 6.6.5: Comparison of calculated PD with MC in Multiplan and measured using CK for all tumor size in NC region




Fig 6.6.6: Comparison of calculated PD with FSPB in Multiplan and measured using CK for all tumor size in NC region

6.7. Comparison of Peripheral Doses in Different Modalities

Treatment plannings of different modalities were compared in terms of PD. FSPB has the lowest PD distributions than the others, and both measurement and calculated (AAA and AXB) PD have the same results in V region. Some examples for different tumor sizes and locations was shown Fig 6.7.1, Fig 6.7.2, Fig 6.7.3 and Fig 6.7.4.




Fig 6.7.1 Calculated PD for small 0,6 cm tumor size in V region.




Fig 6.7.2 Measured PD for small 0,6 cm tumor size in V region.




Fig 6.7.3 Calculated PD for larger 2 cm tumor size in V region.




Fig 6.7.4 Measured PD for larger 2 cm tumor size in V region.

Left part of the tumor in NC was entirely filled with air because of its location in the brain. Calculated both AAA and AXB have same results in left and rights side of the tumor. Plans which were calculated with RT algorithms had lower

measurement results in the left side of the tumor. Some PD examples for different tumor sizes and locations was shown Fig 6.7.5, Fig 6.7.6, Fig 6.7.7 and Fig 6.7.8.




Fig 6.7.5 Calculated PD for 1 cm tumor size in NC region.




Fig 6.7.6 Measured PD for 1 cm tumor size in NC region.




Fig 6.7.7 Calculated PD for 2 cm tumor size in NC region.




Fig 6.7.8 Measured PD for 2 cm tumor size in NC region.

7. DISCUSSION

During radiotherapy treatments, the main concern is applying prescribe dose to target while sparing the dose constraints for organ at risk. The dose which is deposited outside the radiation field is very important in patient especially for whom receives SRS with long live expectancy(43)(56). The effect of absorbing lower doses has long been recognized and still is an issue of research to minimize long-term risk of secondary cancer(57)(58)

There are many studies defining PD resources in the literature (6,33,43,56). The main component of PD are:

- Leakage from treatment machine head
- Scattered radiation from machine components
- Scattered radiation within the patient

The first and second contributors depend on treatment devices and choice of treatment technique, but the last contributor depends on many parameters such as field size, radiation energy, etc. In this study, we concentrate on PD which is related to machine type and treatment technique.

Suresh R. et al reported the dose calculation accuracy of AXB and their results indicated that AXB was superior to AAA in predicting doses beyond air gaps as compared with measured data(59). Several studies showing that AXB was more accurate to use for dose computations in heterogeneous media than AAA by comparison against measurements or MC calculations. (31)(60)(61).

In this study, we calculate same treatment plans with both AAA and AXB for different tumor size and location in the brain. Both algorithms were in good agreement in V and PF locations which has homogeneous media inside the tumor. However, dose differences were quite large in NC location which has air cavity inside the tumor. These measurements were taken 2cm away from field edges as

shown in Table.6.2.1. The same comparison was done for RT and MC calculation algorithm in Multiplan. Large differences were observed between two algorithms as shown in Table.6.2.2. Especially in NC region at 2cm away from field edge, some differences were higher than 10% which is so high for clinical tolerances.

Ning W. et al showed the precision of the film dosimetry for SRS and radiotherapy using Gafchromic EBT3(62). In our study, Gafchromic film was used for verification of plans. Gamma index analysis was done by using $\pm 3\%$ DTA criteria and all results were shown in Table.6.2.1 and Table.6.2.2. All measurement Gamma index values were in good agreement both in AXB and AAA algorithm for V and NC regions. RT and FSPB algorithms showed lower Gamma index values in comparison with MC in NC region but a good agreement was observed in V region.

G Kragl. et al reported that PD can not be calculated easily with a high degree of accuracy(63). With the absence of FF, the major source of scattered radiation is removed, resulting in low dose exposure to the patient's normal structures for linear accelerator treatments(63). Compare to the conformal radiotherapy, SRS treatments utilize IMRT or VMAT techniques which have multiple beam entrance to achieve highly conformal dose distributions. These techniques include beamlets inside the beam which were increased the number of MU and beam on times. Both two factors increase normal tissue low dose values. At the edge of treatment field, PD is predominantly created by scattered photons. With increasing distance from field edge patient scatter dominates and leakage through to the linac head starts to contribute PD. The long distance from field edge leakage radiation becomes the major source of PD(63)(64)(65). If beam energies above 10MV, neutrons start to contribute PD(66).In our study, we calculated and measured TB PD for AAA and AXB algorithms in V and NC regions which are shown in Fig.6.5.2-Fig.6.5.5. AAA and AXB were in good agreement in both homogeneous(V) and heterogeneous(NC) regions. Our study was concentrated on PD at the field edge. PD were related directly to the field size. Large tumor volume irradiations increase PD. In CK treatments commonly arbitrary beam directions are used. Because of these pattern of the machine, entrance or exit dose from beams oriented from inferior or superior

oblique angles with respect to the patient can also contribute to PD near to the field edge(43). In our study, calculated and measured PDs at the edge of the field border were in good agreement both for MC and FSPB algorithms in V region. Same results were observed in NC region. We found negligible differences between measured and calculated doses in V region for RT algorithm. Differences became significant in NC region for RT algorithm. Film measurements were lower than calculated in both V and NC region for RT algorithm. We observed that increasing field size also increases PD for CK treatment as we mentioned before for TB.

Petti showed that the PD is mainly proportional to delivered MU rather than delivered dose(43). Total MU values were nearly same in both AAA and AXB algorithms in Eclipse TPS as shown Table 6.4.1. MLC plans calculated with FSPB algorithm were completely separate plan from RT and MC plans. The total MU values were very low in FSPB algorithm causing much lower PD in CK as shown 6.4.2. Reducing MU values, a number of beam entrance and gradient index value always decrease PD both in CK and TB treatments.

Petti's data suggested that CK PD is larger than both gamma Knife and IMRT treatments(43). In our study, we obtained same results in plans which were calculated using RT algorithm as shown in Fig.6.7.1-8. The plans calculated with RT had larger PD than the plans calculated with MC and FSPB in V and NC regions. PD ratios were directly related to treatment techniques. Non-coplanar arbitrary beam treatments had higher PD in all plans.

8. CONCLUSION

Many factors affect PD in modern radiotherapy techniques. Machine configuration is one of the major factors which effects directly PD in the patient. PD related to machine configuration cannot be blocked but changing some other factor such as treatment technique, algorithm, heterogeneity in the field, the number of MU value etc. could change the value of PD. In this study, we have represented that the changes in treatment technique and used various algorithm may cause large differences in PD.



9. REFERENCES

1. Slater JM. Ion Beam Therapy. 2012;320:3–17. Available from: <http://link.springer.com/10.1007/978-3-642-21414-1>
2. Bernier J, Hall EJ, Giaccia A. Radiation oncology: a century of achievements. *Nat Rev Cancer*. 2004;4(9):737–47.
3. Frankel RI. Centennial of Röntgen's discovery of x-rays. *West J Med* [Internet]. 1996;164(6):497–501. Available from: <http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=1303625&tool=pmcentrez&rendertype=abstract>
4. Baskar R, Lee KA, Yeo R, Yeoh KW. Cancer and radiation therapy: current advances and future directions. *Int J Med Sci* [Internet]. 2012;9(3):193–9. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/22408567>
5. Lomax NJ, Scheib SG. Quantifying the degree of conformity in radiosurgery treatment planning. *Int J Radiat Oncol Biol Phys*. 2003;55(5):1409–19.
6. Fraass B, van de Geijn J. Peripheral dose from megavolt beams. *Med Phys* [Internet]. 1983;10(6):809–18. Available from: <http://scitation.aip.org/content/aapm/journal/medphys/10/6/10.1118/1.595359>
7. Modalities T. Introduction to Stereotactic Radiation (Radiosurgery / SRS). 2017;
8. Kirkpatrick CJ, Wagner M, Köhler H, Bittinger F, Otto M, Klein CL. The cell and molecular biological approach to biomaterial research: A perspective. *J Mater Sci Mater Med*. 1997;8(3):131–41.

9. Song CW, Park H, Griffin RJ, Levitt SH. Radiobiology of Stereotactic Radiosurgery and Stereotactic Body Radiation Therapy. 2012;51–61.
10. Brenner DJ. The Linear-Quadratic Model Is an Appropriate Methodology for Determining Isoeffective Doses at Large Doses Per Fraction. *Semin Radiat Oncol.* 2008;18(4):234–9.
11. Tryggestad E, Christian M, Ford E, Kut C, Le Y, Sanguineti G, et al. Inter- and intrafraction patient positioning uncertainties for intracranial radiotherapy: A study of four frameless, thermoplastic mask-based immobilization strategies using daily cone-beam CT. *Int J Radiat Oncol Biol Phys.* 2011;80(1):281–90.
12. Cerviño LI, Pawlicki T, Lawson JD, Jiang SB. Frame-less and mask-less cranial stereotactic radiosurgery: a feasibility study. *Phys Med Biol.* 2010;55(7):1863–73.
13. Brown JM, Giaccia AJ. The Unique Physiology of Solid Tumors : Opportunities (and Problems) for Cancer Therapy The Unique Physiology of Solid Tumors : Opportunities (and Problems) for Cancer Therapy1. *Cancer Res.* 1998;58(9):1408–16.
14. Paddick I, Lippitz B. A simple dose gradient measurement tool to complement the conformity index. *J Neurosurg.* 2006;105 Suppl:194–201.
15. Feuvret L, Noël G, Mazeron JJ, Bey P. Conformity index: A review. Vol. 64, *International Journal of Radiation Oncology Biology Physics.* 2006. p. 333–42.
16. Wagner TH, Bova FJ, Friedman WA, Buatti JM, Bouchet LG, Meeks SL. A simple and reliable index for scoring rival stereotactic radiosurgery plans. *Int J Radiat Oncol Biol Phys.* 2003;57(4):1141–9.

17. Kataria T, Sharma K, Subramani V, Karrthick KP, Bisht SS. Homogeneity Index: An objective tool for assessment of conformal radiation treatments. *J Med Phys* [Internet]. 2012;37(4):207–13. Available from: <http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=3532749&tool=pmcentrez&rendertype=abstract>
18. Yin F-F, Wong J, Balter J, Benedict S, Bissonnette J-P, Craig T, et al. The Role of In-Room kV X-Ray Imaging for Patient Setup and Target Localization Report of AAPM Task Group 104. *Data Management*. 2009.
19. Bissonnette J-P, Balter P a., Dong L, Langen KM, Lovelock DM, Miften M, et al. Quality assurance for image-guided radiation therapy utilizing CT-based technologies: A report of the AAPM TG-179. *Med Phys*. 2012;39(4):1946.
20. Srivivasan K, Mohammadi M, Shepherd J. Applications of linac-mounted kilovoltage Cone-beam Computed Tomography in modern radiation therapy: A review. *Polish J Radiol*. 2014;79:181–93.
21. Benedict SH, Yenice KM, Followill D, Galvin JM, Hinson W, Kavanagh B, et al. Stereotactic body radiation therapy: The report of AAPM Task Group 101. *Med Phys* [Internet]. 2010;37(8):4078. Available from: <http://scitation.aip.org/content/aapm/journal/medphys/37/8/10.1118/1.3438081>
22. Shirato H, Shimizu S, Kunieda T, Kitamura K, Van Herk M, Kagei K, et al. Physical aspects of a real-time tumor-tracking system for gated radiotherapy. *Int J Radiat Oncol Biol Phys*. 2000;48(4):1187–95.
23. Shende R, Gupta G, Patel G, Kumar S. Assessment and performance evaluation of photon optimizer (PO) vs . dose volume optimizer (DVO) for IMRT and progressive resolution optimizer (PRO) for RapidArc planning using a virtual phantom. 2016;4(3).

24. Varian Medical Systems OS. Eclipse Algorithms Reference Guide Eclipse. 2011;(December):1–368.
25. Vanetti E, Nicolini G, Nord J, Peltola J, Clivio A, Fogliata A, et al. On the role of the optimization algorithm of RapidArc® volumetric modulated arc therapy on plan quality and efficiency. *Med Phys*. 2011;38(2011):5844.
26. Cozzi L, Dinshaw KA, Shrivastava SK, Mahantshetty U, Engineer R, Deshpande DD, et al. A treatment planning study comparing volumetric arc modulation with RapidArc and fixed field IMRT for cervix uteri radiotherapy. *Radiother Oncol* [Internet]. 2008;89(2):180–91. Available from: <http://dx.doi.org/10.1016/j.radonc.2008.06.013>
27. Sievinen J, Ulmer W, Kaissl W. AAA photon dose calculation model in Eclipse. Palo Alto Varian Med Syst [Internet]. 2005;Varian doc:1–23. Available from: http://www.rtsalon.cn/upload/RTsalon_p_3218_2.pdf
28. Tillikainen L, Helminen H, Torsti T, Siljamäki S, Alakuijala J, Pyyry J, et al. A 3D pencil-beam-based superposition algorithm for photon dose calculation in heterogeneous media. *Phys Med Biol*. 2008;53(14):3821–39.
29. Kielar KN, Mok E, Hsu A, Wang L, Luxton G. Verification of dosimetric accuracy on the TrueBeam STx: rounded leaf effect of the high definition MLC. *Med Phys* [Internet]. 2012;39(10):6360–71. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/23039672>
30. Hrbacek J, Lang S, Klöck S. Commissioning of photon beams of a flattening filter-free linear accelerator and the accuracy of beam modeling using an anisotropic analytical algorithm. *Int J Radiat Oncol Biol Phys*. 2011;80(4):1228–37.

31. Vassiliev ON, Wareing TA, McGhee J, Failla G, Salehpour MR, Mourtada F. Validation of a new grid-based Boltzmann equation solver for dose calculation in radiotherapy with photon beams. *Phys Med Biol* [Internet]. 2010;55(3):581–98. Available from: <http://stacks.iop.org/0031-9155/55/i=3/a=002?key=crossref.32a9435d5666fc4f24d83ce44bba291e%5Cnhttp://www.ncbi.nlm.nih.gov/pubmed/20057008%5Cnhttp://iopscience.iop.org/0031-9155/55/3/002>

32. Failla GA, Wareing T, Archambault Y, Thompson S. Acuros ® XB advanced dose calculation for the Eclipse TM treatment planning system. 2010;32. Available from: https://www.varian.com/sites/default/files/resource_attachments/AcurosXBClinicalPerspectives_0.pdf

33. Vlachopoulou V, Antypas C, Delis H, Tzouras A, Salvaras N, Kardamakis D, et al. Peripheral doses in patients undergoing Cyberknife treatment for intracranial lesions. A single centre experience. *Radiat Oncol* [Internet]. 2011;6(1):157. Available from: <http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=3228676&tool=pmcentrez&rendertype=abstract%5Cnhttp://www.royalsocietypublishing.org/journal/ro/content/6/1/157>

34. Accuray. *Physics Essentials Guide*. 2015;1–678.

35. Accuracy Incorporated. *CyberKnife Equipment Specification*. 2009;1–36. Available from: <papers2://publication/uuid/D9B2A0ED-A64D-4356-882C-FB482F6798AC>

36. Bailey S, Calculations B. *Treatment Planning – Manual Calculations*. 2015;301:8–11.

37. Glide-Hurst C, Bellon M, Foster R, Altunbas C, Speiser M, Altman M, et al. Commissioning of the Varian TrueBeam linear accelerator: a multi-institutional study. *Med Phys* [Internet]. 2013;40(3):31719. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/23464314>
38. Chang Z, Wu Q, Adamson J, Ren L, Bowsher J, Yan H, et al. Commissioning and dosimetric characteristics of TrueBeam system: Composite data of three TrueBeam machines. *Med Phys* [Internet]. 2012;39(11):6981. Available from: <http://scitation.aip.org/docserver/fulltext/aapm/journal/medphys/39/11/1.4762682.pdf?expires=1447859013&id=id&acname=2112784&checksum=0F210486D1EBC2F1AF5C5A8DD81DFEF6>
39. Korreman SS, Juhler-Nøttrup T, Persson GF, Navrsted Pedersen A, Enmark M, Nystrom H, et al. The role of image guidance in respiratory gated radiotherapy. *Acta Oncol*. 2008;47(7):1390–6.
40. Zytkovicz a, Daftari I, Phillips TL, Chuang CF, Verhey L, Petti PL. Peripheral dose in ocular treatments with CyberKnife and Gamma Knife radiosurgery compared to proton radiotherapy. *Phys Med Biol*. 2007;52:5957–71.
41. Kilby W, Dooley JR, Kuduvalli G, Sayeh S, Maurer CR. The CyberKnife Robotic Radiosurgery System in 2010. *Technol Cancer Res Treat* [Internet]. 2010;9(5):433–52. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/20815415>
42. McGuinness CM, Gottschalk AR, Lessard E, Nakamura JL, Pinnaduwege D, Pouliot J, et al. Investigating the clinical advantages of a robotic linac equipped with a multileaf collimator in the treatment of brain and prostate cancer patients. *J Appl Clin Med Phys*. 2015;16(5):284–95.

43. Petti PL, Chuang CF, Smith V, Larson D a. Peripheral doses in CyberKnife radiosurgery. *Med Phys.* 2006;33(6):1770–9.
44. Blackwell CR, Coursey BM, Gall KP, Galvin JM, Mclaughlin WL, Meigooni AS, et al. Radiochromic Film Dosimetry. 1998;25(63).
45. Butson MJ, Yu PKN, Cheung T, Metcalfe P. Radiochromic film for medical radiation dosimetry. 2003;41:61–120.
46. Williams MJ, Metcalfe PE. Radiochromic film dosimetry and its applications in radiotherapy. 2011;75–99.
47. Jordan K. Review of recent advances in non 3D dosimeters.
48. Papaconstadopoulos P, Hegyi G, Seuntjens J, Devic S, Papaconstadopoulos P, Hegyi G, et al. A protocol for EBT3 radiochromic film dosimetry using reflection scanning A protocol for EBT3 radiochromic film dosimetry using reflection scanning. 2014;122101.
49. Media D. Dosimetry media, type ebt-3.
50. Search H, Journals C, Contact A, Iopscience M, Address IP. Determination of the recombination correction factor k_s for some specific plane-parallel and cylindrical ionization chambers in pulsed photon and electron. 2007;35.
51. RADIATION PHYSICS | Thimble Ionization Chambers Model N 30013 Waterproof PTW Farmer Ionization Chamber. :37217.
52. Slab Phantom. 2011;4182011.
53. Cagnon C. from *Computed Tomography Advances in Technology ...* 2008;

54. Christopher G. Radiochromic Film.
55. Lewis D. Correcting lateral response artifacts from flatbed scanners for radiochromic film dosimetry. 2015;42(1):416–29. Available from: <http://dx.doi.org/10.1118/1.4903758>
56. Betta E Di, Fariselli L, Vecchio A Del, Broggi S. Evaluation of the peripheral dose in stereotactic radiotherapy and radiosurgery treatments a 2010;3587–94.
57. Risk THE, Health TO, Exposure F, Low TO, Radiation OFI. Guest Editorial.
58. Brenner DJ, Doll R, Goodhead DT, Hall EJ, Land CE, Little JB, et al. Cancer risks attributable to low doses of ionizing radiation : Assessing what we really know. 2003;100(24):3–8.
59. Dosimetric evaluation of Acuros XB dose calculation algorithm with measurements in predicting doses beyond different air gap thickness for smaller and larger field sizes.
60. Bush K, Gagne IM, Zavgorodni S, Ansbacher W, Beckham W, Alto P. Dosimetric validation of Acuros ® XB with Monte Carlo methods for photon dose calculations. 2011;2208–21.
61. Fogliata A, Nicolini G, Clivio A, Vanetti E, Cozzi L. Dosimetric evaluation of Acuros XB Advanced Dose Calculation algorithm in heterogeneous media. 2011;1–15.
62. Wen N, Lu S, Kim J, Qin Y, Huang Y, Zhao B, et al. Precise film dosimetry for stereotactic radiosurgery and stereotactic body radiotherapy quality assurance using Gafchromic TM EBT3 films. *Radiat Oncol* [Internet]. 2016;1–11. Available from: <http://dx.doi.org/10.1186/s13014-016-0709-4>

63. Kragl G, Baier F, Lutz S, Albrich D, Dalaryd M, Kroupa B, et al. Flattening filter free beams in SBRT and IMRT : Dosimetric assessment of peripheral doses. 2011;21:91–101.
64. Kry SF, Vassiliev ON, Mohan R. Out-of-field photon dose following removal of the. 2155.
65. 1.597525.pdf.
66. Ry STFK, Alehpour MOS, Ollowill DASF, Osen ISIR, Physics R, Oncology R, et al. OUT-OF-FIELD PHOTON AND NEUTRON DOSE EQUIVALENTS FROM STEP-AND-SHOOT INTENSITY-MODULATED RADIATION THERAPY M ARILYN S TOVALL , P H . D .,* D EBORAH A . K UBAN , M . D ., † R . A LLEN W HITE , P H . D ., ‡. 2005;62(4):1204–16.

10. CURRICULUM VITAE

Kişisel Bilgiler

ADI	Hanieh	Soyadı	MABHOUTİ
-----	--------	--------	----------

Eğitim Düzeyi

	Mezun Olduğu Kurumun Adı	Mezuniyet Yılı
Doktora/Uzmanlık		
Yüksek Lisans	Tabriz Özel Bilim ve Araştırma Üniversitesi	2013
Lisans	Mahabad Özel Üniversitesi Fen Edebiyat Fakültesi Fizik Bölümü	2011
Lise	Niyayeş Lisesi	2006

İş Deneyimi (Sondan geçmişe doğru sıralayın)

Görevi	Kurum	Süre (Yıl-Yıl)
Medikal Fizik Stajı	Medipol Mega Üniversite Hastanesi Radyasyon Onkolojisi Bölümü	2014-2017

Yabancı Dilleri	Okuduğunu Anlama	Konuşma	Yazma	KPDS Puanı	YDS Puanı
İngilizce	iyi	iyi	iyi		
Almanca	Orta	Orta	Orta		

Bilgisayar Bilgisi

Microsoft Office Programları	Orta