

FOR REFERENCE

COMPUTER AIDED NEUTRON ACTIVATION ANALYSIS TAKEN FROM THIS ROOM

by

HALUK BÖGRÜN

B.S. in Mechanical Engineering
Middle East Technical University
1978

Bogazici University Library



Submitted to the Faculty of the School of Engineering
in Partial Fulfillment of
the Requirements for the Degree of
Master of Science
in

Nuclear Engineering


BOĞAZIÇI UNIVERSITY


1982


APPROVAL OF READERS


This is to certify that I have read and approved the thesis entitled "A Computer Aided Neutron Activation Analysis", submitted by Halûk Böğrün in partial fulfillment of the requirements for the degree of Master of Science in Nuclear Engineering, School of Engineering, Boğaziçi University.

Date :

1- Dr. Vural ALTIN 

2- Doç. Dr. Fahir I. BORAK (Thesis Advisor) 

3- Prof. Dr. Turan ENGİNOĞLU (Department Chairman) 

4- Doç. Dr. İlsen ÖNSAN 

ABSTRACT

Neutron Activation Analysis (NAA) has become a widely used and important analytical technique in the determination of the minor constituents of many samples. The composition of these samples may range from macro-quantities (greater than 1%) down to micro-quantities (parts per million). The micro qualitative and quantitative analyses can be performed by NAA. For this purpose the basic concepts and the theory of NAA is introduced. The detection methods and various types of detectors are discussed. The auxiliary electronic equipments and the statistical methods used in the evaluation of the results are also given. A pre-developed computer code, namely CORGAM, is presented in order to unfold the gamma-ray spectra obtained from the activation analysis of mixtures of materials. The code which was previously programmed for the IBM 360/50 computer is updated for use in the BOĞAZIÇI University computer UNIVAC-1106. The data from a sample problem are used as input to run the program in B.U. computer. The output obtained from this run is observed to fit exactly to the results previously obtained. Finally the program and the input library are recorded onto a magnetic tape, and the necessary information required to run the program is given.

ÖZET

Günümüzde Nötron Aktivasyon Analizi (NAA) yolu ile, çeşitli malzeme kompozisyonlarında, malzemeyi teşkil eden elementlerin cins ve miktarlarının tayini oldukça önem kazanmıştır. Bu çalışmada, NAA'nın dayandığı prensipler, genel teorisi, analiz sırasında kullanılan detektör sistemleri incelenmiş ve neticelerin istatistikî değerlendirilmesi için kullanılan teknikler tanıtılmıştır. Diğer analitik yöntemlere oranla daha hassas sonuç verme üstünlüğüne sahip olan nötron aktivasyon analizinde sonuçların değerlendirilmesi sırasında karşılaşılan uzun zaman alan hesaplamalar, bilgisayarların hızlı hesaplama yeteneğinden yararlanılarak yapılmaya çalışılmıştır. Bu amaçla, daha önce hazırlanmış olan bir bilgisayar programından yararlanılmış ve bu program (CORGAM) üzerinde gerekli değişiklikler yapılarak sadece tasarlandığı üzere IBM 360/50 tipi bilgisayarlarda değil UNIVAC-1106 tipi bilgisayarlarda da kullanılabilmesi sağlanmıştır. Program değiştirilmiş hali ile manyetik bir banda kaydedilmiş ve gelecekteki uygulamalar için gerekli bilgiler ayrıca sunulmuştur.

ACKNOWLEDGMENT

The author would like to express his gratitude to Doç.Dr. Fahir I. Borak who, as his thesis advisor, offered invaluable guidance and assistance. Special thanks are extended to Prof. Dr. Turan Enginol, Head of Nuclear Engineering Department, and to Dr. Vural Altın for their constructive criticisms during this work. A very special expression of love and gratitude is extended to the author's wife because of the deep devotion, encouragement, and uncomplaining sufferings throughout this work.

TABLE OF CONTENTS

| | Page |
|----------------------------------------------------------------------------------------------|------|
| ABSTRACT | i |
| ÖZET | ii |
| ACKNOWLEDGMENT | iii |
| LIST OF TABLES | v |
| LIST OF FIGURES | vi |
| 1. INTRODUCTION | |
| 1.1 GENERAL | 1 |
| 1.2 NEUTRON ACTIVATION ANALYSIS | 3 |
| 1.3 BASIC PRINCIPLES OF GAMMA-RAY SPECTROMETRY | 5 |
| 1.4 OUTLINE OF THE THESIS | 7 |
| 2. GENERAL PRINCIPLES | |
| 2.1 FORMATION OF RADIOACTIVE NUCLEI THROUGH THE INTERACTION OF NEUTRONS WITH MATTER | 8 |
| 2.2 INTERACTION OF GAMMA-RAYS WITH MATTER AND GAMMA-RAY SPECTRUM | 12 |
| 2.3 SCINTILLATION DETECTORS | 21 |
| 2.4 AUXILIARY ELECTRONIC EQUIPMENT | 28 |
| 2.5 STATISTICAL INTERPRETATION OF THE RESULTS | 30 |
| 2.6 USE OF COMPUTERS IN NEUTRON ACTIVATION ANALYSIS | 35 |
| 3. CORGAM | |
| 3.1 A COMPUTER CODE FOR NEUTRON ACTIVATION ANALYSIS | 37 |
| 3.2 THEORY OF CORGAM | 39 |
| 3.3 EXPERIMENTAL PROCEDURE OF CORGAM | 46 |
| 4. DISCUSSION AND CONCLUSION | |
| REFERENCES | 52 |
| APPENDICES | |
| A. PROPERTIES OF ESTIMATORS | 55 |
| B. THE LEAST SQUARES METHOD | 57 |
| C. CURVE FITTING METHODS | 59 |
| D. LAGRANGIAN INTERPOLATION | 62 |
| E. DETERMINANT EVALUATION BY CHIO'S METHOD | 64 |
| F. DESCRIPTION OF THE PROGRAMS | 65 |

| | | |
|----|--------------------------------------------------------|-----|
| G. | LISTING OF THE INPUT DATA AND INPUT OPTIONS | 77 |
| H. | LISTING OF THE UNIVAC VERSION OF THE CODE (CORGAM) ... | 84 |
| J. | FLOW DIAGRAM OF THE CODE | 101 |
| K. | LISTING OF THE INPUT DATA FOR SAMPLE PROBLEM | 104 |
| L. | LISTING OF THE OUTPUT OF THE SAMPLE PROBLEM | 115 |
| M. | LISTING OF THE UPDATING PROGRAMS | 148 |
| N. | COMPUTER CODES LISTING | 152 |

LIST OF TABLES

| Table | Page |
|--------------------------------------------------------------------------------------------|------|
| 2.1a Calibration sources for gamma-ray spectrometry..... | 19 |
| 2.1b Energies and relative intensities in the decay of 1200 y $^{166}\text{Ho}^m$ | 20 |
| 2.2 Properties of some commonly used phosphors | 22 |
| 2.3 Characteristics of a few commercial photomultipliers . | 26 |
| 2.4 Probability levels | 31 |
| 3.1 Compilation of nuclear properties of several isotopes. | 40 |
| 3.2 Neutron flux levels for KSU TRIGA Mark II reactor..... | 47 |
| 3.3 Reference spectra information | 48 |

LIST OF FIGURES

| Figure | Page |
|------------------------------------------------------------------------------------------|------|
| 1.1 Comparison of a theoretical and actual response of a scintillation spectrometer..... | 6 |
| 2.1 Decay scheme of ^{60}Co | 9 |
| 2.2 Decay scheme of ^{137}Ce | 10 |
| 2.3 Photoelectric cross section as a function of photon energy for lead | 13 |
| 2.4 Schematic diagram of Compton scattering | 14 |
| 2.5 Absorption coefficient for NaI(Tl) as a function of gamma-ray energy | 17 |
| 2.6 Gamma spectrum of ^{24}Na measured with a NaI(Tl) | 18 |
| 2.7 PM tube (Venetian-blind dynode structure) | 24 |
| 2.8 PM tube (Box type dynode structure) | 24 |
| 2.9 PM tube (Circularly focussed dynode structure) | 25 |
| 2.10 Pulse-height spectrum of ^{137}Cs | 25 |
| 2.11 Spectrometer system | 29 |
| 2.12 Pulse-height to time conversion | 29 |
| 2.13 The normal distribution function | 32 |
| 2.14 Signal detection and critical level | 33 |
| 3.1 Outline of correlation algorithm | 45 |
| 3.2 Shielded counting cavity dimensions | 49 |
| 3.3 Block diagram of the NAA instrumentation system | 47 |

1. INTRODUCTION

1.1- GENERAL

The nuclei of atoms are stable only when they contain a certain number of neutrons and protons. The number of neutrons in the nucleus of an atom determines whether that atom is radioactive or stable. (There are exceptions. For a few elements there are no stable nuclei. In some cases, there are other differences that make certain nuclei radioactive.) Thus, while all sodium atoms contain 11 protons, only those sodium atoms that contain 12 neutrons are stable. For other elements there may be more than one number of neutrons that results in stability; for instance there are 10 stable atoms (isotopes) of tin; each containing a different number of neutrons in their nuclei.

The fact that nuclei can absorb additional neutrons, which, in many cases, results in the conversion of a stable nucleus to a radioactive one, makes Neutron Activation Analysis (NAA) possible. Because radioactive nuclei decay in unique ways and yield radiations that are often distinct, and can be measured even in very small amounts, measurements of these radiations can determine the kind and the number of radioactive atoms that are present.

In order to perform NAA, the first condition is that the element of interest gives by some nuclear reaction a radioactive isotope. The probability of the reaction (cross-section), the isotopic abundance of the target nucleus, and the half-life of the formed isotope should be large enough to allow measurements of the emitted radiation. The type and the energy of this radiation is also of interest, mainly in overcoming nuclear interferences induced by the matrix and the other impurities.

Through its nuclear character, the neutron activation technique allows the analysis of a wide variety of elements in widely differing matrices, making it a powerful analytical tool not only in trace analysis but also in determining alloy elements and minor constituents. Therefore the field of application of the method is very wide. It has been used to determine traces in ultra pure elements, in biological materials and minerals, in analysis where only small samples are available, in identification of industrial products by trace characterization, in industrial analysis of steel, in geochemistry for the determination of isotopic abundancies, etc. (Al-61, Ch-60, Pi-62)

Although NAA was proposed by von Hevesey and Levi (He-36) followed in 1938 by charged particle activation analysis applied by Seaborg et al. (Se-38) the method was considered as rather a curiosity until mid-forties. Micro constituents analysis by NAA, however became possible only after the development of high flux reactors, accelerators and reliable scintillation spectrometers. Basic principles of NAA are discussed in detail by many authors (So-69, Ba-64, Ev-55, Pr-64, OK-72). In addition some application of digital computers to NAA were presented at a symposium in 1962 (OK-62).

1.2- NEUTRON ACTIVATION ANALYSIS

The irradiation of most atomic nuclei by thermal neutrons causes these nuclei to become radioactive. The subsequent emission of beta or gamma radiations from these nuclei is characteristic of the particular isotope. If the energy of emission of the radiations are known, the isotope can be identified. The magnitude of the radioactivity is directly related to the amount of the isotope present in the specimen. The detection of radioactivity requires sensitive techniques. However neutron activation analysis permits the accurate appraisal of very small amounts of an element. Neutrons are able to pass through a considerable thickness of most substances. It is possible, therefore to analyze materials in solid, liquid or gaseous phases (Wa-63).

Besides the high sensitivity of thermal neutron activation analysis, this technique offers a number of other advantages such as the simplicity of the nuclear reaction, transparency of most materials to thermal neutrons and non-destruction of the sample, etc. In principle this means that in most cases only neutron capture reactions have to be considered and that consequently there is no change in the nature of the irradiated elements.

Unfortunately thermal NAA can not solve a number of important trace analysis problems, as for instance those of light elements, e.g. carbon, nitrogen and oxygen. These problems can sometimes be solved by applying an other method of activating particles: photons or charged particles such as protons, deuterons, helium ions, etc. If neutron analysis is possible for the whole mass of the sample, this is not true for charged particle analysis as these have a very short range. While all analytical techniques require some chemical treatment before the analysis, the NAA may require the chemical treatment only after irradiation and often after removal of surface contamination. One further advantage of NAA is in many cases the selectivity or even specificity of its nature. It should be remembered that the activation is a reaction of the nucleus, and thus independent of the chemical state of the considered element (So-69).

By applying gamma scintillation spectrometry to activation analysis it becomes possible to analyze simultaneously a number of different elements. A careful statistical treatment can demonstrate that NAA can be as precise and even better than other analytical methods.

The procedure can be summarized as follows: irradiate a sample and a standard under the same conditions; count both under the same conditions; the ratio of these rates will be equal to the ratio of the masses of the element of interest. If one takes activation in a thermal reactor, one should remember among many other parameters that the neutron energy varies from a thermal energy of 0.025 eV up to about 20 MeV. Consequently activation not only occurs with thermal neutrons, but also with epithermal and fast neutrons. Many other complexities may be met during activation, such as unsuspected reactions, production of the isotope of interest by second order reactions or radioactive growth of a daughter isotope, irreproducible geometry, absorption of radiation, instability of the instrumentation, non-linearity in detector response, etc. From this short list of errors it can be concluded that NAA is a complex analytical technique and that experience is required to obtain satisfactory results.

1.3- BASIC PRINCIPLES OF GAMMA-RAY SPECTROMETRY

The basic principle of NAA is the determination of the activity of the product nuclei Y from the activation of the parent nuclei X by means of neutrons. The nuclear reaction of this type can be represented in short notation, as follows



where n refers to neutron and k to the emitted particles or radiation during the reaction. If the activation process of nuclei X, is assumed to be carried out in constant neutron flux, and if the initial number of atoms of the parent nuclei is also assumed to be constant, then the well known differential equation governing the rate of formation of Y will have the following solution, as given in detail in section 3.2.

$$A_y \equiv N_y \lambda_y = \sigma \phi_0 N_x (1 - e^{-\lambda_y t}) \quad (1.2)$$

where, A_y is the activity of Y, atoms/sec.,

N_y is the number of atoms of Y,

N_x is the number of atoms of X,

λ_y is the decay constant for product Y, sec.⁻¹

ϕ_0 is the neutron flux, neutrons/cm² sec,

σ is the activation cross-section for production of Y from X, cm²,

t is the time of irradiation, sec.

When a sample is removed from the neutron flux, the activity begins to decrease exponentially. Then the equation 1.2 becomes,

$$A_y \equiv N_y \lambda_y = \sigma \phi_0 N_x (1 - e^{-\lambda_y t}) e^{-\lambda_y t^0} \quad (1.3)$$

where, t^0 is the corrected decay time to account for decay of the sample during a finite counting time, and is given by;

$$t^0 = t_r - \frac{t_{1/2}}{0.693} \left[\frac{t_{1/2} (1 - \exp(-0.693 t_s / t_{1/2}))}{0.693 t_s} \right] \quad (1.4)$$

where, t_r is the time from removal from neutron flux to the beginning of the count, t_s is the counting time, $t_{1/2}$ is the half-life of the radioactive product.

The activated atoms will decay to lower energy levels, in order to reach more stable states. If during this decay process, the emitted excess energy is in the form of gamma-rays, both qualitative and quantitative analysis are possible.

The interaction principle of gamma-rays with matter can be used in NAA to identify the gamma emitting isotope and its quantity in the decaying material. In this work main reference is made to scintillating crystals as interacting matter with gamma-rays.

A scintillation detector converts the incoming gamma-rays into scintillations, and into electrical signals through a photomultiplier (PM) tube. These electrical signals are then amplified, measured and registered. The measurements are performed by the use of multi-channel analyzers (MCA). A MCA is an electronic equipment which stores the electrical signals in different locations (channels) according to their amplitudes. Since the amplitude of the signal is proportional to the incident gamma-ray energy a pulse height versus energy distribution can be obtained. Figure 1.1 is a plot of the expected theoretical response of a monoenergetic gamma emitter of energy E_0 and the actual spectrometer response, i.e., gamma-ray spectrum (He-64).




Figure 1.1- Comparison of a theoretical and actual response of a scintillation spectrometer to a monoenergetic gamma emitter of energy E_0 (He-64).

1.4- OUTLINE OF THE THESIS

The NAA technique is defined, its use and a brief history is given in section 1. The theory of NAA is briefly stated in the subsection 1.2. Before going into a more detailed study of the general principles of NAA, a short summary of the principles behind NAA are stated and the gamma-ray spectrometry is introduced in the same section. Section 2. is a more detailed discussion of the general principles of NAA. In this section, the interaction of neutrons with matter, formation of radioactive nuclei, their decay schemes are studied. The following subsection defines the types of interaction of gamma radiation with matter. The formation of gamma spectrum is also introduced in more detail in this subsection. The detectors used in common NAA work, their principle of operation, the scintillation process and scintillation detectors are studied in subsection 2.3. Information on the electronic equipment used in NAA and the statistics of radiation can be found in subsections 2.4 and 2.5 respectively. The use of computers in NAA being an important event, some predeveloped computer codes are named in the last subsection of the second section. Among the computer codes which were available in NEA Computer Library (NEA-CPL), CORGAM is analyzed in detail. It is in section 3 that CORGAM is introduced and the theory used in this computer program is given in some detail in subsection 3.2. One can find the experimental procedures considered in this code in subsection 3.3. In this computer program, there exist some special techniques. These techniques are given in the appendices A, B, C, D and E. Finally a listing of the program, its flow diagram, the input data and the input options, the updating programs and a sample problem run are given in the appendices F, G, H, J, K, L and M.

2. GENERAL PRINCIPLES

2.1- FORMATION OF RADIOACTIVE NUCLEI THROUGH THE INTERACTION OF NEUTRONS WITH MATTER

If a moving particle collides with another one, kinetic energy is exchanged between them in agreement with the laws of conservation of energy and momentum. If the potential energy of the system remains unchanged, kinetic energy being conserved during the collision, the phenomenon is called elastic scattering. The scattering is inelastic if one of the particles is left in an excited state after collision. The symbols (n,n) and (n,n') are often used to denote these processes in neutron elastic and inelastic scattering. Neutrons may disappear in an interaction as the result of absorption reactions, the most important of which, is the (n,γ) reaction. This process is also known as radiative capture, since one of the products of the reaction is instantaneous gamma radiation. Neutrons also disappear in charged-particle reactions such as (n,p) or (n,α) reactions. Sometimes, two or more neutrons are emitted when a nucleus is struck by a high energy neutron. The process here are of the $(n,2n)$ or $(n,3n)$ type. A closely related process is the (n,pn) reaction, which also occurs with highly energetic incident neutrons. Finally, when a neutron collides with certain heavy nuclei, the nucleus splits into two large fragments with the release of considerable energy. This is the fission process. Detailed discussion of these interactions are given in many text books (Ev-55, La-72, Pi-64, Fr-67).

In neutron activation analysis the most important reaction is of the (n,γ) type. This type of reaction is the most common of all reactions, since thermal neutrons induce this reaction in nearly all nuclei. It also occurs with a very high probability for a number of nuclei at particular energies in the epithermal range. This latter phenomenon is known as resonance capture. The gamma-rays which are emitted in the (n,γ) reaction usually have energies of several MeV.

The emission of this gamma radiation is due to the release of excess energy of the excited nucleus. The energy of this gamma-ray is equal to the difference in the energy levels of the excited and stable states of the nucleus. In some cases two or more gamma-rays of discrete energies are emitted in cascade. These gamma-rays are usually emitted within 10^{-13} sec. after the beta-particle, electron capture or alpha-particle decay occurs. Sometimes, this excited state has a characteristic half-life of its own.

Whenever this half-life is of the order of 10^{-7} to 10^{-6} sec. or longer, the excited state is called an isomer of the product atom. And the subsequent decay of this isomer by the emission of a gamma-ray is called an isomeric transition (Ga-67). The usual case is illustrated by ^{60}Co (Figure: 2.1).




Figure: 2.1- Decay scheme of ^{60}Co (Le-68).

The detection of these gamma radiations by suitable detectors, and the analysis of the spectrum formed by these radiations, may be used to identify the radiating nuclide, its amount, etc. The spectrum formed by these gamma-rays, due to the reasons discussed in the following subsections is unique. That is a given radionuclide has its own spectrum. Once the gamma spectrum is obtained from an unknown nuclide, qualitative and quantitative analysis is possible by using special techniques. These techniques are studied in detail in many text books (So-69, He-64, Cu-68). A brief description of these methods can be found in the following subsections, and a computer application of a specially developed method, namely the correlation algorithm is presented in section 3 (Ec-68).

2.2- INTERACTIONS OF GAMMA-RAYS WITH MATTER AND GAMMA-RAY SPECTRUM

Gamma-rays, like other electromagnetic radiation, have the properties of both waves and of discrete particles, or quanta. At the relatively high energies of isotopic gamma radiation, the quantum characteristics are the most important. Thus, for the purposes of considering their interaction with matter, gamma photons may be regarded as particles having a given energy. There are four kinds of interaction of photons with matter. Fano (Fa-53) in his article states that these are:

- 1- Basic interaction with atomic electrons
- 2- Basic interaction with nuclear particles (nucleons)
- 3- Interaction with electric field surrounding charged particles (nuclei and electrons)
- 4- Interaction with meson field surrounding nucleons.

And there are three effects that the interaction may have upon a photon,

- A- Outright absorption
- B- Elastic (coherent) scattering
- C- Inelastic (incoherent) scattering.

Different types of processes may arise from each kind of interaction (1,2,3,4) leading to each of the end results (A,B,C). Thus there are twelve types of processes in all. Many of these processes are quite infrequent, some have not yet been observed. The following three processes are the most important ones. These are;

- a) Photoelectric effects, which is the basic interaction with atomic electrons leading to an outright absorption,
- b) Compton scattering, which is the basic interaction with atomic electrons leading to inelastic (incoherent) scattering,
- c) Pair production, which is the interaction with electric field surrounding charged particles (nuclei and electrons) leading to an outright absorption.

PHOTOELECTRIC EFFECT is most important for photons of low energy. In this process, an incident gamma photon impinges on a bound orbital electron and transfers all of its energy to the electron, ejecting it from the atom. If E_b is the binding energy of the electron and the recoil energy of the atom is neglected, then $E = hv - E_b$, where hv is the energy of the electromagnetic quantum and E is the kinetic energy imparted to the electron. The quantum of radiation is possible only because the remainder of the atom can receive some momentum. The photoelectric effect is frequently used to

determine gamma-ray energies. This may be accomplished by measurement of the total ionization due to the photoelectrons in a proportional or scintillator counter.

The photoelectric effect has its highest probability for electrons that are more tightly bound. The incident photon energy must exceed the binding energy of the electron, of course, before this reaction can occur. For incident photon energies above the K-shell binding energy, experiments show that 80 % of the photoelectric effect interactions occur with K-shell electrons. The probability for this interactions, increases rapidly with decreasing photon energies. Figure 2.3 shows how the photoelectric effect probability varies with incident photon energy.




Figure: 2.3- Photoelectric cross-section as a function of photon energy for lead (Ga-67).

The discontinuities occurring at the K and L shell binding energies are called the K and L absorption edges. The probability of photoelectric interaction is approximately proportional to $Z^{4.4}$ of the absorber and to $E^{-3.5}$ of the interacting photon (Ga-67). Thus the photoelectric effect occurs with high Z absorbers and low photon energies.

Instead of giving up its entire energy to a bound electron, a photon may transfer only a part of its energy to an electron, which in this case may be either bound or free; the photon is not only degraded in energy but also deflected from its original path. This process is called the Compton effect or COMPTON SCATTERING named after A.H. Compton who first reported it in 1922 (Ko-61). This effect is shown schematically in figure 2.4.




Figure: 2.4- Schematic diagram of Compton scattering of a gamma-ray by an electron (Fr-67).

Consider the expression which relates the total energy of a particle E to its momentum p .

$$E = (E_0^2 + c^2 p^2)^{1/2} \quad (2.1)$$

The quantity E_0 is the total energy of the particle when it is at rest and is given by $m_0 c^2$ where m_0 is the rest mass of the particle. It must also be recalled that the rest mass of a photon is zero.

If E_γ is the energy of the initial gamma-ray, E'_γ is the energy of the gamma-ray after being scattered by an angle θ , E_0 is the rest energy of the electron (0.511 MeV), p and φ are the momentum of the electron and the scattering angle of the electron respectively. Then the conservation of total energy requires that;

$$E_\gamma + E_0 = E'_\gamma + (E_0^2 + c^2 p^2)^{1/2} \quad (2.2)$$

The conservation of momentum along the X direction and Y direction are

$$\frac{E_\gamma}{c} = \frac{E'_\gamma}{c} \cos \theta + p \cos \varphi \quad \text{X-direction} \quad (2.3)$$

$$\frac{E'_\gamma}{c} \sin \theta = p \sin \varphi \quad \text{Y-direction} \quad (2.4)$$

Combining equations 2.3 and 2.4 the following equations are obtained.

$$E_{\gamma}^2 - 2E_{\gamma}E'_{\gamma} \cos \theta + E'_{\gamma}{}^2 = c^2 p^2 \quad (2.5)$$

$$\frac{1}{E'_{\gamma}} - \frac{1}{E_{\gamma}} = \frac{1 - \cos \theta}{E_0} \quad (2.6)$$

Since $E = \frac{hc}{\lambda}$ where c is the velocity of light, equation 2.6 becomes;

$$\lambda' - \lambda = \frac{h}{m_e c} (1 - \cos \theta) \quad (2.7)$$

where m_e is the rest mass of the electron. The quantity $\frac{h}{m_e c} = 2.426 \times 10^{-10}$ cm is the well known Compton wavelength of the electron. Equation 2.7 shows that for a given incident energy there is a minimum energy (maximum wavelength) for the scattered gamma-ray and that this occurs for scattering in the backward direction ($\cos \theta = -1$).

$$(E'_{\gamma})_{\min} = \frac{E_0}{2} \frac{1}{1 + E_0/2E_{\gamma}} \quad (2.8)$$

For large incident gamma-ray energies ($E_{\gamma} \gg \frac{1}{2} E_0$) the minimum energy of the scattered gamma-rays approaches $\frac{1}{2} E_0 = 250$ KeV. For this reason scintillation spectra of high-energy gamma-rays always show a "back-scattering" peak at ≤ 250 KeV which is caused by Compton scattering and a valley between photopeak and Compton continuum whose width corresponds to the minimum energy (≤ 250 KeV) carried off by gamma-rays Compton scattered in the crystal (Fr-67).

The Compton effect is important for photons in the energy range of 0.6 to 2.5 MeV. The probability of Compton interaction increases with increasing Z of the absorber and with decreasing energy of the gamma-ray.

The third mechanism by which electromagnetic radiation can be absorbed is the PAIR PRODUCTION process. The process is so termed because it results in the production of a pair of particles of opposite sign, the electron and the positron. As has been pointed out above, the positron is a particle having the mass and dimensions of an electron, but a positive rather than a negative charge. Pair production is possible only for a photon having an energy of at least 1.02 MeV. When such an energetic photon traverses a powerful electric field, such as that surrounding a nucleus, the photon can be annihilated with the production of a pair of particles of equal and opposite charge. In this way, the electrical neutrality of the original photon is conserved.

The energy equivalent of the mass of each particle is 0.51 MeV and the sum of these masses, 1.02 MeV, is the minimum photon energy required for production of particles. If the photon has an energy in excess of 1.02 MeV, that excess energy is dissipated as the kinetic energy of the electron and positron created. As in the case of photoelectric effect, the incident photon disappears in this process. The positron so formed is then annihilated by collision with orbital electrons, both positron and electron being destroyed in the collision. The energy resulting from the annihilation process is emitted as electromagnetic or gamma radiation. It can be seen that positron annihilation is the reverse of pair production, in that the positron is destroyed with the consequent production of two 0.51 MeV gamma photons. The absorption of a gamma photon by matter through pair production is always associated with the emission of this relatively low energy gamma-radiation.

From the figure 2.5 it can be seen that the photoelectric process is the most important process below about 1 MeV. The Compton process is fairly important over the range from 0.1 to 10 MeV, and pair production can take place only above 1.02 MeV. Because a gamma-ray can interact with a detector by these three principal processes, a fairly complex response, Figure 2.6 is obtained from the spectrometer system. The photoelectric process which results in the total deposition of the gamma-ray's energy in the scintillation crystal is characterized by the Gaussian-shaped full energy-peak. Compton scattering is characterized by the nearly flat continuum below the four small peaks. The backscatter peak corresponds to the photoelectric process for a source gamma-ray whose initial direction of travel was 180° away from the scintillator but was backscattered into the crystal. The energy of the source gamma-ray in excess of the rest mass energy of the electron-positron pair (1.02 MeV) is transferred to the electron-positron pair as kinetic energy. Therefore, a small peak, called the two photon escape peak (or double escape), will occur at a point equal to the source gamma-ray energy minus 1.02 MeV. If one of the 0.51 MeV annihilation gamma-rays interact in the crystal by the photoelectric process and the other escapes interaction, there will be a small peak, called the one photon escape peak (or single escape), located at a point equal to the source gamma-ray energy minus 0.51 MeV. The third peak, called the annihilation peak, is a result of photoelectric interaction of 0.51 MeV gamma-rays produced by the pair production process occurring in materials outside the scintillator detector. Because of these interactions, a gamma-ray interacting in a particular scintillation crystal produces a unique gamma-ray spectrum.




Figure: 2.5- Absorption coefficient for NaI(Tl) as a function of gamma-ray energy (Ec-68).

The detection of gamma-rays is mostly performed by means of a NaI(Tl) scintillator or with a Ge(Li) junction. As the output of these detectors is proportional to the energy of the incident gamma-ray energy, spectrometry is possible. According to the size and the absorbing properties of the detector and to the energy of the measured gamma-ray, total or partial absorption can occur. Figure 2.6 shows the gamma spectra ²⁴Na (E_1 2.76 MeV, E_2 1.38 MeV) obtained with a 2"X2" NaI(Tl) detector. In this figure (a), (b), (c) represents the total absorption while (d), (e), (f) denote partial absorption. All processes giving rise to total absorption contribute to the formation of the peak at the highest energy side of the spectrum: the so called photo-peak. The coincidence of two or more gamma-rays in the detector, gives rise




Figure: 2.6- Gamma spectrum of ^{24}Na measured with a NaI(Tl) (1) and with a Ge(Li) (2) detector (So-69).

to a sumpeak. The relative intensity of the sumpeak to the single peak is a function of the source detector geometry.

The energy of the unknown gamma-ray can be determined from the pulse-height of the photopeak maximum. The calibration of the spectrometer can be performed by means of isotopes emitting gamma-rays of known energies. catalogues of gamma-ray spectra measured with NaI(Tl) and Ge(Li) detectors are also available for the identification of gamma spectra (Cr-64 , He-64 , Le-68). In table 2.1a and 2.1b some calibration sources and respective energies are given.

Quantitative data of the gamma-ray intensity can be obtained from the area or the height of the photopeak. Since radiation is a statistical phenomenon, the photopeak can be described as a Gaussian error curve. The qualitative and quantitative of complex spectra, can be simplified by computer techniques. The qualitative analysis of the sample by comparison with the standards (reference spectra) can be performed by means of a least squares method, which also allows the calculation of the statistical error. The use of computers can be appreciated in activation analysis not only when dealing with a large number of routine analysis, but also for various tasks, e.g. gamma spectrum analysis, decay curve analysis, statistical treatment of data

optimizing, decay and counting conditions, etc. In section 3 a computer code "CORGAM" for complex gamma-ray analysis, is presented in detail. The method used in this code is given in the subsections of the same section, together with a complete listing of the program and a sample problem, in the appendices H, K and L.

| Isotope | Energy (keV) | Isotope | Energy (keV) |
|-------------------|----------------|--------------------------------|----------------|
| ²⁴¹ Am | 59.568 ± 0.017 | ⁹⁵ Nb | 765.83 ± 0.07 |
| ¹³¹ I | 80.166 ± 0.009 | ⁵⁴ Mn | 834.84 ± 0.07 |
| ¹⁵³ Gd | 97.43 ± 0.02 | ⁸⁸ Y | 898.01 ± 0.07 |
| ¹⁵³ Gd | 103.18 ± 0.02 | ²⁰⁷ Pb | 1063.82 ± 0.28 |
| ¹⁷⁷ Lu | 112.97 | ⁶⁰ Co | 1173.13 ± 0.04 |
| ¹⁴¹ Ce | 145.44 ± 0.05 | ⁶⁰ Co | 1332.39 ± 0.05 |
| ¹³⁷ Cs | 165.84 ± 0.03 | ²⁴ Na | 1368.40 ± 0.04 |
| | | ⁵⁶ Co(D.E.) | 1576.9 ± 0.32 |
| ¹⁷⁷ Lu | 208.36 | ²⁰⁸ Tl(ThC'')(D.E.) | 1592.3 ± 0.13 |
| ²⁰³ Hg | 279.12 ± 0.05 | ²⁴ Na(D.E.) | 1731.6 ± 0.16 |
| ¹³¹ I | 364.47 ± 0.005 | ⁸⁸ Y | 1836.1 ± 0.07 |
| | | ⁵⁶ Co | 2035.2 ± 0.50 |
| | | ⁵⁶ Co(D.E.) | 2232.2 ± 0.67 |
| ¹⁹⁸ Au | 411.776 ± 0.01 | ²⁴ Na(S.E.) | 2242.6 ± 0.14 |
| | | ⁵⁶ Co | 2598.9 ± 0.30 |
| Annihilation | 511.006 ± 0.02 | ²⁰⁸ Tl(ThC'') | 2614.3 ± 0.09 |
| ²⁰⁷ Pb | 569.65 ± 0.10 | ²⁴ Na | 2753.6 ± 0.12 |
| ¹³⁷ Cs | 661.59 ± 0.07 | ⁵⁶ Co | 3202.4 ± 0.65 |
| | | ⁵⁶ Co | 3254.2 ± 0.65 |
| | | ⁵⁶ Co | 3273.6 ± 0.4 |
| | | ⁵⁶ Co | 3452.5 ± 0.75 |

S.E. = single escape D.E. = double escape.

Table: 2.1a- Calibration sources for gamma-ray spectrometry (So-69).

| Energies (keV) | Relative energies | | | |
|-------------------|-------------------|--------------|--------|-------------|
| | This work | Ref. 5 | Ref. 4 | Ref. 2 |
| 80.57 ± 0.20 | 14.48 ± 0.48 | 14.55 ± 0.45 | 13.33 | 14.50 ± 2.9 |
| 94.55 ± 0.15 | 0.3 | | | 0.16 ± 0.03 |
| 121.60 ± 0.15 | 0.78 ± 0.18 | | | 0.7 ± 0.50 |
| 160.70 ± 0.20 | 0.36 ± 0.15 | | | 0.35 ± 0.1 |
| 184.50 ± 0.25 | 100 | 100 ± 1 | 100 | 100 |
| 215.60 ± 0.25 | 3.94 ± 0.09 | 4.15 ± 0.06 | 4.44 | 3.8 ± 0.4 |
| 231.20 ± 0.30 | 0.36 ± 0.03 | 0.32 ± 0.05 | | 0.3 ± 0.2 |
| 260.15 ± 0.30 | 1.77 ± 0.12 | 1.42 ± 0.10 | 1.11 | 1.8 ± 0.5 |
| 280.30 ± 0.30 | 38.61 ± 0.46 | 43.6 ± 0.4 | 33.33 | 39.5 ± 2.8 |
| 300.45 ± 0.40 | 4.77 ± 0.09 | 5.45 ± 0.05 | 4.44 | 4.8 ± 0.4 |
| 365.80 ± 0.30 | 2.93 ± 0.06 | 3.72 ± 0.07 | 2.22 | 2.9 ± 3 |
| 410.95 ± 0.30 | 15.50 ± 0.19 | 16.8 ± 0.2 | 13.33 | 15.8 ± 1.2 |
| 451.55 ± 0.25 | 3.48 ± 0.07 | 4.30 ± 0.08 | 3.32 | 3.5 ± 0.7 |
| 464.90 ± 0.25 | 2.00 ± 0.07 | 1.66 ± 0.08 | 2.22 | 2 ± 0.4 |
| 529.95 ± 0.30 | 10.16 ± 0.32 | 13 ± 0.40 | 13.33 | 10.3 ± 1 |
| 570.50 ± 0.30 | 6.77 ± 0.14 | 7.08 ± 0.14 | 7.77 | 6.8 ± 0.7 |
| 594.90 ± 0.20 | 1.28 ± 0.18 | 0.74 ± 0.10 | | 1.2 ± 0.4 |
| 610.95 ± 0.15 | 1.48 ± 0.27 | 1.59 ± 0.32 | | 1.4 ± 1.0 |
| 670.48 ± 0.2 | 7.01 ± 0.25 | 7.35 ± 0.29 | 6.66 | 7 ± 0.7 |
| 691.15 ± 0.15 | 1.85 ± 0.09 | 1.62 ± 0.08 | 1.11 | 1.9 ± 0.4 |
| 711.70 ± 0.20 | 71.65 ± 0.68 | 71.5 ± 0.7 | 64.44 | 72.5 ± 6 |
| 725.50 ± 0.2 | 0.2 | | | 0.19 ± 0.1 |
| 736.70 ± 0.20 | 0.46 ± 0.04 | 0.50 ± 0.05 | | 0.45 ± 0.15 |
| 752.50 ± 0.25 | 16.06 ± 0.40 | 15.2 ± 0.30 | 15.55 | 16.1 ± 1.2 |
| 779.20 ± 0.45 | 3.72 ± 0.07 | 3.88 ± 0.06 | 4.44 | 3.8 ± 0.3 |
| 793.10 ± 0.20 | 0.2 | | | 0.17 ± 0.1 |
| 810.80 ± 0.40 | 76.38 ± 0.82 | 76.4 ± 0.80 | 66.66 | 76 ± 8 |
| 830.80 ± 0.30 | 12.07 ± 0.25 | 12.9 ± 0.30 | 12.22 | 12.5 ± 1 |
| 875.90 ± 0.20 | 1.14 ± 0.07 | 0.91 ± 0.04 | 0.88 | 1.15 ± 0.15 |
| 889.20 ± 0.15 | 0.7 ± 0.15 | | | 0.73 |
| 896.50 ± 0.10 | 0.7 ± 0.15 | | | 0.73 |
| 951.50 ± 0.20 | 3.5 ± 0.14 | 3.16 ± 0.12 | 3.66 | 3.6 ± 0.6 |
| 995.35 ± 0.25 | 0.68 ± 0.18 | | | 0.73 ± 0.24 |
| 1118.10 ± 0.20 | 0.30 | 0.26 ± 0.02 | 0.33 | |
| 1147.20 ± 0.20 | 0.38 ± 0.05 | 0.26 ± 0.02 | 0.22 | 0.38 ± 0.06 |
| 1241.80 ± 0.30 | 1.22 ± 0.05 | 1.06 ± 0.04 | 1.11 | 1.25 ± 0.25 |
| 1275.00 ± 0.25 | 0.75 ± 0.20 | | | 0.8 ± 0.15 |
| 1283.10 ± 0.20 | 0.38 ± 0.04 | 0.22 ± 0.02 | 0.33 | 0.8 ± 0.15 |
| 1401.30 ± 0.15 | 0.86 ± 0.05 | 0.72 ± 0.02 | 0.77 | 0.93 ± 0.09 |
| 1427.50 ± 0.20 | 0.65 ± 0.03 | 0.69 ± 0.02 | 0.66 | 0.69 ± 0.07 |

^a The intensity 184.50 keV is taken as 100 units.

Table: 2.1b- Energies and relative intensities^a in the decay of 1200 y ¹⁶⁶Ho^m. The errors quoted are one standard deviation.

2.3- SCINTILLATION DETECTORS

One of the oldest methods of nuclear-radiation detection makes use of light scintillations. It is well known that Rutherford and his collaborators used this technique for detecting alpha particles in the famous alpha-scattering experiments (Pr-64). The principle of these detectors is the classical one, long used in the spinthariscopes also (P1-62). The same principle is applied in the detectors of radioactivity which are called scintillation counters; these consist of two parts, a scintillator in which the luminous effects are produced, and a photomultiplier (PM) to convert the light pulses into high amplitude electrical signals. The chief advantages of scintillation counters are that they provide pulses proportional to the energy of the particles, and they have high efficiency for gamma radiation as well having large sensitive areas. Before discussing the scintillation detection mechanism, it is necessary to discuss the scintillation phenomenon itself. Many materials are excited by radiation and subsequently lose part of this excitation energy by the emission of low energy photons of radiation. However very few materials emit radiation of a wavelength that easily interacts with photomultiplier phosphor material and yet easily penetrate common glasses. Even fewer are transparent to their own emitted radiation (Ga-67). Luminescence is the general term describing the absorption of energy of any type by a substance and the reemission of some fraction of it as photons of visible or near visible light. Radiation energy is but one type of energy that can excite materials to initiate luminescence. Two types of luminescent emission can occur. By absorption of energy, electrons can be raised to higher excited states and can return to the ground state either directly or via a metastable state. This deexcitation which occurs in a time shorter than 10^{-8} s is called fluorescence (immediate emission). The energy of the emitted light is less than that of the absorbed one to reach the excited state. The electron occupying a metastable state, can be raised to a higher excited state and subsequently return to the ground state by emission of a photon. This process called phosphorescence, occurs after a time interval of about 10^{-8} sec after excitation (delayed emission). Scintillation theory is discussed in more detail by many authors (Bi-53, Cu-53).

The following groups of materials are known to possess desired properties as scintillators:

1. Crystalline inorganic substances
2. Crystalline organic substances

3. Solid organic solutions (plastics)
4. Liquid organic solutions
5. Rare gases

Among the desirable properties for good scintillators are high yield of fluorescent light, transparency for this emitted light, rapid light-emission in combination with the lowest possible phosphorescence and finally, a spectral distribution of the light suitable with regard to the sensitivity of existing photomultipliers. The demands in other respects, (density, size, geometrical shape, etc.) depend on the special situation in which the detector is to be used. A survey of the characteristics of the most commonly used phosphors is given in table 2.2.

| Materials | Group | Density gm/cm ³ | Wavelength of maximum emission A° | Relative pulse height % | Alpha beta ratio % | ev/photon | Decay time nsec |
|---------------------------|----------------------------|-------------------------------|--------------------------------------------|----------------------------------|-----------------------------|-----------|-----------------------|
| Anthracene crystal | Single cryst. organic | 1.25 | 4470 | 100 | 9 | 65 | 30 |
| Trans-stilbene crystal | Single cryst. organic | 1.15 | 4110 | 60 | 9 | 100 | 5 |
| Liquid phosphors | Organic solution | 0.86 | ~4000 | 40-60 | 9 | 150 | 5 |
| Plastic phosphors | Organic solution | 1.06 | ~4000 | 28-48 | 9 | 150 | 4 |
| NaI(Tl) | Single cryst. inorganic | 3.67 | 4130 | 210 | 44 | 30 | 250 |
| ZnS(Ag) | Single cryst. inorganic | 4.10 | 4500 | 200-300 | 100 | | 200 |
| LiI(Eu) | Single cryst. inorganic | 4.06 | 4700 | 70 | 95 | | 1200 |

Table: 2.2- Properties of some commonly used phosphors (Pr-64, So-69).

The inorganic scintillators are crystals of inorganic salts, primarily the alkali-halides, containing small amounts of impurities as activators for the luminescent process. The yield of fluorescent light is connected with the existence of luminescent centers in the crystal lattice. For these, the transition from the excited state, to the ground level is accompanied by the emission of a photon. Luminescent centers are created by the addition of certain activating substances. Sodium Iodide, activated by the addition of thallium, is the most widely used scintillator for gamma-ray spectrometry.

The important properties, and the theory of various other scintillators are given in many text books (Pr-64, Bi-53, Cu-53, So-69).

Radiation can be measured by the fluorescence phenomenon if a means of measuring the amount of light produced is available. The photomultiplier tube serves this purpose. A PM tube is a phototube in which the currents of photo-electrons from the tubes cathode is multiplied many times through the process of electron multiplication. The principle elements in each of these tubes are the photocathode for releasing the photoelectrons, the dynodes for producing the electron multiplication, and the anode for collecting the current pulse. An incident photon interacts with the photocathode material (usually a mixture of Antimony and Cesium) deposited on the inside of the glass cylinder and produces a photo-electron. A positive voltage applied to the first dynode with respect to the photocathode draws this photo-electron to the dynode. The photo-electron is accelerated sufficiently to cause electron multiplication upon its interception by the first dynode. Several electrons are typically emitted in the interaction of the photo-electron with the first dynode. These, in turn, are drawn to the second dynode which is maintained at a positive potential higher than the first one. The electron multiplication continues in this fashion at each dynode until all the electrons so produced are collected at the anode, which is located at the end of the dynode chain. The amount of electron multiplication at each dynode is a strong function of the voltage between dynodes. An approximate relation for the multiplication factor, M , at each stage is,

$$M = KV^7 \quad (2.9)$$

where V is the voltage in volts and K is a constant, and the total multiplication is given as :

$$A = M^n \quad (2.10)$$

where M is the multiplication at each stage and n is the number of stages. A more detailed discussion on PM tubes is given by Price (Pr-64). Figures 2.7, 2.8 and 2.9 are the schematic diagrams of three type of dynode arrangement for PM tubes.

Two important properties of PM tubes are the relative constancy of the electron multiplication factor M and of the transit time of the electrons through the multiplier structure. The former property effects the pulse-height resolution of the scintillation detector and the latter the resolving time. The resolution of a detector system is a measure of the extent

to which monoenergetic particles produce pulse heights of a single value. In gamma-ray spectrometry, the spread in the total-energy peak is a measure of the resolution, which is normally defined as,

$$R = \frac{W_n}{E_n} \quad (100) \quad (2.11)$$

where R is the resolution in %, W_n is the full width at half maximum (MeV), E_n is the gamma-photon energy (MeV).

^{137}Cs generally is used as a comparison standard for resolution of a crystal. Literature values for typical resolutions range from 5-13 % (Pr-64).




Figure: 2.7- ASCOP photomultiplier tube with Venetian-blind dynode structure. (Pr-64)




Figure: 2.8- Du Mont 6292 photomultiplier tube; diagram and photograph Box type dynode structure. (Pr-64)




Figure: 2.9- RCA 6342 photomultiplier tube; diagram and photograph circularly focussed dynode structure. (Pr-64)

Figure 2.10 is the pulse-height spectrum obtained with a NaI(Tl) counter for a ^{137}Cs source.




Figure: 2.10- Pulse-height spectrum obtained with a NaI(Tl) scintillation counter for a ^{137}Cs source. (P1-62)

The resolving time for counting measurements is the minimum time which must exist between successive events if they are to be counted as two events. Table 2.3 is a compilation of the important properties of various PM tubes.

| Type | Maximum outside dimensions, in. | | Cathode | | | Dynode system | | | | Over-all voltage, volts |
|----------|---------------------------------|-------------------|-----------------------|---------------------|--------------------------------------|----------------|----------------|--------------------|--------------------|-------------------------|
| | Diam. | Length | Area, cm ² | Response | Av. sensitivity, $\mu\text{a/lumen}$ | No. of dynodes | Type | Material | Gain | |
| | | | | | | | | | | |
| RCA: | | | | | | | | | | |
| 5819 | 2 $\frac{1}{4}$ | 5 $\frac{13}{16}$ | 14.2 | S-11 | 50 | 10 | Circ. focused | Cs ₃ Sb | 2.3×10^6 | 1,250 |
| 6342 | 2 $\frac{1}{4}$ | 5 $\frac{13}{16}$ | 14.2 | S-11 | 60 | 10 | Circ. focused | AgMg | 0.55×10^6 | 1,500 |
| 6655 | 2 $\frac{1}{4}$ | 5 $\frac{13}{16}$ | 14.2 | S-11 | 50 | 10 | Circ. focused | Cs ₃ Sb | 2.3×10^6 | 1,250 |
| 6199 | 1 $\frac{9}{16}$ | 4 $\frac{9}{16}$ | 7.75 | S-11 | 45 | 10 | Circ. focused | Cs ₃ Sb | 2.8×10^6 | 1,250 |
| 6372 | 2 $\frac{9}{16}$ | 7 $\frac{3}{4}$ | 80 | S-11 | 33 | 10 | Circ. focused | Cs ₃ Sb | 2.5×10^6 | 1,200 |
| 1P28 | 1 $\frac{9}{16}$ | 3 $\frac{11}{16}$ | 1.9 | S-5 | 40 | 9 | Circ. focused | Cs ₃ Sb | 5.3×10^6 | 1,250 |
| 6810A | 2 $\frac{1}{4}$ | 5 $\frac{13}{16}$ | 14.2 | S-11 | ... | 14 | Circ. focused | AgMg | 12.5×10^6 | 2,000 |
| 8055 | 5 $\frac{1}{2}$ | 6 $\frac{3}{4}$ | 97 | S-11 | 75 | 10 | Venetian blind | CuBe | 1.5×10^6 | 2,000 |
| Du Mont: | | | | | | | | | | |
| 6291 | 1 $\frac{1}{2}$ | 4 $\frac{1}{4}$ | 6.4 | S-11 | 60 | 10 | Box | AgMg | 2×10^6 | 2,100 |
| 6292 | 2 $\frac{1}{16}$ | 5 $\frac{5}{8}$ | 13.4 | S-11 | 60 | 10 | Box | AgMg | 2×10^6 | 2,100 |
| 6363 | ... | 6 $\frac{1}{8}$ | 31.4 | S-11 | 60 | 10 | Box | AgMg | 2×10^6 | 1,800 |
| 6364 | 5 $\frac{1}{4}$ | 7 $\frac{1}{2}$ | 88.8 | S-11 | 60 | 10 | Box | AgMg | 2×10^6 | 1,800 |
| 6365 | $\frac{3}{4}$ | 2 $\frac{3}{4}$ | 1.26 | S-11 | 50 | 6 | Box | AgMg | 3×10^6 | 1,200 |
| EMI-US: | | | | | | | | | | |
| 9536B | 2 | 4 $\frac{3}{4}$ | 14.6 | S-11 | 70 | 10 | Venetian blind | Cs ₃ Sb | 3×10^6 | 1,700 |
| 9545B | 12 $\frac{1}{2}$ | 14 $\frac{1}{8}$ | 470 | S-11 | 50 | 11 | Venetian blind | Cs ₃ Sb | 4×10^6 | 2,000 |
| CRS: | | | | | | | | | | |
| CL1012 | 1 $\frac{1}{2}$ | 4 $\frac{3}{4}$ | 8 | S-11 | 60 | 10 | Venetian blind | AgMg | 2×10^6 | 1,750 |
| ASCOP: | | | | | | | | | | |
| 511A | 1 $\frac{1}{4}$ | 4 $\frac{1}{4}$ | 5.1 | (Depends on window) | 70 | 11 | Venetian blind | AgMg | 10^6 | 2,100 |
| 513A | 2 | 7 $\frac{1}{2}$ | 14.6 | S-11 | 60 | 11 | Venetian blind | AgMg | 10^6 | 2,100 |

Table: 2.3- Characteristics of a few commercial photomultipliers.(

Besides scintillation counters a wide variety of detectors are used for radiation detection purposes. A detailed study of these detectors may be found in references (Pr-64, Cr-64, Bi-53, Cu-53, Ko-61, Ga-67, Ev-55).

Gas filled detectors:

Three of the oldest but still very useful nuclear radiation detector types are the ionization chambers, the proportional counter, and the Geiger-Müller (GM) tubes. The difference in the three systems can be seen through the pulse-height versus applied-voltage curves of figure 2.11. The system consists of a gas-filled chamber with a central electrode well insulated from the chamber walls. These detectors use the principle of electrical conductivity of the gas resulting from the ionization produced by the nuclear radiation.

Semiconductor radiation detectors:

In this detector, an electric field is set up across a semiconducting medium of low electrical conductivity. When a charged particle passes through the medium, electron hole pairs are produced in it. The charges are separated by the field and collected at electrodes, producing an electrical signal which can be amplified, recorded and used to give information about the particles detected.

Several other types of detectors are in use and are quite important for specialized detection problems. These detectors are Cloud Chambers, Nuclear Track Plates, Spark Chambers, Cherenkov Counters, Chemical Detectors, Calorimetric Methods, and various types of Neutron Detectors. (Pr-64).


2.4- AUXILIARY ELECTRONIC EQUIPMENT

As it is discussed in the previous subsections the output of a scintillation detector is an electrical signal. In order to analyse this signal some auxiliary equipment is needed. A complete instrument for the measurement of nuclear radiation consists of a combination of a detector one or more electronic units which perform the functions required. The main functions of these electronic devices are supplying the necessary voltage to the detector amplifying the detector output signal to a suitable level; registering, presenting and analysing the amplified signal (Er-64).

Like all others, scintillation detectors also, need a very well stabilized high voltage supply, which is used to create a potential difference between the successive dynodes of the PM tube. Preamplifiers are used where signals are so small that a normal amplifier is not able to amplify them sufficiently or where some amplification or impedance matching is required to offset cable losses. Finally there exists a large variety of electronic equipment which is used as the analysing apparatus (Pr-64).

One of the important properties of a scintillation detector is its ability to give an indication of the energy of the incident radiation. Because of this property, the pulse-height distribution curves from a scintillation detector can be used to determine the energy spectrum of incident radiation. The object in pulse-height analysis is to determine the total number of pulses whose heights fall in selected intervals or channels. If this analysis is made by sampling one channel at a time, the device is referred to as a single-channel pulse-height analyzer (SCA). However, in many applications it is essential to accumulate data without any loss, so multi-channel analyzers (MCA) are used for this purpose. These instruments all consist of some or all of the following elements: 1) a section that associates each input signal with a specific amplitude channel, 2) a memory, or data storage device that contains the information on this number of pulses in each channel until this information is recorded otherwise, 3) a data display device, 4) auxiliary readout devices such as curve plotters, paper-tape or card-punches, magnetic tape, typewriters, etc.

Below is a block diagram of an electronic arrangement used in pulse-height spectrometry. An over all discussion of the techniques of MCA's as well as other methods for pulse-height-to-digital conversion is given by Chase (Ch-61).



The most widely used method employs pulse-height-to-time converters. Figure 2.12 illustrates this process. First the signal pulse "a" is converted to a flat-topped one "b" by means of a pulse-stretching circuit. At the start of the pulse, a linearly rising signal "c" is switched on, and a series of equally spaced pulses "d" is started. When "c" reaches "b", the pulse train "d" is stopped. In this manner the pulse heights are converted into trains of clock pulses, and different channels now correspond to trains with different numbers of clock pulses. These trains of pulses can be sorted, stored, displayed and read out by various techniques that have been developed for digital computers (Ch-61).




Figure: 2.12- Schematic representation of pulse-height-to-time conversion.

The MCA's are a logical extension of SCA's. By employing n differential discriminators set at successive window levels and each having its own counting system, an n -channel analysis is obtained. The principle advantage of this approach is that its time resolution can be made less than 1 μ sec. However there are two serious disadvantages. First, this arrangement, having n separate channels consequently n circuits, is costly to construct and maintain. Second, it is a difficult design to maintain equal window widths. In general the MCA's are very sophisticated instruments, with good accuracy and versatility.

2.5- STATISTICAL INTERPRETATION OF THE RESULTS

The measurement of nuclear radiation by the detection systems introduced in previous subsections, involve phenomena which are statistical in nature. By the use of statistical analysis, estimates can be made of the accuracy of the measurements, and procedures can be devised which will minimize errors due to random nature of the process. Obviously, the statistical analysis cannot obviate the need for sound analytical techniques, but unnecessary loss of operator time in the laboratory can be eliminated and the reliability of the results increased by the use of statistical methods.

Two kinds of error can be met in radiation detection; 1) systematic errors, which determine the accuracy of a result, 2) random errors, which determine the reproducibility or precision. By systematic errors the experimental results deviate from the true value, because of a bias. Statistics, on the other hand, deal with random errors. It makes no sense to define an error as the deviation of a measurement from the true value, as the latter is usually unknown. Thus "error" can only be defined after evaluation of the most probable value of the quantity measured. For a number "n" of observations x_1, x_2, \dots, x_n , made under the same conditions and all equally reliable, the arithmetic mean (\bar{x}) is the best value for this set of observations.

$$\bar{x} = \sum_{i=1}^n x_i / n \quad (2.12)$$

If all observations are not equally reliable a statistical weight W_i can be introduced. The most probable value is then the "weighted mean", sometimes called "general or probable mean" (Gr-74):

$$\bar{x} = \frac{\sum_{i=1}^n W_i x_i}{\sum_{i=1}^n W_i} \quad (2.13)$$

The choice on W_i , the weighing parameters, depends upon the experimental conditions and the scatter of data. The reproducibility of an analytical method is usually expressed by the standard deviation $\sigma(x)$ on a single determination. The standard deviation is defined as the square root of the average value of the square of the individual deviations from the most probable value (mean), taken for a large number of observations. Thus,

$$\sigma^2 = \frac{1}{n-1} \sum_{i=1}^n (\bar{x} - x_i)^2 \quad (2.14)$$

where $n-1$ indicates the degrees of freedom of the observation set.

As the arithmetic mean for an infinite number of observations is usually not accessible, it is only possible to calculate the probability that this value is contained within the limits $\bar{x} \pm \Delta x$. The number Δx can be chosen arbitrarily small. It is common practice to express Δx in terms of the standard deviation, e.g. $\pm \Delta x = \pm t \sigma$, where t depends on the probability level and on the number of degrees of freedom. If the latter is large, t is replaced by u . Where "u" is the number of standard deviations for large degrees of freedom. The degree of confidence increases, as the "confidence interval" increases (see table 2.4).

| Constant u | Probability P | P' | Error |
|--------------|-----------------|---------|--------------------------------------------------|
| 0.6745 | 0.500 | 0.500 | "Probable error" |
| 1.000 | 0.683 | 0.317 | "Standard deviation" (σ) (rms error) |
| 1.177 | 0.761 | 0.239 | |
| 1.645 | 0.900 | 0.100 | "nine-tenths error" or "reliable error" |
| 1.960 | 0.950 | 0.050 | "ninety-five hundredths error" |
| 2.576 | 0.990 | 0.010 | "ninety-nine hundredths error" |
| 2.807 | 0.995 | 0.005 | |
| 3.000 | 0.9973 | 0.0027 | 3σ |
| 3.291 | 0.999 | 0.001 | |
| 4.000 | 0.99994 | 0.00006 | 4σ |

Table: 2.4- Probability levels (So-69).

The radioactive decay has a statistical character. If a radionuclide is counted several times in identical experimental conditions, different counting rates (number of counts per unit time) will be observed, even for very long lived species. It can be shown that the probability $P(n)$ of obtaining n disintegrations in a time Δt from N_0 original radioactive atoms is given by the Binomial distribution (Fr-67).

$$P(n) = \frac{N_0!}{(N_0 - n)! n!} (1 - e^{-\lambda \Delta t})^n (e^{-\lambda \Delta t})^{N_0 - n} \quad (2.15)$$

where $(1 - e^{-\lambda \Delta t})$ is the probability of an atom disintegrating in a time Δt and $e^{-\lambda \Delta t}$ the probability of surviving through the time Δt . The Binomial distribution may be replaced by the simpler Poisson distribution (Fr-67), the probability $P(n)$ of finding n disintegration is then given by equation 2.16, where m is the average to be expected

$$P(n) = \frac{m^n e^{-m}}{n!} \quad (2.16)$$

on condition that $\lambda \Delta t \leq 0.01$ and $N_0 \gg 100$.

For small disintegration rates, the distribution is not symmetrical around $n=m$. For m (or n) > 100 and $m \cong n$ the Poisson distribution can be replaced by the normal or Gaussian distribution:

$$P(n) = \frac{(2\pi)^{-1/2}}{\sigma} \exp\left(-\frac{(m-n)^2}{2\sigma^2}\right) \quad (2.17)$$

where, σ and m designate the standard deviation and the true mean, respectively. Figure 2.13 illustrates the form of the distribution along with its most important properties.




Figure: 2.13- Illustration the normal distribution function (Pr-64).

It will be noticed that this distribution is symmetrical about $n=m$. For the Poisson and Gaussian distributions one may derive $\sigma = \sqrt{m}$, or for large n , $\sigma \cong \sqrt{n}$.

An important concept of statistical interpretation of radioactive disintegration is the confidence intervals, discussed in detail by Curie (Cu-68). If the observed number of counts from a sample + background and from the background alone, for an equal time of measurement Δt , are respectively C_{s+b} and C_b , then the count difference $C_s = C_{s+b} - C_b$ is a measure of the net activity of the sample. However, owing to the statistical fluctuations, this count difference at the averages of the background \bar{C}_b and the sample \bar{C}_s may obtain various values. At a relatively small value of C_s , the presence of the activity in the sample becomes doubtful. Therefore a so-called "critical value" L_c is introduced, that is often called "decision limit" (Cu-68). For $C_s > L_c$ one assumes that the signal is present, for $C_s < L_c$ the signal will not be detected. Due to this assumption, two kinds of errors may occur:

1-The measured value $C_s > L_c$, one concludes that the activity > 0 when in fact it is equal to zero (type 1 error).

2-The measured value $C_s < L_c$, one concludes that the activity is equal to zero when in fact the activity > 0 (type 2 error).

As appears from figure 2.14, the critical level is mathematically given by:

$$L_c = u_1 \sqrt{0} \quad (2.18)$$

The parameter u_1 is a multiple of the standard deviation, determining the probability of making the type 1 error. Usual values for u_1 are 1.645, 1.960, 3.000. Contrary to table 2.4, this does not correspond to a confidence level of 90, 95 and 99.73 % respectively but to a confidence level of 95, 97.5 and 99.865 % respectively. Indeed the probability of making the type 1 error is respectively 5, 2.5, 0.135 % only (one side of the Gauss curve). For that reason, the symbol u has been replaced by u_1 .




Figure: 2.14- Signal detection. Definition of critical level L_c (So-69).

A number of tests have been devised to answer the question as to whether or not a given set of data fits a certain distribution function. The concept of Poisson's index of dispersion is particularly useful in this regard. To make this test, the quantity χ^2 , defined as

$$\chi^2 = \frac{\sum_{i=1}^n (\bar{x} - x_i)^2}{\bar{x}} \quad (2.19)$$

is calculated. Here n is the number of times that a given counting determination is made, x_i is the value obtained in the i^{th} measurement, and \bar{x} is the

average value. If the value of χ^2 is sufficiently large so that there is only a small probability of occurrence of a value this high, it is an indication of trouble in the experiment. Likewise if χ^2 is small, the corresponding data should be questioned because the spread in it is so small. Thus it is seen that the Poisson index-of-dispersion test can be helpful in determining whether or not the counting equipment is functioning satisfactorily.

Chauvenet's criterion is frequently used for the rejection of outlying data. The criterion states that any one in a series of n readings may be rejected when its deviation from the mean of the series is such that the probability of occurrence of all deviations from the mean as large or larger is less than $n/2$ (So-69, Gr-74).

An overall application of statistical methods briefly mentioned here can be found in the computer code CORGAM, which is introduced in section 3.

2.6- USE OF COMPUTERS IN NEUTRON ACTIVATION ANALYSIS

When experimental data resulting from an activity counting and thus subject to statistical fluctuations, have to be treated, a smoothing routine can be desirable. This is the case when dealing with gamma spectra or decay curves, which have to be analyzed by means of a computer. The simplest smoothing routine is the graphical smoothing, but it can be applied only when data are graphically evaluated.

Mathematical data convolution techniques largely correct for statistical scatter. The simplest convolution sets the convolution integer equal to one, giving rise to a "rectangular" convolution, i.e., a normal moving average. One also can use a symmetrical triangular convolution, one-sided exponential or symmetrical exponential convolution. But a more exact procedure can be achieved by using a least squares method. In this technique one tries to fit the best curve through the set of points under consideration. In this way it is possible to obtain a set of equations of the n^{th} order with $2n+1$ coefficients as unknowns. It is obvious that the calculations needed for one smoothed point are time consuming.

After data smoothing, the analysis of the gamma spectra can be performed in several ways. With the development of high resolution detectors, using a large number of channels per spectrum, the computer becomes practically unavoidable for analysis purposes.

The location of the photo-peak in a gamma spectrum, and the calculation of the corresponding area is of extreme importance in NAA. This can be done by a computer program, and qualitative and quantitative analysis can be made possible in this way. The technique employed in photo-peak positioning is the use of the first derivative of the smoothed spectrum. The exact location of the peak can be found by means of the second derivative. Both, first and second derivative can be calculated according to the least squares method.

Once the photo-peaks have been located, the activity under the peaks can be determined by subtracting the detector background. As the net photo-peak activities are obtained, one could conceive a complete activation analysis program. By treating the spectra of standards and samples in the same way, and correcting for irradiation time, decay time, and counting time, the results of the analysis can be computed. When the components or constituents of a complex gamma spectrum are qualitatively known, as is often the case in NAA, the only problem is to determine "how much" of each is present

in the sample. From a complex spectrum and from the spectra of composing standards one can compute the best value for the ratios of each component to its respective standard and also the standard deviation (So-69).

A great number of computer programs are available, for the analysis of complex gamma-ray spectra (NEA-CPL). Among them the ones recently prepared are given in the appendix N.

3. C O R G A M

3.1- CORGAM-A COMPUTER CODE FOR NEUTRON ACTIVATION ANALYSIS

The importance of computers in NAA was mentioned in previous sections, and a list of computer programs developed for the analysis of gamma-ray spectra is also given in appendix N.

Each of these codes have their own techniques in unfolding complex gamma-ray spectra. All these codes are prepared in order to estimate the intensity coefficients of the general linear hypothesis, which will be discussed in the subsection 3.2. There are many methods for estimating these coefficients. These methods, stripping, least squares, incremental, iterative, vector analysis, and stepwise regression, have been used for estimating these coefficients. All of these procedures may use either least squares or approximate least squares techniques. "Regression analysis" has been used to describe the least squares procedure for fitting a dependent variable, e.g., a complex spectrum, to a set of independent variables, e.g., a set of reference spectra. Because a model is used in which only significant variables are retained in the final solution, a more appropriate description for this work is "Correlation Analysis". (NEA-CPL, Ec-68).

Among the codes available by NEA-CPL in the experimental data processing field, "CORGAM" (A Correlation Algorithm for Gamma-Ray Spectra) is the one which uses correlation analysis in unfolding gamma-ray spectra.

In the present section of the thesis, a detailed description of CORGAM will be given. The listing of the program, flow diagram, input data for sample problem, instruction for users and the output of a sample problem are given in the appendices.

In CORGAM, a correlation algorithm is coded to allow the unfolding of complex gamma-ray spectra, typically collected in a NAA procedure.

CORGAM,

- a) Compensates for electronic shifts in the data,
- b) Corrects for the background,
- c) Normalizes the data to a fixed neutron flux level,
- d) Allows a choice of weighting factors,
- e) Allows a choice of methods for calculation of standard deviations.

The code requires a matrix of reference gamma-ray spectra. These spectra

can be in a raw data form. All of the modifications available to the complex gamma-ray spectra are also available to the reference gamma-ray spectra. In addition, a decay correction is available for the reference gamma-ray spectra. Only the reference gamma-ray spectra that have intensity coefficients which have significance at a prescribed level of significance are retained in the final solution. The intermediate ones, i.e., those solutions that contain reference gamma-ray spectra which have nonsignificant intensity coefficients at the prescribed level, are printed out. Therefore, several solutions are imbedded in the final solution.

A backward elimination least squares method is employed. The elimination is one in which the student's t value of each intensity is compared to a t value at a prescribed level of significance (Gr-74). The calculation of the standard deviation for each intensity coefficients accounts for randomness in both the complex and reference gamma-ray spectra.

Currently the program is restricted to a problem of 400 channels and 15 reference spectra. There would be no difficulty in expanding this to as many as 4096 channels. However, the number of reference spectra is limited by the accuracy desired in the inversion procedure. The inversion program currently being used would suffer (from the accuracy standpoint) if the number of reference spectra were greater than 25.

CORGAM is a novel method of unfolding complex gamma-ray spectra because the variance calculations include contributions from both the complex and reference spectra. Also, only reference spectra that have intensity coefficients which are significant at a preselected level are retained in the output. The intermediate results can be printed, and the level of significance of each of the intensity coefficients can be observed.

A typical analysis, involving 256 channels and 14 reference spectra and using reference spectra which were suitably aligned, requires less than 1 minute of execution time on IBM 360/50.

CORGAM was coded for IBM 360/50 using Fortran IV, level G. An intermediate disk storage is required. The source deck contains 1303 cards. The code is updated by the author for UNIVAC-1106 using Fortran IV and Cobol. The updated form of the code is recorded onto a magnetic tape for further use in BU-NE Department.

The tape contains in the,

1. File: IBM 360/50 version of CORGAM sample data,
2. File: UNIVAC-1106 version of CORGAM,
3. File: Sample problem input data,
4. File: Sample problem output data.

3.2- THEORY OF CORGAM

If a nuclear reaction, such that



is forced to proceed, the radioactivity of the product nuclei Y, is directly proportional to the number of atoms of the parent nuclei X. The equation governing the rate of formation of Y is given by the well known differential equation.

$$\frac{dN}{dt} = \dot{N} = \sigma N_0 \phi - \lambda N \quad (3.2)$$

where,

N is the time dependent nuclide (Y) concentration,

\dot{N} is the derivative with respect to time,

N_0 is the time dependent parent nuclide (X) concentration,

σ is the microscopic activation cross-section,

ϕ is the time dependent neutron flux.

The parent nuclide concentration can be given in terms of N and the initial parent nuclide concentration, N_0^0 , as

$$N_0 = N_0^0 - N \quad (3.3)$$

Therefore equation 3.2 can be written as,

$$\frac{dN}{dt} = \dot{N} = \sigma \phi N_0^0 - (\lambda + \sigma \phi) N \quad (3.4)$$

If the neutron flux is a constant $\phi = \phi_0$ then equation 3.4 has the following solution,

$$N = \frac{\sigma \phi_0 N_0^0}{(\lambda + \sigma \phi_0)} (1 - e^{-(\lambda + \sigma \phi_0)t}) \quad (3.5)$$

It can be seen from equation 3.5 that the time dependency of N_0 can be ignored if $\sigma \phi_0 \ll \lambda$ or $\sigma \phi_0 / \lambda \ll 1.0$.

The following table is a compilation of λ , σ and $\sigma \phi_0$ values for many of the elements used in common NAA work. The value for ϕ_0 used in this table

of the Nuclides" by D.T.Goldman 1965.

| isotope | $t_{1/2}$ sec $\times 10^{-3}$ | λ sec $^{-1}\times 10^4$ | $\sigma^a\phi_0$ sec $^{-1}\times 10^{11}$ | $\sigma\phi_0/\lambda$ |
|------------------|-----------------------------------|-------------------------------------|-----------------------------------------------|------------------------|
| Na ²⁴ | 54.00 | 0.128 | 0.13 | 1.01×10^{-7} |
| Mg ²⁷ | 0.57 | 12.160 | 0.03 | 2.47×10^{-10} |
| Al ²⁸ | 0.14 | 50.220 | 0.23 | 4.55×10^{-10} |
| S ³⁷ | 0.31 | 22.650 | 0.14 | 6.16×10^{-10} |
| Cl ³⁸ | 2.24 | 3.097 | 0.43 | 1.39×10^{-8} |
| K ⁴² | 44.64 | 0.155 | 1.10 | 7.09×10^{-7} |
| Ca ⁴⁹ | 0.53 | 13.120 | 1.10 | 8.38×10^{-9} |
| Ti ⁵¹ | 0.35 | 19.910 | 0.14 | 7.03×10^{-10} |
| V ⁵² | 0.23 | 30.640 | 4.90 | 1.60×10^{-8} |
| Mn ⁵⁶ | 9.29 | 0.746 | 13.30 | 1.76×10^{-6} |
| Cu ⁶⁴ | 46.44 | 0.144 | 4.50 | 3.12×10^{-7} |
| Br ⁸⁰ | 1.08 | 6.416 | 11.40 | 1.78×10^{-7} |
| Br ⁸² | 127.08 | 0.055 | 3.20 | 5.87×10^{-6} |
| I ¹²⁸ | 1.50 | 4.620 | 6.40 | 1.39×10^{-7} |

Table: 3.1- Compilation of nuclear properties of several isotopes (Ec-68).

The values for $\sigma\phi_0/\lambda$ given in table 3.1 been much less than unity, the time dependency of the equation 3.5 can be ignored. With these assumptions equation 3.2 can be solved to obtain A, an expression for activity,

$$A \equiv \lambda N = \sigma N_0 \phi_0 (1 - e^{-\lambda t}) e^{-\lambda t} \quad (3.6)$$

where t is the time of neutron flux exposure, t' is the time after removal from the neutron flux, ϕ_0 is the neutron flux, a constant.

When the activity is predominantly gamma radiation, both a qualitative and quantitative analysis can be performed through the use of a gamma-ray spectrometer system and the proper data analysis procedure.

The gamma-ray spectrum of a complex sample containing several isotopes (m in number) may be represented by:

$$y_i = \sum_{j=1}^m B_{ij} \quad , \quad i = 1, 2, 3, \dots, n \quad (3.7)$$

where, y_i is the counts registered in channel i of a MCA,
 B_{ij} is the counts registered in channel i from j^{th} constituent isotope

Reference gamma-ray spectra are collected for each of the m constituent isotopes using known quantities of radioisotopes. The reference gamma-ray spectra of each of the m constituent isotopes are related to the complex gamma-ray spectrum by,

$$y_i = \sum_{j=1}^m x_{ij} B_j + e_i \quad , \quad i = 1, 2, 3, \dots, m \quad (3.8)$$

where, B_j is the intensity coefficient for isotope j , i.e., a coefficient proportional to the quantity of isotope j in the complex sample,
 x_{ij} is the counts registered in channel i from the j^{th} isotope of known quantity,
 e_i is the deviation at point i of the estimated value of y_i ,
 $\sum_{j=1}^m x_{ij} B_j$, from the observed value, y_i ,

Since negative quantities of radioactivity are physically impossible, the following restriction is necessary:

$$B_j \geq 0 \quad , \quad j = 1, 2, 3, \dots, m \quad (3.9)$$

Equation 3.8 and inequality 3.9 form a constrained general linear hypothesis. Therefore, the problem is to estimate the intensity coefficients. The minimization of the sum of the squared error, i.e., least squares, generally yield estimators for the intensity coefficients with desirable statistical properties, e.g., least squares estimators are unbiased if the data are free of systematic errors (Ec-68).

The linear hypothesis, in matrix form, is

$$\underline{Y} = \underline{X}\underline{B} + \underline{e} \quad (3.10)$$

where,

- Y: (nx1) vector of dependent variable values,
- X: (nxm) matrix of independent variable values,
- B: (mx1) vector of solution values,
- e: (nx1) vector of errors (deviation values),
- n: number of data points (channels),
- m: number of reference isotopes.

The sum of the squared errors is given by:

$$e'e = (Y - XB)'(Y - XB) , \quad (3.11)$$

which reduces to,

$$e'e = Y'Y - 2B'X'Y + B'X'XB \quad (3.12)$$

The estimator, b for B results from minimizing e'e. If the elements of B are linearly independent, e'e is differentiated with respect to each element of B, term by term, the result is equated to zero, and arranged in matrix form after B is replaced by its estimator, \hat{B} . This results in;

$$(X'X) b = X'Y \quad (3.13)$$

The estimate, b, is obtained by solving equation 3.13 for b:

$$b = (X'X)^{-1}(X'Y) \quad (3.14)$$

Irrespective of the error probability distribution function, the estimates, b, for the true parameters, B, calculated from equation 3.14 minimize the sum of the squared error. Rewriting equations 3.13 and 3.14 and including the error variance-covariance matrix, the following equations are obtained (Gr-74)

$$(X'V^{-1}X) b = (X'V^{-1}Y) \quad (3.15)$$

and
$$b = (X'V^{-1}X)^{-1}(X'V^{-1}Y) \quad (3.16)$$

where V^{-1} is the inverse of the error variance-covariance matrix and is often denoted as a weighting matrix, W.

$$W \equiv V^{-1} \quad (3.17)$$

If the errors are assumed to be statistically independent, the covariance terms for e will be zero and estimate for the variance of e_i , $\sigma^2(e_i)$, can be obtained from,

$$\sigma^2(e_i) = \sigma^2(y_i) + \sum_{j=1}^m [B_j^2 \sigma^2(x_{ij}) + x_{ij}^2 \sigma^2(B_j)] \quad (3.18)$$

where $\sigma^2(x_{ij})$: the variance of the counts in channel i from isotope j
 $\sigma^2(B_j)$: variance of the estimator for the j^{th} isotope.

An iterative approach is used to solve for b , as follows:

1. Estimate $\sigma_1^2(e_i)$, given in equation 3.18 using only $\sigma^2(y_i)$,
2. Calculate the first estimate, b_1 , of b given in equation 3.16 using the estimated $\sigma_1^2(e_i)$ from step 1.,
3. Estimate $\sigma_2^2(e_i)$ by using b_1 for B and $\sigma^2(b_1)$ for $\sigma^2(b)$ in equation 3.18,
4. Calculate b_2 using the estimated $\sigma_2^2(e_i)$
5. Continue this procedure until a termination criterion is met.

This iterative procedure can be mathematically stated as,

$$b_n = (X^0 W_n X)^{-1} (X^0 W_n Y) \quad (3.19)$$

$$W_{in} = [\sigma_n^2(e_i)]^{-1} = \left\{ \sigma^2(y_i) + \sum_{j=1}^m [b_{j,n-1}^2 \sigma^2(x_{ij}) + x_{ij}^2 \sigma_{n-1}^2(b_j)] \right\} \quad (3.20)$$

$$b_{j_0} = \sigma_0^2(b_j) = 0 \quad (3.21)$$

The termination criterion used is:

$$\chi^2/DF \leq 1 + U_p \sqrt{(2/DF)} \quad (3.22)$$

where,
 χ^2 : Chi squared, i.e., $\sum_{i=1}^n \frac{[y_i - \sum_{j=1}^m x_{ij} b_j]^2}{\sigma^2(e_i)}$

U_p : Confidence interval factor,

DF: Degrees of freedom, i.e., $n-m-1$,

$/DF$; is called the quality of fit number.

U_p factor corresponds to a desired level of significance. Since the numerator of χ^2 is the minimization criterion used in this method, it is not reasonable to assume that the iteration procedure will increase χ^2/DF . Therefore, only a one-sided region is included in equation 3.22 (Ec-68). The mathematical model of this approach can be found in some detail in

The least squares is a logical method to determine the intensity coefficients for NAA problem. (Ak-77). However, one of the difficulties in applying least squares technique to the resolution of gamma spectra is that the calculated intensity coefficients can be negative or it may be difficult to determine the reference spectra that comprise the complex spectrum. For example, the least squares technique may yield positive intensity coefficients for all the reference spectra, but some may not be significant. To overcome those difficulties, the t value for the null hypothesis is calculated for each intensity coefficient (Gr-74). The t values then are used in the correlation algorithm to eliminate reference spectra for which the calculated intensity coefficients are negative or not significant. An outline of the correlation algorithm used in the code is given in figure 3.1. Some of the statistical techniques used in CORGAM, namely, properties of estimators, least squares technique, curve fitting method, interpolation and determinant evaluation methods are given in the appendices. And a more detailed logic diagram is also given in the appendix J. If needed a complete description of the methods used in CORGAM can be found in its manual.

This code, CORGAM has been used to analyze many gamma-ray spectra from several cereal and feed grain samples (Ec-68). Each sample was irradiated in the rotary specimen rack of a KSU TRIGA Mark II nuclear reactor at 100 kW power for 15 minutes. A small (0.02 gm) copper powder sample was irradiated with each grain sample. The area of the 0.511 MeV photopeak from ^{64}Cu was used for neutron fluence normalization. Each sample was transferred to a nonirradiated container and was counted at 10 cm from 3" x 3" NaI(Tl) scintillation crystal, which was interfaced to a TMC-4096 multiparameter system. The reference spectra required by CORGAM had been collected in an identical manner (Ec-68).

Due to the special purpose of cereal analysis problems, the reference spectra contained in the library of the code is formed from those spectra of the elements contained in agricultural samples. These reference isotopes are: ^{24}Na , ^{42}K , ^{56}Mn , ^{38}Cl , ^{80}Br , ^{82}Br , ^{49}Ca , ^{28}Al , ^{27}Mg , ^{64}Cu , ^{128}I , and ^{37}S . But it is possible to replace the reference spectra by the ones needed for the analysis of other mixtures. It should be remembered that the reference data collection and that of the complex one should be done in identical conditions




Figure: 3.1- Outline of correlation algorithm for resolution of gamma spectra.

3.3- EXPERIMENTAL PROCEDURE OF CORGAM

Previous use of CORGAM was for the routine analysis of agricultural samples. In order to perform this analysis, the experimental procedure used is given in this section (Ec-68). The input data of the same experiment is used to run the program. The listing of the input data and that of the output obtained from the sample problem run are given in the appendix L. During the data collection of the sample problem the following equipment were used.

1. Shielded counting arrangement: A lead shielded counting cavity with interior dimensions of 24" x 24" x 30" was constructed of 2" x 4" x 8" lead bricks, with exterior walls 4" thick, and top and bottom 2" thick. The interior of the cavity was lined with 0.025" thick cadmium sheet followed by 0.025" thick copper sheet. Shielded counting cavity and wall construction details are illustrated in figure 3.2. By using this cavity background is reduced by a factor of 20 % .
2. Detector assembly: The scintillation crystal used is Harshaw 3"x3" NaI(Tl), solid type crystal, Model 12512, serial CC362 integral line. It uses an RCA 8054 photomultiplier tube, the resolution of the crystal is 8.07 %.
3. Preamplifiers: Two main purposes of a preamplifier are to serve an impedance matching and an amplification function. The preamplifier used in this experiment was TMC model DS-13 with a gain of approximately 2.0.
4. Multichannel analyzer: The multichannel analyzer is used to sort and store voltage pulses according to their respective heights. The pulses are sorted into differential pulse height units, i.e., from v_0 to $v_0 + v$. That function is carried out by converting the pulse height to time and measuring the elapsed time from pulse base to pulse height. The analyzer used in this experiment was the TMC 4096 multiparameter system.
5. Neutron source: KSU TRIGA Mark II nuclear reactor was used as a source of thermal neutrons. A typical irradiation for neutron activation analysis is made in either the pneumatic transfer unit or in the rotary specimen rack. In table 3.2 the approximate neutron flux values per watt of power at those two locations are compared.

The block diagram of the NAA instrumentation system is shown in figure

3.3.

6. Reference spectra: Spectroscopically pure compounds were irradiated to collect the reference spectra used in this experiment. Table 3.3 is a compilation of information about each of these compounds.

| Neutron Energy | Pneumatic transfer unit n/cm ² -sec-watt | Rotary specimen rack n/cm ² -sec-watt |
|----------------|--------------------------------------------------------|-----------------------------------------------------|
| Thermal | 1.7×10^7 | 7.2×10^6 |
| Fast | 1.4×10^7 | 6.0×10^6 |

Table: 3.2- Neutron flux levels for KSU TRIGA Mark II nuclear reactor (EG-68)




Figure: 3.3- Block Diagram of the NAA instrumentation system.

Table: 3.3- Reference spectra information.(Ec-68)

| Element | Chemical Form | %H ₂ O Hydration | Reference isotope | Half Life (min) | Irradiation Policy " | Weight (gm) |
|---------|---------------------------------|-----------------------------|-------------------|-----------------|----------------------|-------------|
| Na | Na ₂ CO ₃ | 6.8 | Na ²⁴ | 900.0 | 1-4-2-10 | 0.229 |
| K | K ₂ CO ₃ | 5.1 | K ⁴² | 744.0 | 4-4-5-10 | 0.299 |
| Mn | metal | 0.0 | Mn ⁵⁶ | 154.8 | 2-4-2-0.1 | 0.0806 |
| Cl | NH ₄ Cl | 0.0 | Cl ³⁸ | 37.3 | 2-5-2-50 | 0.0130 |
| Br | NH ₄ Br | 0.0 | Br ⁸⁰ | 18.0 | 2-65-2-90 | |
| | | | Br ⁸² | 2118.0 | 2-1707-10 | 0.0304 |
| Ca | CaCO ₃ | 1.4 | Ca ⁴⁹ | 8.8 | 2-5-2-10 | 0.194 |
| Al | Al ₂ O ₃ | 3.9 | Al ²⁸ | 2.3 | 1-12-2-10 | 0.0628 |
| Mg | MgO | 5.5 | Mg ²⁷ | 9.5 | 2-5-2-10 | 0.0788 |
| Cu | powder | 0.0 | Cu ⁶⁴ | 774.0 | 1-112-3-10 | 0.0321 |
| I | KI | 0.3 | I ¹²⁸ | 25.0 | 2-5-2-50 | 0.0252 |
| S | powder | 0.0 | S ³⁷ | 5.1 | 2-5-2-50 | 0.2064 |

" The first number is the irradiation time in minutes. The second number is the decay time in minutes. The third number is the counting time in minutes. The fourth number is the reactor power level in kilowatts.

The reference spectra were collected using the following procedure; A small amount (≤ 1.0 gram) of each spectroscopically pure reference compound, and a small amount of copper powder were irradiated in the rotary specimen rack (RSR) of the KSU TRIGA Mark II nuclear reactor for the times given in table 3.3. After the respective times had elapsed, the reference spectra were collected for the counting time given in table 3.3, using the solid " x3" NaI(Tl) scintillation crystal.




Figure: 3.2- Shielded counting cavity dimensions and details (Ec-68).

Instabilities in the spectrometer system used to collect the gamma-ray spectra can cause shifts to occur in these spectra. To compensate for these inconsistencies, the gamma-ray photopeaks of the spectra were fit to a Gaussian function using a linearized least squares method. The mid-point corresponding to the gamma-ray energy was calculated from the relation, $X_0 = E/G$ where E is the gamma-ray energy in MeV. If a reference spectrum contains more than one gamma-ray photopeak, these mid-channel values were used to calculate the slope and intercept of the relation, $E = a_0 + GX_0$, where E is the gamma-ray energy in MeV and X_0 is the mid-channel value of the photopeak. However, if the reference spectrum contains only one photopeak, the intercept was assumed to be zero. Using these parameters, each spectrum was shifted to a slope and intercept of 14.5 keV/ch and zero, respectively. Each of the gamma-ray photopeaks in the shifted spectra was fitted to a Gaussian function (Ec-68).

A more detailed description of the procedure and various applications performed by using CORGAM can be found in the reference Ec-68.

4. DISCUSSION AND CONCLUSION

The importance of neutron activation analysis was mentioned in the previous subsections. It is also mentioned that the determination of the micro constituents of a complex sample was possible down to the ppm (parts per million) level. In the study of the principles of NAA, it is noticed that the interaction of gamma radiation with the detector material was very important. The gamma-rays resulting from the decay of activated nuclei, interact with the detector. If the material of the detector is a phosphor, scintillations are produced in the phosphor resulting from the interaction of the incoming photons with the detector material. These scintillations are transformed into electrical signals by the photomultiplier tube. The signals are then amplified measured and recorded by the use of proper electronic equipments. The most widely used electronic equipment is a multichannel analyzer, which gives a distribution of the incoming photons with respect to their energies. This amplitude versus energy distribution of the photons is called the gamma-ray spectrum of the activated nuclei. Due to the character of the interactions of gamma-rays with a particular scintillation crystal, produces a unique gamma-ray spectrum. The analysis of this spectrum by the methods described in this work, gives the necessary information about the nature and the concentration of the constituents of the complex sample. Among the methods used in common NAA work, the correlation algorithm is found to be the best method for the resolution of gamma-ray spectra.

A predesigned computer code which uses the correlation algorithm in the unfolding procedure of a complex gamma-ray spectrum is also analyzed. In this analysis it is observed that the program is restricted to a problem of 400 channels and 15 reference spectra. There would be no difficulty in expanding this to as many as 4096 channels. However due to the accuracy desired in the inversion procedure, the number of reference isotopes is limited. CORGAM is a novel method of unfolding complex gamma-ray spectra because only reference spectra that have intensity coefficients which are significant at a preselected level are retained in the output and the variance calculations include contributions from both the reference and complex spectra. The large computing time required by the NAA can be eliminated by the use of a computer program. A typical analysis, involving 256 channels and 14 reference spectra, requires less than 1 minute of execution time on an IBM 360/50. Although the computer memory requirements are large for the correlation

algorithm as it was programmed for the IBM 360/50, it is possible to segment the program for adaptation to a smaller computer. In order to use this program in NAA, changes should be made on the reference spectra and on the flux reference spectra according to the laboratory conditions. The choice of the spectrometr system affects the results in various ways. When a NaI(Tl) crystal is used the interferences of the various reference spectra became significant. Use of a Ge(Li) semiconductor detector eliminates these interferences to a desired level.

Presently the code contains eleven reference spectra in its library, if needed these spectra can be changed by the ones required for a specific analysis. Since the number of reference spectra is limited by 15 isotopes, at each special problem, need may arise to change the reference library. It is concluded that this procedure was very time consuming when the spectrometer system is not branched to the computer where the program will be run.

REFERENCES

- Ad-69 Adams F, Dams R, "Radiochim. Acta", 10, 1 (1968), "Radioanalyst. Chem", 3, 99 (1969).
- Al-61 Alimarin I.P., Yakovlev Y.V. "Industrial Laboratory", USSR, 26/8 (1961).
- Ak-77 Aktas Z. "Numerical Analysis", ODTÜ, (1977).
- Ay-62 Ayres F.Jr. "Matrices", Shaum Pub.Co. New York (1962).
- Ba-63 Baker P.S., Fucellio D.A., Gerrard M.W. "Radioisotopes in industry" (1963).
- Bi-53 Birks J.B. "Scintillation Counters", Mc Graw Hill, (1953).
- Br-65 Brownlee K.A. "Scintillation Theory and Methodology in Science and Engineering", 2.ed., J.W. and sons, (1965).
- Ch-60 Cheronis N.D. "Microchemistry", 5.vol. (1960).
- Ch-61 Chase R.L. "Nuclear Pulse Spectrometry", Mc Graw Hill book Co. New York (1961).
- Cr-64 Crouthamel C.E. "Applied gamma-ray Spectrometry", Pergamon Press inc. New York (1964).
- Cu-53 Curran S.C. "Luminescence and Scintillation Counters", Academic New York (1953).
- Cu-68 Curie L.A. "Analytical chemistry", 40, 587. (1968).
- Ec-68 Eckhoff N.D. "A Correlation Algorithm for Gamma Ray Spectra", Kansas state Universty, (1968).
- Ec-74 Eckhoff N.D. "A correlation Algorithm for Gamma ray Spectra" Nuclear Instruments and Methods, North-Holland Pub. Co. (1974).
- Er-64 Erwall L.G., Fdrsborg H.G., Junggren K. "Industrial Isotopes Techniques", (1964).
- Ev-55 Evans R.D. "The Atomic Nucleus", Mc Graw Hill, (1955).
- Fa-53 Fano U. "Interactions of Gamma Rays", Nucleonics, Vol. 8, 9, 10 August-October (1953).
- Fr-67 Friedman A.K. "Nuclear Physics", Kogakusha LTD. (1967).
- Ga-67 Gardner R.P., Ely R.L.Jr. "Radioisotope Measurement Applications in Engineering", Reinhold Pub.Co. (1967).
- Gr-74 Groot M.H.De "Probability and Statistics", Addison Wesley Pub.Co. (1974).
- Gr-61 Graybill F.A. "An Introduction to Linear Statistical Models", Mc Graw Hill, New York, (1961).
- He-64 Heat R.L. "Scintillation Spectrometry-Gamma Ray Spectrum Catalog" 2.ed. Vols 1, 2 IDO, 16880, (1964).
- He-36 Hevesey G., Levi H., "Mathematics, Physics, Medicine", 14(5). (1936).

- Ho-77 Holman J.P. "Experimental Methods for Engineers", Mc Graw Hill Kogakusha LTD. (1977).
- IBM IBM, manuel H20-0166-2 System 360 Scientific Subroutine Package Application Description (360 A-CM-03X).
- Jo-74 Jordan C. "Calculus of Finite Differences", New York Chelsea (1974).
- Ke-72 Keisch B. "Neutron Activation Analysis", ERDA, (1972).
- Ko-61 Kohl J., Zentner R.D., Lukens H.R. "Radioisotope Applications Engineering", D. Van Nostrand Co. (1961).
- La-72 Lamarsh J.R. "Introduction to Nuclear Reactor Theory", Addison Wesley Pub. Co. (1972).
- Le-68 Lederer C.M., Hollander J.M., Perlman I. "Table of Isotopes", /6.ed. John Willey Co. (1968).
- Mo-63 Mood A.M., Graybill F.A. "Introduction to the Theory of Statistics", Mc Graw Hill Co. (1963).
- NEA-CPL NEA-Computer Program Library, Argonne Code Center ANL-7411, Ispra-Italy.
- OK-62 O'Kelley G.D. "Applications of Computers to Nuclear and Radiochemistry", Proc. Symp. (Gatlinburg, Tennessee), (1962).
- Pi-64 Piraux H. "Radioisotopes and their Industrial Applications", Philips Technical Lib. (1964).
- Pi-62 Pinta M. "Recherche et Dosage des Elements Traces", Dunod, Paris (1962).
- Po-62 Pontryagin L.S., et al., "The Mathematical Theory of Optimal Processes", John Wiley Co. (1962).
- Pr-64 Price W.J. "Nuclear Radiation Detection", Mc Graw Hill 2.ed. (1964).
- Se-38 Seaborg G.T., Livingood J.J., "J. Am. Chem. Soc.", 1784-786, (1938).
- Sm-66 Smith L.H. "Statistical Technics in Neutron Activation Analysis", Analytica Chimica Acta, 36 (1966).
- So-69 Soete D. De, Gijbels R., Hoste J. "Neutron Activation Analysis", Wiley Interscience, Chemical Analysis Vol: 34, (1969).
- Tr-67 Trombka J.I. "A Method for the Analysis of Pulse-Height Spectra Containing Gain Shift and Zero Drift Compensation", Transactions of the ANS, Vol. 10, no: 1, (1967).
- Wa-71 Warren J. Mc G., Gordon, Breach, "Non Destructive Testing", 2.ed. Science Pub., (1971).
- Yo-73 Young D.M., Gregory R.T., "A Survey of Numerical Mathematics", Addison Wesley Pub. Co. (1973).

A P P E N D I C E S

APPENDIX A

-PROPERTIES OF ESTIMATORS

The estimate of the intensity coefficients (concentrations) of the equation 3.14 have the following properties;

- a) Unbiased (Gr-61): An estimator, a , for the parameter α is defined as an unbiased estimator of α if

$$E(a) = \alpha \quad (A-1)$$

where E represents the mathematical expectation. The mathematical expectation is defined by the following theorem (Ho-63).

Theorem: Let \underline{x} be a random variable with density $f(x)$. The expected value of a function u of the random variable \underline{x} is,

$$E(u(x)) = \sum_x u(x_1) f(x_1) \quad (A-2)$$

if \underline{x} is discrete, and

$$E(u(\underline{x})) = \int_{-\infty}^{\infty} u(x) f(x) dx \quad (A-3)$$

if \underline{x} is continuous.

Therefore, the mathematical expectation is just the mean value of the function u of the random variable \underline{x} .

- b) Minimum variance (Gr-61): An estimator, a , is defined to be the minimum variance estimator of the parameter α if

$$E((a - E(a))^2) \leq E((a' - E(a'))^2) \quad (A-4)$$

The term $E((a - E(a))^2)$ is defined to be the variance of the parameter α .

- c) Sufficiency (Ho-63): Let $f(\underline{x}, \alpha)$ be a joint frequency function involving the parameter α and the random variable \underline{x} , n in number. Let $a_1, h_1(\underline{x}), a_2, h_2(\underline{x}), \dots, a_k, h_k(\underline{x})$ be a set of statistics and $a'_1, h'_1(\underline{x}), a'_2, h'_2(\underline{x}), \dots, a'_k, h'_k(\underline{x})$ be any other set of estimators which are not functions of a_1, a_2, a_k . If $g(a'/a)$, the conditional density function of a' given a , does not involve α , the set a is called a set of sufficient statistics. This means that "all the information" contained in the n values of \underline{x} is given

by the k parameters \underline{a} where $k < n$.

- d) Consistency (Mo-63): If \underline{a} is a sequence of estimators of α , they are simple consistent estimators, if for every $\xi > 0$

$$\lim_{n \rightarrow \infty} P(|\bar{a}_n - \alpha| \leq \xi) = 1 \quad (\text{A-5})$$

Further, the estimators, \underline{a} , are squared error consistent estimators of α , if

$$\lim_{n \rightarrow \infty} \left\{ E((\bar{a}_n - \alpha)^2) - \frac{1}{n} \sum_1^k (a_i - \alpha)^2 \right\} = 0 \quad (\text{A-6})$$

where

$$\bar{a}_n = \frac{1}{n} \sum_1^k a_i$$

Squared error consistency implies that both the bias and variance of \underline{a} approach zero as the number of sample points increase.

- e) Efficiency (Br-65): efficiency of an estimator is defined to be the ratio of the variance of the minimum variance estimators to the variance of the estimator in question. Thus, minimum variance estimators are said to be efficient estimators .

APPENDIX B-THE LEAST SQUARES METHOD

The method is used to determine the intensity coefficients for the NAA problem. To illustrate the method, consider a set of n points, (x, y) . The points are assumed to be related by the linear hypothesis (Gr-74):

$$y_1 = B_0 + B_1 x_1 + e_1, \quad 1, 2, 3, \dots, n \quad (\text{B-1})$$

where e_1 is the deviation at point 1 of the estimated value, from its observed value, y_1 , B_0, B_1 are the true values of the intercept and slope, respectively. The estimators, b_0 and b_1 , for B_0 and B_1 can be obtained by minimizing the sum of the squared errors;

$$S = \sum_{i=1}^n e_1^2 = \sum_{i=1}^n (y_1 - (B_0 + B_1 x_1))^2 \quad (\text{B-2})$$

S is minimized by differentiating with respect to B_0 , and B_1 and equating the results to zero. Also, b_0 and b_1 , the estimators, are substituted for B_0 and B_1 . The resulting equations are;

$$\sum_{i=1}^n (y_1 - (b_0 + b_1 x_1)) = 0 \quad (\text{B-3})$$

$$\sum_{i=1}^n \{x_1 (y_1 - (b_0 + b_1 x_1))\} = 0 \quad (\text{B-4})$$

Equations B-3 and B-4 can be rearranged into a form called the normal equations.

$$b_0 n + b_1 \sum_{i=1}^n x_1 = \sum_{i=1}^n y_1 \quad (\text{B-5})$$

$$b_0 \sum_{i=1}^n x_1 + b_1 \sum_{i=1}^n x_1^2 = \sum_{i=1}^n y_1 x_1 \quad (\text{B-6})$$

Solving equations B-5 and B-6 for b_0 and b_1 , results in

$$b_1 = \frac{\sum_{i=1}^n x_1 y_1 - ((\sum_{i=1}^n x_1)(\sum_{i=1}^n y_1))/n}{\sum_{i=1}^n x_1^2 - (\sum_{i=1}^n x_1)^2/n} = \frac{\sum_{i=1}^n (x_1 - \bar{x})(y_1 - \bar{y})}{\sum_{i=1}^n (x_1 - \bar{x})^2} \quad (\text{B-7})$$

$$b_0 = \bar{y} - b_1 \bar{x} \quad (\text{B-8})$$

where,

$$\bar{y} = \left(\sum_{i=1}^n y_i \right) / n, \quad \bar{x} = \left(\sum_{i=1}^n x_i \right) / n.$$

Equation B-1 can be rewritten in terms of the estimator, b_1 , and the mean value \bar{x} and \bar{y} , through the use of equation B-8.

$$y_i = \bar{y} + b_1(x_i - \bar{x}), \quad i = 1, 2, 3, \dots, n \quad (\text{B-9})$$

Equation B-9 can be rewritten in matrix form as,

$$\underline{Y} = \underline{X}b_1, \quad (\text{B-10})$$

where $\underline{Y} = (y - \bar{y}) = (n \times 1)$ vector of the dependent values

$\underline{X} = (x - \bar{x}) = (n \times m)$ matrix of the independent variables.

In matrix form, the linear hypothesis model, equation B-1, is more appropriately written as,

$$\underline{Y} = \underline{B}^0 \underline{X}^0 + \underline{e} \quad (\text{B-11})$$

where \underline{X}^0 is the transpose of the $(n \times m)$ matrix X

\underline{B}^0 is the transpose of the $(m \times 1)$ solution vector B

\underline{e} is the $(n \times 1)$ vector of deviations.

More often equation B-11 is written in the following form:

$$\underline{Y} = \underline{X}B + \underline{e} \quad (\text{B-12})$$

APPENDIX C-CURVE FITTING METHODS

The following method was used to fit the photopeaks to a Gaussian function (Ec-68).

First the photopeaks were assumed to be pure Gaussian, i.e.,

$$y_1 = y_0 \exp(-(x_1 - x_0)^2 / b_0) , \quad (C-1)$$

y_1 is the counts in channel 1 of a photopeak,

y_0 is the count corresponding to channel x_0 ,

x_0 is the mean channel of the Gaussian function

x_1 is the channel number

$$b_0 = 2\sigma^2,$$

σ^2 is the variance of the Gaussian function.

Equation C-1 is non linear and must be linearized to apply the principle of least squares. The logarithm of equation C-1 yields a linear form.

$$\ln y_1 = \ln y_0 - \frac{x_0^2}{b_0} + \frac{2x_0}{b_0} x_1 - \frac{1}{b_0} x_1^2 \quad (C-2)$$

Equation C-2 can be rewritten in terms of new variables,

$$z_1 = a_1 + a_2 x_1 + a_3 x_1^2 \quad (C-3)$$

where,

$$z_1 \equiv \ln y_1 ,$$

$$a_1 \equiv \ln y_0 - \frac{x_0^2}{b_0} , \quad a_2 = \frac{2x_0}{b_0} .$$

$$a_3 \equiv -\frac{1}{b_0}$$

The principle of least squares applied to equation C-3 yields the following matrix equation

$$a = (X^0 W X)^{-1} (X^0 W Z) , \quad (C-4)$$

where a is a (3×1) vector of the coefficients a_1, a_2 and a_3 ,

X is a $(n \times 3)$ matrix whose first column is composed of all unity values, the second column is composed of the x_1 values, and the third column is

composed of the x_2 values

X' is the transpose of X ,

Z is a $(n \times 1)$ vector of the z_1 values,

W is a $(n \times n)$ diagonal weighting matrix whose elements are the inverse of the variances of the z_1 values,

n is the number of data points.

The variance of z_1 can be obtained by

$$\sigma^2(z_i) = \frac{1}{y_i^2} \sigma^2(y_i) = \frac{1}{y_i} \quad (C-5)$$

It should be noted that by applying the principle of least squares to equation C-3 instead of to equation C-1 the minimum total squared deviation is,

$$S_D = \sum_{i=1}^n \left\{ y_i - y_0 \exp \left[-(\lambda_i - \lambda_0)^2 / b_0 \right] \right\}^2 \quad (C-6)$$

The desired minimum total deviation is,

$$S_L = \sum_{i=1}^n \left[\ln y_i - \left(\ln y_0 - \frac{\lambda_0^2}{b_0} + \frac{2\lambda_0}{b_0} \lambda_i - \frac{1}{b_0} \lambda_i^2 \right) \right]^2 \quad (C-7)$$

Minimizing S_L does not guarantee a minimum for S_D . To minimize S_D requires an iterative technique which employs a Taylor's series expansion and can result in unstable solutions for cases with variations in the data. In addition this method requires a significantly longer amount of computation time.

Variance values of y_0 , x_0 , and b_0 , could be obtained by

$$\sigma^2(y_0) = \left[C_{4,4} + \frac{a_2^2}{4a_3^2} \left(C_{2,2} + \frac{a_2^2}{4a_3^2} C_{3,3} \right) \right] \exp \left[2 \left(a_1 - \frac{a_2^2}{4a_3} \right) \right] \quad (C-8)$$

$$\sigma^2(x_0) = \frac{1}{4a_3^2} \left[C_{2,2} + \left(\frac{a_2}{a_3} \right)^2 C_{3,3} \right] \quad (C-9)$$

$$\sigma^2(b_0) = C_{3,3} / a_3^4, \quad (C-10)$$

where,

a_1, a_2, a_3 , are the elements of the vector a , and $C_{4,4}, C_{2,2}, C_{3,3}$ are the diagonal elements of the $(X'WX)^{-1}$ matrix. Since the y_0 value is just a normalization term, it could also be calculated by minimizing S_D after x_0 and b_0 were obtained by minimizing S_L . This method yielded lower variance values for y_0 (Gr-74). Hence, y_0 and (y_0) were calculated from,

$$y_0 = \frac{\sum_{i=1}^n w_i y_i k_i}{\sum_{i=1}^n w_i k_i^2}, \quad (C-11)$$

$$\sigma^2(y_0) = \left[\sum_{i=1}^n w_i k_i^2 \right]^{-1} \quad (C-12)$$

where

n is the number of data points used in the fit,

$$k_i = \exp(-(x_i - x_0)^2 / b_0)$$

In deriving equation C-12 it was assumed that $w_i = \frac{1}{y_i}$.

APPENDIX D-LAGRANGIAN INTERPOLATION

This interpolation method is used in the subprogram REFIN, to interpolate between the decay spectra of a reference isotope.

The method is based on polynomial interpolating functions (Ak-77). Assume that a function $y = f(x)$ is given at $(n-1)$ points, say $x_0, x_1, x_2, \dots, x_n$, and determine a polynomial $F(x)$ of degree n or less such that

$$F(x_i) = f(x_i) \quad i = 0, 1, 2, \dots, n \quad (D-1)$$

It may be proved that the polynomial exists and it is unique. For this purpose two alternative proofs will be given.

Proof 1:

Let $F(x)$ function is written as

$$F(x) = L_0 y_0 + L_1 y_1 + \dots + L_n y_n \quad (D-2)$$

where

$$y_0 = f(x_0), \quad y_1 = f(x_1), \quad \dots, \quad y_n = f(x_n) \quad (D-3)$$

and

$$L_i(x) = \prod_{\substack{j=0 \\ j \neq i}}^n \left(\frac{x - x_j}{x_i - x_j} \right) \quad i = 0, 1, 2, \dots, n \quad (D-4)$$

It is clear that each $L_i(x)$ term is a polynomial of degree n or less and rewriting it explicitly as

$$L_i(x) = \frac{(x-x_0)(x-x_1)(x-x_2)\dots(x-x_{i-1})(x-x_{i+1})\dots(x-x_n)}{(x_i-x_0)(x_i-x_1)(x_i-x_2)\dots(x_i-x_{i-1})(x_i-x_{i+1})\dots(x_i-x_n)} \quad (D-5)$$

it is obvious that $f(x_j) = F(x_j)$ because for $j \neq i$, where $j = 0, 1, \dots, n$, all $L_i(x_j)$ terms are zero. For $j=i$ case $L_i(x_j) = 1$ and therefore,

$$E(x_i) = L_0(x_i)y_0 + L_1(x_i)y_1 + L_2(x_i)y_2 + \dots + L_n(x_i)y_n \quad (D-6)$$

$$\text{becomes } F(x_i) = y_i \quad \text{where } y_i = f(x_i) \quad (D-7)$$

Proof 2:

An alternative method to determine the interpolating polynomial will be given using the method of undetermined coefficients.

(a) being any convenient value, then

$$F(x_1) = C_0(x-a)^n + C_1(x-a)^{n-1} + \dots + C_n \quad (D-8)$$

The coefficients C_0, C_1, \dots, C_n will be such that

$$F(x_i) = f(x_i) \quad \text{or} \quad F(x_i) = y_i \quad (D-9)$$

is satisfied. Obviously, C_k coefficients must satisfy the following system of linear algebraic equations,

$$\begin{aligned} C_0(x_0-a)^n + C_1(x_0-a)^{n-1} + \dots + C_n &= y_0 \\ C_0(x_1-a)^n + C_1(x_1-a)^{n-1} + \dots + C_n &= y_1 \\ \vdots & \\ C_0(x_n-a)^n + C_1(x_n-a)^{n-1} + \dots + C_n &= y_n \end{aligned} \quad (D-10)$$

This system of $n+1$ equations in terms of $n+1$ unknowns will have a solution if the coefficient matrix is nonsingular. The determinant of the coefficient matrix, which is known as Vandermonde determinant, is written as

$$\Delta = \begin{vmatrix} (x_0-a)^n & (x_0-a)^{n-1} & \dots & (x_0-a) & 1 \\ (x_1-a)^n & (x_1-a)^{n-1} & \dots & (x_1-a) & 1 \\ \vdots & \vdots & \ddots & \vdots & \vdots \\ (x_n-a)^n & (x_n-a)^{n-1} & \dots & (x_n-a) & 1 \end{vmatrix} \quad (D-11)$$

and its value is given as

$$\Delta = \prod_{\substack{i,j=0 \\ i < j}}^n (x_i - x_j) \quad (D-12)$$

Since x_i are distinct, it follows that $\Delta \neq 0$ and C_k can be uniquely determined. Quite often in computations $a=0$ is taken ($J_0=74$).

APPENDIX E

-DETERMINANT EVALUATION BY CHIO'S METHOD (PIVOTAL CONDENSATION)

This method is used in the subprogram DET in order to evaluate the determinant of the (X^0WX) matrix.

This method continuously reduces the order of the determinant by one, using a series of subtractions (Ak-77). For the explanation of the method suppose a $(4,4)$ matrix A is given.

$$|A| = a_{11} \begin{vmatrix} 1 & a'_{12} & a'_{13} & a'_{14} \\ a_{21} & a_{22} & a_{23} & a_{24} \\ a_{31} & a_{32} & a_{33} & a_{34} \\ a_{41} & a_{42} & a_{43} & a_{44} \end{vmatrix} = a_{11} \begin{vmatrix} 1 & 0 & 0 & 0 \\ a_{21} & a''_{22} & a''_{23} & a''_{24} \\ a_{31} & a_{32} & a_{33} & a_{34} \\ a_{41} & a_{42} & a_{43} & a_{44} \end{vmatrix}$$

1 2 3

where

$$a'_{12} = \frac{a_{12}}{a_{11}}, \quad a'_{13} = \frac{a_{13}}{a_{11}}, \quad a'_{14} = \frac{a_{14}}{a_{11}}$$

in the second step, performing the following operations,

- 1: $(2^{\text{nd}} \text{ col.}) - a'_{12} \times (1^{\text{st}} \text{ col.})$
- 2: $(3^{\text{rd}} \text{ col.}) - a'_{13} \times (1^{\text{st}} \text{ col.})$
- 3: $(4^{\text{th}} \text{ col.}) - a'_{14} \times (1^{\text{st}} \text{ col.})$

yields a reduction in the

order of the determinant;

$$|A| = a_{11} \begin{vmatrix} a''_{22} & a''_{23} & a''_{24} \\ a''_{32} & a''_{33} & a''_{34} \\ a''_{42} & a''_{43} & a''_{44} \end{vmatrix} = a_{11} a''_{22} \begin{vmatrix} 1 & a'''_{23} & a'''_{24} \\ a''_{32} & a''_{33} & a''_{34} \\ a''_{42} & a''_{43} & a''_{44} \end{vmatrix} = a_{11} a''_{22} \begin{vmatrix} 1 & 0 & 0 \\ a''_{32} & a'''_{33} & a'''_{34} \\ a''_{42} & a'''_{43} & a'''_{44} \end{vmatrix}$$

- 1: $(2^{\text{nd}} \text{ col.}) - a'''_{23} \times (1^{\text{st}} \text{ col.})$
- 2: $(3^{\text{rd}} \text{ col.}) - a'''_{24} \times (1^{\text{st}} \text{ col.})$

$$|A| = a_{11} a''_{22} \begin{vmatrix} a'''_{33} & a'''_{34} \\ a'''_{43} & a'''_{44} \end{vmatrix} = a_{11} \cdot a''_{22} \cdot \left[(a'''_{33} \times a'''_{44}) - (a'''_{43} \times a'''_{34}) \right]$$

The procedure is repeated till a $(2,2)$ determinant is finally reached.

(Yc-73, Ay-62).

APPENDIX F-DESCRIPTION OF THE PROGRAMS

CORGAM is the main program through which the flow of the program is controlled. This program is used to:

1. Read the program flow parameters.
2. Define many of the data parameters.
3. Read the input data which describes the experimental conditions, such as irradiation time, counting times, etc.
4. Read the complex spectrum, i.e., the spectrum to be analyzed.
5. Subtract background from the complex spectrum, normalize it, and shift it to compensate for electronic instabilities in the MCA.
6. Calculate the irradiation and decay time correction factors.
7. Form the $(X^T W X)$ matrix.
8. Calculate the concentration estimators.
9. Delete any and all non-significant reference isotopes.
10. Print or punch the output information.

The following 19 subprograms are called by CORGAM, their description are given in the following pages, and they are listed in order called.

- | | |
|------------|------------|
| 1. INIT | 17. CHI |
| 2. PHOTO | 18. VARNCE |
| 3. SMINV | 19. WRIT1 |
| 4. INDEK | |
| 5. SHIFT | |
| 6. SHIFTS | |
| 7. WATE | |
| 8. REFS | |
| 9. REFIN | |
| 10. XFORM | |
| 11. PEAKS | |
| 12. SAVE | |
| 13. DET | |
| 14. DELETE | |
| 15. STAT | |
| 16. WRIT | |

The following is the listing of the subprograms called by CORGAM. The description of the subprograms are also given together with the variables used and their options.

1. INIT:

This subprogram is used to initialize the irradiation and decay correction factors and many of the program control parameters.

```
SUBROUTINE INIT (NITER,RAT,SAV,C,ISOCHN,MODE,NBKGD,NUBKG,NSHFT,NPK,LPP,NORM,
NCTIN,NAREA,NCHCU,NCULL,NCUUL,NRHO,NUMVAR,NVAR,NREFS,ITER,NR,NW,LOCHN,NHICHN,
NCHAN,MC,MVAR,T,AREA,BKGTIM,TIMIRR,PHRL,UP,XNCOM)
```

```
DIMENSION C(15)
```

NITER: Number of iterations used in variance calculations

RAT: Correction factor

SAV: Correction factor

C(I): Correction factor for the Ith reference isotope, if -1.0 is entered for C(I) value, the Ith isotope is deleted from the analysis.

ISOCHN: Number of channels in the complex spectrum.

MODE: 1 Read by 5E14.8 format

2 Read by 6X,12F6.0,2X format

NBKGD: 0 no background subtraction

Other will subtract background

NUBKG: 0 no new background will be read

Other, a new background will be read

NSHFT: 0, no spectrum shifting

Other, will shift spectrum

NPK: Number of peaks in complex spectrum used in shifting

LPP: Index of peak used in shifting ratio

NORM: 0, no normalization of complex spectrum is needed

Other, normalize the complex spectrum

NCTIN: 0, read precalculated correction factors

Other, read unity (or less than zero for deletion)

NAREA: 0, read photopeak area for copper reference

1, calculate pp area for Cu ref., read Cu spectrum (read any)

2, use same area as previous spectrum.

NCHCU: Number of channels in Cu reference spectrum.
 NCULL: Lower fit limit on Cu reference photopeak.
 NCUUL: Upper fit limit on Cu reference photopeak.
 NRHO: 0, will not punch shifted output
 Other, will punch shifted output
 NUMVAR: Number of variances to be averaged
 NVAR: 1, read variance by 5E14.8 format
 2, variance: 1.0
 3, variance: RHO
 4, use variance previously read
 NREFS: 0, new set of reference spectra are read
 1, no new reference spectra are to be read
 ITER: Maximum number of iterations on variance
 NR: 0, calculate average weight variance from the long way
 Other, ratio: YAVE/YAVEAEVE
 NW: 0, will not punch shifted output for reference spectra
 Other, will punch shifted output for reference spectra
 LOCHN: Lower channel of least squares fit
 NHICHN: Upper channel of least squares fit
 NCHAN: Number of channel in variance spectrum if variance are to be read
 MC: Number of reference spectra which have a set of decay spectra. These
 spectra must be in order and they must be the first of the reference
 spectra to be read in.
 MVAR: 0, variance of reference spectrum: reference spectrum
 Other, read variance of reference spectrum
 T: Student's t value for null hypothesis test
 AREA: Copper reference photopeak normalized area
 BKGTIM: Counting time of background (minute)
 TIMIRR: Complex sample irradiation time (minute)
 FRHL: Flux reference half-life (minute)
 UP: U_p factor of Chi-squared test (confidence interval factor)
 XNCDM: Complex spectrum normalization constant.

This subprogram is read only for the complex spectrum whose data are not given in the beginning of the main program. Generally the first complex spectrum data is given in DATA statement and in NAMELIST. The complex spectra variables following the first one are initialized by this subprogram.

2. PHOTO:

This program is used to fit the photopeak data of a spectrum to a linearized form of the Gaussian function. The fitting technique is discussed in the Appendix C.

```
SUBROUTINE PHOTO (RES,XO,IPPLL,IPPUL,EPP,AREA)
DIMENSION B(3),Q(15,15),RES(400),RE(3),BZ(60),C(10)
```

RES(I): Counts registered at the Ith channel of the flux reference spectrum.

XO: x_0 (x sub o) of the equation C-1 of appendix C.

IPPLL: NCULL(defined in INIT)

IPPUL: NCUUL(defined in INIT)

EPP: Midpoint energy of the photopeak of interest (MeV).

AREA: defined in INIT

3. SMINV:

This subprogram is the scientific subroutine designed for the IBM 360 to invert symmetric matrices. It is modified by the author by its UNIVAC version in order to use in Boğaziçi University computer.

```
SUBROUTINE SMINV (R,N)
DIMENSION A(230),R(15,15)
```

R(I,J): Ith row ,Jth coloumn element of the matrix

N: Order of the matrix

The subprogram is used to obtain the inverse of a symmetry matrix A of order N. Only the upper triangular part of A is stored. On exit the inverse replaces A in storage. The inverse matrix is found by successive congruent transformations to obtain the diagonal. Inverting the diagonal elements and performing the congruent transformations in reverse gives the inverse matrix. In this code A is the (X^0WX) matrix.

4. INDEX:

This subprogram is apart of the subroutine SMINV and is used to store a two-dimensional matrix as a row vector. It is applicable only for symmetrical matrices. It allows reference to a two-dimensional matrix with subscripts (I,J) while the matrix is stored as a vector.

FUNCTION INDEX (I8,J8,N)

I8: Row of matrix

J8: Coloumn of matrix

N: Order of matrix (maximum value 15)

While a matrix of order N has $N \times N$ elements, the same matrix stored as a vector will have $(N^2 + N)/2$ elements, provided that the matrix is symmetric. In this way storage space is reduced.

5. SHIFT:

This program is used to calculate the input parameters to SHIFTS.

SUBROUTINE SHIFT (NCHN,NPK,LPP,NWRITE,Y,PO)

NCHN: ISOCHN (defined in INIT)

NPK: defined in INIT

LPP: defined in INIT

NWRITE: NRHO (defined in INIT)

Y(I): RHO(I), counts registered at the I^{th} channel of the complex spectrum
($I_{\text{max}}: 400$, i.e., up to 400 channels can be used).

PO: Identification number of the spectrum to be shifted.

6. SHIFTS:

This program is used to shift a spectrum to a new gain and baseline, hence it will compensate for electronic instabilities in the MCA.

SUBROUTINE SHIFTS (SPECT,S,RATIO,NCHN,NCHN2,BLSHFT)

DIMENSION SPECT(400),S(400)

SPECT(I): Counts registered at the I^{th} channel of the unshifted spectrum.

S(I): Counts registered at the I^{th} channel of the shifted spectrum by using
SPECT(I).

RATIO: Shifting ratio

NCHN: ISOCHN (defined in INIT), unshifted channel number.

NCHN2: The shifted new channel number.

BLSHFT: Baseline shift.

7. WATE:

This program is used to establish the original form of the weighting matrix. It calculates or reads the weight factors and reads or sets variances using (A) as storage.

```
SUBROUTINE WATE (NVAR, NUMVAR, NCHAN, RHO, XNCOM, III, OMG, A, RES, ISOCHN)
DIMENSION RES(400), RHO(400), OMG(400), A(400)
```

NVAR: Defined in INIT
 NUMVAR: Defined in INIT
 NCHAN: Defined in INIT
 RHO(I): Counts registered at the Ith channel of the complex spectrum.
 XNCOM: Defined in INIT
 III: Identification number of the complex spectrum (if several complex spectra are to be analyzed).
 OMG(I): Ith element of the variance vector.
 A(I): Ith element of the vector X in (X⁰WX) matrix.
 RES(I): Defined in PHOTO
 ISOCHN: Defined in INIT.

8. REFS:

This program is used to read all pertinent reference spectra concerning each spectrum. Options are provided for each spectrum to have a background subtraction, to be normalized, and to be corrected for decay. In addition, a set of reference spectra (up to five) for an element which upon irradiation produces several non separable isotopes, such as Br Br⁸⁰, Br⁸², may be read in. This program uses these data by interpolation to obtain the reference spectrum at the decay time of the complex spectrum.

```
SUBROUTINE REFS (NISO, NHICHN, ISO, XNORM, O, AT, BKG, NW, DECAY, COUNT, III, NFAC,
FRRHL, D, H, MC, NO, RHO, A)
DIMENSION ISO(15), H(15), D(15), XNORM(15), O(15), AT(400, 15), BKG(400), RES(400),
P(15), CALC(400), NFAC(15), CT(400, 15), T(5), C(5), RHO(400), A(400)
```

NISO: Number of reference isotopes used in the analysis of the complex spectrum
 NHICHN: Defined in INIT
 ISO(I): Number of channels in the Ith reference spectrum
 XNORM: Normalization constant
 O(I): Identification name of the reference spectrum

AT(I,J): Counts registered at the I^{th} channel of the J^{th} reference spectrum
 BKG(I): Counts registered at the I^{th} channel of the background spectrum.
 NW: Defined in INIT
 DECAY: Decay time for complex spectrum (minute)
 COUNT: Counting time for complex spectrum (minute)
 III: Defined in WATE
 NFAC(I): Decay time specification for the I^{th} reference isotope
 FRRHL: Reference spectrum, flux reference isotope half-life (minute)
 D(I): Irradiation time, in minutes, of reference sample spectrum
 H(I): Half-life, in minutes, of reference isotope spectrum
 MC: Defined in INIT
 NO: Number of decay spectra of the reference isotope (max 5)
 RHO(I): Defined in WATE
 A(I): Defined in WATE

9. REFIN:

This program is used to interpolate, channel by channel between those spectra discussed in REFS, by Lagrangian method. The method is discussed in the appendix D.

SUBROUTINE REFIN (AT,III,T,C,CT,NO,M,DECAY,COUNT,ISO,RHO,A)
 DIMENSION AT(400,15),CT(400,15),T(5),C(5),ISO(15),IS(15),RHO(400),A(400)

AT(I,J): Defined in REFS
 III: Defined in REFS
 T(K): Decay time in minutes of the K^{th} decay spectra of the J^{th} reference
 C(K): Counting time in minutes of the K^{th} decay spectra of the J^{th} reference
 CT(I,K): Counts registered at channel I of the K^{th} decay spectra for J^{th} reference isotope
 NO: Defined in REFS
 M: Tag name of the reference isotope
 DECAY: Defined in REFS
 COUNT: Defined in REFS
 ISO(J): Defined in REFS
 RHO(I): Defined in WATE
 A(I): Defined in WATE

10. XFORM:

This program is used to form the (X^0WX) matrix, its inverse, and the matrix correlation coefficients. A provision is included to evaluate the determinant of the matrix of correlation coefficients which is used as a measure of the efficiency for extracting the maximum information from the analysis problem. Further, the calculated identity matrix $(X^0WX)^{-1}(X^0WX)$ is formed and the error of average element can be found. Finally, the accuracy of the inverse matrix can be improved by an iteration method.

```
SUBROUTINE XFORM (NFORM,NISO,R,RI,MAXCON,P,O,LOCHN,NHICHN,A,AT,NCORR,NWRIT,
NIDEN,NINTF,RII,BETA,THETA,OI,RHO,III)
DIMENSION R(15,15),RI(15,15),RII(15,15),P(15),O(15),A(400),AT(400),OI(15),
Q(15,15),S(15,15),BETA(15),THETA(15),RHO(400)
```

NFORM: Specification for the $(X^0WX)^{-1}$ matrix to be read or formed

NISO: Defined in REFS

R(I,J): The inverse cross product matrix $(X^0WX)^{-1}$

RI(I,J): The cross product matrix (X^0WX)

MAXCON: The remaining number of reference isotopes

P(I): Tag name of the Ith reference isotope

O(I): Defined in REFS

LOCHN: Defined in INIT

NHICHN: Defined in INIT

A(I): Defined in WATE

AT(I): Defined in REFS

NCORR: Control variable which specifies formation of correlation coefficients

NWRIT: Control variable which specifies writing options for the (X^0WX) and its inverse

NIDEN: Control variable which specifies formation of identity matrix

NINTF: Control variable for interference calculation

RII(I,J): R(I,J), the $(X^0WX)^{-1}$ matrix

BETA(J): Concentration estimators (ppm) of the Jth reference isotope

THETA(J): Interferences in percent

OI(J): Intermediate vector

RHO(I): Defined in WATE

III: Defined in REFS

11. PEAKS:

This program is used to select the data near each photo-peak which is to be used in the area-of-interest unfolding procedure.

SUBROUTINE PEAKS (AT,A,RHO,LOCHN,NHICHN,NISO,III)

AT(I,J): Defined in REFS

A(I): Defined in WATE

RHO(I): Defined in WATE

LOCHN: Defined in INIT

NHICHN: Defined in INIT

NISO: Defined in REFS

III: Defined in REFS

12. SAVE:

This program is used to store the original $(X^0WX)^{-1}$ matrix and the corresponding reference spectra identification tags in order to allow stacking of the complex spectra

SUBROUTINE SAVE (NISO,O,OI,RI,RII)

NISO: Defined in REFS

O(J): The tag name of the reference spectrum

OI(J): O(J)

RI(I,J): The $(X^0WX)^{-1}$ matrix

RII(I,J): RI(I,J)

13. DET:

This program is used to evaluate the determinant by pivotal condensation (Chio's method) of the matrix composed of the correlation coefficients of the $(X^0WX)^{-1}$ matrix. The method is discussed in the appendix E.

SUBROUTINE DET (R,N,DETERM)

R: The matrix to be used $(X^0WX)^{-1}$

N: Order of the matrix

DETERM: The determinant of the matrix R.

14. DELETE:

This program is used to delete from the analysis the reference spectra which are not significant at a pre-selected level

SUBROUTINE DELETE (NISO,REMOVE,ISOCHN,BETA,R,C,O,THETA,L)

NISO: Defined in REFS

REMOVE: Intermediate value for the reference isotope

ISOCHN: Defined in INIT

BETA(I): New concentration estimator for Ith reference isotope

R(I,J): The $(X^0WX)^{-1}$ matrix

C(J): Defined in REFIN

O(J): Defined in REFS

THETA(J): Defined in XFORM

L: Tag name of the deleted isotope

15. STAT:

This program is used to calculate the standard deviations of the concentration estimators. In addition if some average weighting matrix is used in the analysis, this program is used to correct the standard deviations for this calculations.

SUBROUTINE STAT (REMOVE,BETA,RHOSQD,NISO,LOCHN,NHICHN,R,C,AT,A,RES,RHO,NVAR, NR,RAT,ST,NITER,F,RCORR,P,MAXCON,RHOSUM,XNORM,MVAR,ISO,III)

REMOVE: Defined in DELETE

BETA: Defined in DELETE

RHOSQD: $\sum (RHO(I))^2$

NISO: Defined in REFS

LOCHN: Defined in INIT

NHICHN: Defined in INIT

R(I,J): Defined in DELETE

O(J): Defined in REFS

AT(J): Defined in REFS

A(I): Defined in WATE

RES(I): Defined in PHOTO

RHO(I): Defined in WATE
 NVAR: Defined in INIT
 NR: Defined in INIT
 RAT: Defined in INIT
 ST: T -the student's t value
 NITER: Number of iterations to be performed
 F: F value for null hypothesis
 RCORR: Correlation coefficient
 P: Tag name for the reference isotope
 MAXCON: Remaining number of isotopes after deletions
 RHOSUM: $\sum (RHO(I))$
 XNORM: Defined in REFS
 MVAR: Defined in INIT
 ISO: Defined in REFS
 III: Defined in REFS

16. WRIT:

This program is used to print the concentration estimators, standard deviations, t values, etc.

SUBROUTINE WRIT (NITER, ITER, O, BETA, STDEV, ST, NISO, F, RCORR)

NITER: Iteration number
 ITER: Maximum iterations
 O(I): Tag number of the reference isotope
 BETA: Concentration estimators (ppm)
 STDEV: Standard deviations of BETA
 ST: Calculated t value
 NISO: Order of the isotope whose values are going to be printed
 F: Calculated F value
 RCORR: Correlation coefficients

This program is composed only of WRITE and FORMAT statements.

17. CHI:

This program is used to calculate the value of Chi-squared divided by the degrees of freedom for each analysis

78
SUBROUTINE CHI (CALC,RES,MAXCON,P,O,LOCHN,NHICHN,BETA,AT,C,RHO,NISO,QUAL,A)

CALC: Intermediate value for QUAL
RES(I): Defined in PHOTO
MAXCON: Defined in STAT
P(J): Defined in STAT
O(J): Defined in REFS
LOCHN: Defined in INIT
NHICHN: Defined in INIT
BETA(J): Defined in DELETE
AT(I): Defined in REFS
C(I): Defined in REFIN
RHO(I): Defined in WATE
NISO: Defined in REFS
QUAL: Quality of fit number
A(I): Defined in WATE
UP: Confidence interval of Chi-squared
DF: Degrees of freedom

18. VARNCE:

This program is used to adjust the weighting matrix to account for the variance in the reference spectra and the concentration estimators

SUBROUTINE VARNCE (A,OMG,STDEV,AT,LOCHN,NHICHN,NISO,MAXCON,C,P,O,BETA,XNORM)

All of the variables are previously defined.

19. WRIT1:

This program is used to write the significant concentrations estimators and their standard deviations, onto an intermediate storage.

SUBROUTINE WRIT1 (O,BETA,STDEV,NISO)

All of the variables are previously defined.

APPENDIX G

- LISTING OF THE INPUT DATA AND INPUT OPTIONS

Listing of input data and input options.

The following cards and decks are read by CORGAM:

A.1 Card One (15A4 format): Identification card for the set of data which may contain more than one complex spectrum.

A.2 Card Two (25I3 format):

NISO = number of isotopes used in analyzing all complex spectra of the data set. Any analysis for which a specific isotope is not desired may be removed by entering a negative unity in the proper location of the irradiation and decay times correction factor card.

NCASES = number of complex spectra contained in the data set.

Six other control parameters are carefully defined in the comment cards of the program.

A.3 Card Three (25I3 format):

NCHG1 = 0; no change in NAM1. See the comment cards in CORGAM for the definition of the variables in NAM1

NCHG1 = 1; change in NAM1. Below is a list of the standard values of the variables of NAM1. Read one or more cards with the necessary input data

NCHG2 = 0; no change in NAM2. This namelist contains the irradiation and decay time correction factors

NCHG2 = 1; change in NAM2. The standard value for all of these correction factors is unity. Read one or more cards with the necessary input data

NCHG3 = 0; no change in NAM3. See the comment cards in CORGAM for the definition of the variables in NAM3.

NCHG3 = 1; change in NAM3. This is the standard value for NCHG3 because not all of the necessary input data is defined. Below is a list of the standard values for the variables of NAM3. Read one or more cards with the necessary input data

A.4. Card Four (15A4 format): Identification card of a specific complex spectrum from the group which is to be analyzed.

A.5 Deck One (5E14.8 or 6X,12F6.0,2X format): This is the complex spectrum.

Note: Card Five and Deck Two are not needed if NBKGD = 0 or NUBKG = 0

A.6. Card Five (15A4 format): Identification card for the background spectrum.

A.7 Deck Two (6X,12F6.0,2X format): This is the background spectrum.

Note: Card Six and Deck Three are not needed unless NAREA = 1

A.8 Card Six (15A4 format): Identification card for the flux reference spectrum

A.9 Deck Three (6X,12F6.0,2X format): This is the flux reference spectrum

Note: Card Seven and Cards Eight are not needed if NSHFT = 0

A.10 Card Seven (8F10.0 format): This card contains the gain (MeV/ch) to which the complex spectrum is to be shifted

A.11 Cards Eight (4(I3,7X,F10.0) format): These cards contain the approximate mid point and energy of each photopeak used in shifting the complex spectrum

Note: If NVAR = 1, Deck Four is needed

A.12 Deck Four (5E14.8 format): This deck contains a variance spectrum corresponding to the complex spectrum.

The following cards and decks are read by subprogram REFS.

A disk is needed for intermediate storage.

Note: There will be NISO decks consisting of Card One and Card two plus a deck of control parameters and a reference spectrum

B.1 Card One (A4,9I3 format): This card contains an identification tag and several control parameters for a reference spectrum. See subprogram REFS comment cards for definitions

B.2 Card Two (8F10.0 format): This card contains the half life (in minutes), the irradiation time (in minutes), flux reference isotope half life (in minutes), and the reference spectrum normalization constant. This normalization constant is read as unity if the reference spectrum is to be normalized.

Note: Card Three is needed only if NO, the number spectra composing a set of reference decay spectra, is > unity.

B.3 Card Three (8F10.0 format): This card contains the decay time (in minutes) and counting time (in minutes) for each spectrum in the set of reference decay spectra. There is a maximum of five spectra allowed in each set of reference decay spectra; all the sets of reference decay spectra must come first in reference spectra decks; and each spectrum in a set of reference decay spectra must be in order of the decay times.

Note: Cards Four, Five, Deck One, Card Six, Deck Two, are not needed if NO is zero or unity

B.4 Card Four (25I3 format): This card contains several control parameters for the set of decay reference spectra.

- B.5 Card Five (8F10.0 format): This card contains the background counting time (in minutes), the flux reference decay time (in minutes), the sample weight (in grams), the weight of the flux reference (in grams), and the flux reference counting time (in minutes).
- B.6 Deck One (5E14.8;6X,12F6.0,2X; or 8F10.0 format): This deck contains the set of decay reference spectra.
 Note: If NORF is zero, Card Six, Deck Two, Card Seven, Card Eight, and Card Nine are not needed. If NNBKG is yew, Card Six and Deck Two are not needed.
- B.7 Card Six (15A4 format): Identification card for the background spectrum
- B.8 Deck Two (6X,12F6.0, 2X format): This deck contains the background spectrum
 Note: There will NO decks composed of Card Seve, Card Eight, and Cards Nine
- B.9 Card Seven (25I3 format): This card contains the number of peaks and the pivotal peak used in shifting the reference spectrum.
- B.10 Card Eight: See Card Seven of CORGAM
- B.11 Card Nine: See Card Eight of CORGRAM
 Note: If NO is zero, some or all of the following cards or decks are needed. If NORF is unity, Card Ten and Card Eleven are needed.
- B.12 Card Ten: See Card Four of REFS
- B.13 Card Eleven (8F10.0 format): This card contains the counting time (in minutes) and the information contained on Card Five of REFS for each reference spectrum
- B.14 Deck Three (5E14.8;6X,12F6.0,2X; or 8F10.0 format): This deck contains the reference spectrum.
 Note: If NORF is not equal to unity, Card Twelve, Deck Four, Cards Fourteen, Card Fifteen and Deck Five are not needed. If NUBKG is zero, Card Twelve and Deck Four are not needed.
- B.15 Card Twelve (15A4 format): Identification card for the background spectrum.
- B.16 Deck Four (6X,12F6.0,2X format): This deck contains a background spectrum
- B.17 Card Thirteen: See Card Seven of CORGAM.
- B.18 Card Fourteen: See Card Eight of CORGAM.
- B.19 Card Fifteen (15A4 format): Identification card for the flux reference spectrum
- B.20 Deck Five (6X,12F6.0,2X format): This deck contains a flux reference spectrum

The following cards and decks are read by subprogram XFORM

Note: If NFORM = 1, the following deck is needed:





C.1 Deck One (3E24.18 format): This deck contains the $(X'W X)^{-1}$ matrix which has been formed previously




Note: If NFORM = 2, no other input is required

If NFORM = 3, the following two cards are needed.

C.2 Card One (25I3 format): This card contains the number of peaks needed for the area of interest (AOI) unfolding procedure and the number of channels on either side of the mid-point channel of the photopeak used in the AOI unfolding procedure.

C.3 Card Two (25I3 format): This card contains the channel numbers which correspond to each peak used in the AOI unfolding procedure.





APPENDIX H

- LISTING OF THE UNIVAC-1106 VERSION OF THE CODE (CORGAM)

```

MAIN
DIMENSION RHO(400), AT(400,15), A(400), THETA(15), BETA(15), P(15), R(15), CORG
115), C(15), O(15), RES(400), CALC(400), RI(15,15), H(15), BKG(400), Q(15), ICORG
115), D(15), ISO(15), S(15,15), S1(400), OMG(400), STDEV(400), ST(15), NFAC(CORG
115), XNORM(15), OI(15), RII(15,15)
NAMELIST /NAM1/ISO,CHN,MODE,NBKG,NUBKG,NSHFT,NPK,LPP,NORM,NCTIN,
1 NAREA,NCHCU,NCULL,NCUUL,NRHO,NUMVAR,NVAR,NREFS,ITER,NR,NW,LOCHN,
2 NHICHN,NCHAN,MC,MVAR/NAM2/C/NAM3/T,DECAY,COUNT,WT,UCT,WTCU,AREA,
3 CUDEC,BKGTIM,TIMIRR,FRHL,UP,XNCOM
DATA ITER, RAT, SAV, C/0,1.0,0.0,15*1.0/
DATA ISO,CHN,MODE,NBKG,NUBKG,NSHFT,NPK,LPP,NORM,NCTIN,
1 NAREA,NCHCU,NCULL,NCUUL,NRHO,NUMVAR,NVAR,NREFS,ITER,NR,
2 NW,LOCHN,NHICHN,NCHAN,MC,MVAR/256,2,1,1,1,3,1,1,1,1,48,33,40,1
3,1,3,0,0,0,1,7,250,256,1,0/
DATA T, AREA, BKGTIM, TIMIRR, FRHL, UP, XNCOM/1.645,1.0,60.0,15.0,
1774.0,2.5758,1.0/
9 FORMAT(1H1)
10 FORMAT(5I3)
11 FORMAT(5E14.8)
12 FORMAT(8F10.0)
13 FORMAT(/24HL (ATRANS = * OMEGA * RHO) /)
14 FORMAT(11HT VALUE = ,F7.4)
15 FORMAT(15A4)
16 FORMAT(/35HL DECAY AND IRR CORRECTION FACTORS /)
17 FORMAT(/23HL SOLUTION VECTOR (PPM) /)
18 FORMAT(14,(1X,10E12.5))
19 FORMAT(2(1X,13,1X,E12.5,1X,E12.5,1X,E12.5))
20 FORMAT(1HI,2(1X,3,4CHN,5X,3HEXP,10X,4HCALC,10X,3HRES,4X))
21 FORMAT(1X,A4,(1X,9E12.5))
22 FORMAT(A4,9I3)
23 FORMAT(14,6(2X,A4,3H = ,E12.5))
24 FORMAT(3E24.18)
25 FORMAT(6X,12F6.0,2X)
26 FORMAT(35H NO VARIABLES SIGNIFICANT AT A T = ,F5.1)
27 FORMAT(/15H CU REF AREA = ,E14.8/)
28 FORMAT(/18H ABS DEV PER PT = ,E14.7,13H NO OF PTS = ,I3 /)
29 FORMAT(/6X,9(4X,A4,4X))
30 FORMAT(/17H SAMPLE WT(GM) = ,E12.5,17H IRR TIME(MIN) = ,E12.5,19H D
1 ECAY TIME(MIN) = ,E12.5,14H COUNT TIME = ,E12.5/19H FLUX REF WT(GM)
2 ) = ,F12.5,23H REF DECAY TIME(MIN) = ,E12.5,23H REF COUNT TIME(MIN)
3 ) = ,F12.5,18H REF HLIFE(MIN) = ,E11.4/17H BKG TIME(MIN) = ,E12.5)
50 FORMAT(/19H Y NORM CONSTANT = ,E15.8/)
51 FORMAT(8F10.3)
98 READ(5,15)(P(I),I=1,15)
READ(5,10) NISO,NCASES,NRES,NFORM,NWRIT,NCORR,NINTF,NIDEN
NRES=0 NO RES CALC
=1 RES CALC BUT NOT WRITTEN OUT,JUST UND OUTSIDE ABS(2)
=2 RES CALC AND WRITTEN OUT
NFORM = 1 READ INV COR MATRIX
= 2 FORM INV COR MATRIX
= 3 TAKE ONLY PEAK CHAN VALUES FOR LST SQRS FIT
NWRIT = 0 WRITE R AND RINV MATRICES
= 1 OTHER DO NOT WRITE R AND RINV MATRICES
NCORR = 0 FORM CORRELATION COEFF
= 1 OTHER DO NOT FORM CORRELATION COEFF
NINTF = 0 INTERFERENCES WILL BE CALCULATED
= 1 OTHER INTERFERENCES WILL NOT BE CALCULATED
NIDEN = 0 CALC THE IDENTITY MATRIX FROM THE INVERSE
= 1 OTHER DO NOT CALC THE IDENT
IF(NISO.EQ.0)CALL EXIT
NHI=NHICHN
WRITE(6,9)
WRITE(6,15)(P(I),I=1,15)
ISO,CHN = NO OF CHN IN COMPLEX SPECTRUM
MODE = 1 READ BY 5E14.8
= 2 READ BY 6X,12F6.0,2X
MAXCON=NISO
NBKGD = 0 NO BKGD SUBTRACTION
= 1 OTHER WILL SUBTRACT BKGD
NUBKG = 0 NO NEW BKGD TO BE READ
= 1 OTHER A NEW BKGD WILL BE READ
NSHFT = 0 NO SPECTRUM SHIFTING
= 1 OTHER WILL SHIFT SPECTRUM
NPK = NO OF PEAKS IN SPECT FOR SHIFTING
LPP = INDEX OF PK USED IN SHIFT RATIO
NORM = 0 NO NORMALIZATION OF COMPLEX NEEDED
= 1 OTHER NORMALIZE THE COMPLEX
NCTIN = 0 READ PRECALC CORRECTION FACTORS
= 1 OTHER READ UNITY (OR LT ZERO FOR DELETION)
NAREA = 0 READ SP AREA FOR CU REF
= 1 UNITY CALC PP AREA FOR CU REF,READ CU SPECT(READ ANY)
= 2 USE SAME AREA AS PREVIOUS SPECTRUM

```

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65
66
67
68
69
70
71
72
73
74
75
76
77
78
79
80
81
82
83
84
85
86
87
88
89
90
91
92
93
94
95
96
97
98
99
100

```

NCHCU = NO OF CHAN IN CU REF SPECT
NCULL = LOWER FIT LIMIT ON CU REF DO
NCUJL = UPPER FIT LIMIT ON CU REF DO
NRHO = 0 WILL NOT PUNCH SHIFTED OUTPUT
      = OTHER WILL PUNCH SHIFTED OUTPUT
NUMVAR = NO OF VAR TO AVERAGED
NVAR = 1 READ VAR BY 5E14.8
      = 2 VAR = 1.0
      = 3 VAR = RHO
      = 4 USE VAR PREVIOUSLY READ
NREFS = 0 NEW SET OF REF SPECT ARE READ
      = 1 NO NEW REF SPECT ARE TO BE READ
ITER = NO OF ITERATIONS ON VARIANCE
      = 0 MUST = 0 IF NVAR = 1,2,4
NR = 0 CAL AVE WT VAR THE LONG WAY
      = OTHER RATIO = YAVE/YAVEAVE
NW = 0 WILL NOT PUNCH SHIFTED OPT FOR REF SPECT
      = OTHER WILL PUNCH SHIFTED OPT FOR REF SPECT
LOCHN = LOWER CHN OF LST SORS FIT
NHICHN = UPPER CHN OF LST SORS FIT
NCHAN = NO OF CHN IN VAR SPECT IF VAR ARE TO BE READ
MC = NO OF REF SPECT WHICH HAVE A SET OF DECAY SPECTRA. THESE
      = SPECTRA MUST BE IN ORDER AND THEY MUST BE THE FIRST OF THE
      = SPECTRA TO BE READ IN
MVAR = 0 VAR OF REF SPECT = REF SPECT
      = OTHER READ VAR OF REF SPECT
C(I) = REAL CORRECTION FACTORS (NCTIN = 0)
      = ZERO OR LESS, USED FOR DELETION, (NCTIN GT 0)
T = STUDENTS T VALUE FOR NULL HYPOTHESIS TEST
      = BASE LINE SHIFT (NEG SHIFTS ZERO LEFT)
E = ENERGY LINE SHIFT IN SHIFTING
DECAY = DECAY TIME (MIN) FOR COMPLEX
COUNT = COUNT TIME (MIN) FOR COMPLEX
AREA = CU - EF PD NORMALIZED AREA
WT = COMPLEX WT (GM)
WTCU = WT OF CU REF (GM)
CUCT = CU REF COUNT TIME (MIN)
IF NAREA = 0, WTCU AND CUCT MUST EQUAL TO 1.0
CUDEC = DECAY TIME OF CU REF (MIN)
BKGTIM = COUNT TIME OF BKGD (MIN)
TIMIRR = COMPLEX SAMPLE IRR TIME
FRHL = FLUX REF HALF LIFE (MIN)
UP = U-SUB-P FACTOR CHISQR TEST
XNCOM = COMPLEX SPECTRUM NORM CONST
DO 1002 III=1, NCASES
IF (III.LE.1) GO TO 3011
CALL INIT (ITER, RAT, SAV, C, ISOCHN, MODE, NBKGD, NUBKG, NSHFT, NDK, LPP, NO
1 2, NCTIN, NAREA, NCHCU, NCULL, NCUJL, NRHO, NUMVAR, NVAR, NREFS, ITER, NR, NW
2, LOCHN, NHICHN, NCHAN, MC, MVAR, T, AREA, BKGTIM, TIMIRR, FRHL, UP, XNCOM)
11 READ (5, 10) NCHG1, NCHG2, NCHG3
IF (NCHG1.EQ.0) GO TO 52
READ (5, NAM1)
52 IF (NCHG2.EQ.0) GO TO 53
READ (5, NAM2)
53 IF (NCHG3.EQ.0) GO TO 54
READ (5, NAM3)
54 READ (5, 15) (P(I), I=1, 15)
DO = P(1)
GO TO (111, 134), MODE
11 READ (5, 11) (RHO(I), I=1, ISOCHN)
GO TO 135
34 READ (5, 34) (RHO(I), I=1, ISOCHN)
35 WRITE (6, 15) (P(I), I=1, 15)
WRITE (6, 46) WT, TIMIRR, DECAY, COUNT, WTCU, CUDEC, CUCT, FRHL, BKGTIM
IF (NRHO.EQ.0) GO TO 3333
WRITE (2, 15) (P(I), I=1, 15)
33 IF (NBKGD.EQ.0) GO TO 3170
FACT = COUNT/BKGTIM
IF (NUBKG.EQ.0) GO TO 374
READ (5, 15) (P(I), I=1, 15)
74 READ (5, 34) (BKG(I), I=1, ISOCHN)
DO 318 I=1, ISOCHN
RHO(I) = RHO(I) - BKG(I) * FACT
IF (RHO(I).LT.0.) RHO(I) = 0.0
18 CONTINUE
70 WRITE (1, 0)
WRITE (TIMIRR.EQ.15) GO TO 317
WRITE (1 - EXP(-ALOG(2.) * TIMIRR / FRHL)) / (1 - EXP(-ALOG(2.) * 15. / FRHL))
17 IF (NORM.EQ.0) GO TO 133
IF (NAREA.EQ.2) GO TO 3171
IF (NAREA.EQ.0) GO TO 370
READ (5, 15) (P(I), I=1, 15)

```

```

CORGG 78
CORGG 79
CORGG 80
CORGG 81
CORGG 82
CORGG 83
CORGG 84
CORGG 85
CORGG 86
CORGG 87
CORGG 88
CORGG 89
CORGG 90
CORGG 91
CORGG 92
CORGG 97
CORGG 98
CORGG 99
CORGG 100
CORGG 101
CORGG 102
CORGG 103
CORGG 104
CORGG 105
CORGG 106
CORGG 107
CORGG 108
CORGG 109
CORGG 110
CORGG 111
CORGG 112
CORGG 113
CORGG 114
CORGG 115
CORGG 116
CORGG 117
CORGG 118
CORGG 119
CORGG 120
CORGG 121
CORGG 122
CORGG 123
CORGG 124
CORGG 125
CORGG 126
CORGG 127
CORGG 128
CORGG 129
CORGG 130
CORGG 131
CORGG 132
CORGG 133
CORGG 134
CORGG 135
CORGG 136
CORGG 137
CORGG 138
CORGG 139
CORGG 140
CORGG 141
CORGG 142
CORGG 143
CORGG 144
CORGG 145
CORGG 146
CORGG 147
CORGG 148
CORGG 149
CORGG 150
CORGG 151
CORGG 152
CORGG 153
CORGG 154
CORGG 155
CORGG 156
CORGG 157
CORGG 158
CORGG 159
CORGG 160
CORGG 161
CORGG 162
CORGG 163

```

```

READ(5,34)(RES(I),I=1,NCHCN)
CALL PHOTO(RES,XO,ICOLL,NCUUL,0.511,AREA)
CUDEC=CUDEC+CUCT/2.
AREA=AREA*EXP(CJDEC*ALOG(2.)/FRHL)/(WTCU*CUCT)
570 WRITE(6,36)AREA
171 XNCOM=AREA*WT*1.0E=06*COUNT*SF
DO 131 I=1,ISOCHN
131 RHO(I)=RHO(I)/XNCOM
WRITE(6,50)XNCOM
133 IF(NSHFT.EQ.0) GO TO 3013
ICHN=ISOCHN
CALL SHIFT(ISOCHN,NPK,LPP,NRHO,RHO,PO)
RATIO=(FLOAT(ISOCHN)-0.5)/FLOAT(ICHN)
NHICHN=FLOAT(NHICHN)*RATIO+0.5
LOCHN=FLOAT(LOCHN)*RATIO+0.5
IF(NVAR.EQ.4) GO TO 3190
013 CALL WATE(NVAR,NUMVAR,NCHAN,RHO,XNCOM,III,OMG,A,RES,ISOCHN)
IF(DECAY.EQ.SAV) GO TO 3190
IF(III.GT.1) GO TO 3195
SAV=DECAY
195 IF(NREFS.EQ.1) GO TO 3190
CALL REFS(NISO,NHICHN,ISO,XNORM,O,AT,RKG,NW,DECAY,COUNT,III,NFAC,
1 FRRHL,D,H,MC,NO)
190 DO 978 I=1,MAXCON
IF(NCTIN.EQ.0) GO TO 978
XN=(1.-EXP(-ALOG(2.)*15./FRRHL))
SF=XN/(1.-EXP(-ALOG(2.)*D(I)/FRRHL))
XL=ALOG(2.)/2H(I)
IF(NFAC(I).EQ.0) GO TO 3337
FACT=1.0
GO TO 130
337 TB=DECAY*ALOG((1.-EXP(-XL*COUNT))/(COUNT*XL))/XL
FACT=EXP(XL*TB)
130 C(I)=FACT*C(I)*1.0E+06*SF*(1.-EXP(-XL*D(I)))/(1.-EXP(-XL*15.))
978 P(I)=C(I)
WRITE(6,16)
WRITE(6,27)NISO,(C(I),I=1,NISO)
IF(NVAR.EQ.4) GO TO 319
99 CALL XFORM(NFORM,NISO,R,RI,MAXCON,P,O,LOCHN,NHICHN,A,AT,NCORR,NWRIC
1 I,NIDFN,NINTF,RII,BETA,THETA,OI,RHO,III)
119 RHOSUM=0.
132 RHOSQD=0.
REMOVE=0.
DO 120 I=LOCHN,NHICHN
IF(A(I).EQ.0.0) A(I)=1.0
RHOSQD=RHOSQD+(RHO(I)**2)/A(I)
120 RHOSUM=RHOSUM+RHO(I)/A(I)
266 RAT=RHOSUM/FLOAT(NHICHN-LOCHN+1)
FORM ATRANS * OMEGA * RHO
253 DO 117 K=1,NISO
117 THETA(K)=0.0
L=1
DO 4050 K=1,MAXCON
IF(P(K).NE.0(L)) GO TO 4050
DO 1000 I=LOCHN,NHICHN
000 THETA(L)=THETA(L)+AT(I,K)*RHO(I)/A(I)
L=L+1
IF(L.GT.NISO) GO TO 4051
550 CONTINUE
551 WRITE(6,13)
WRITE(6,18)NISO,(THETA(I),I=1,NISO)
WRITE(6,9)
701 DO 106 J=1,NISO
BETA(J)=0.0
DO 106 I=1,NISO
106 BETA(J)=BETA(J)+R(J,I)*THETA(I)
DO 125 I=1,NISO
25 BETA(I)=BETA(I)*C(I)
DO 500 I=1,NISO
000 REMOVE=REMOVE+BETA(I)*THETA(I)/C(I)
ISOCHN=NISO
L=1
02 IF(C(I)) 601,600,600
01 CALL DELETE(NISO,REMOVE,ISOCHN,BETA,R,C,O,THETA,L)
L=L-1
00 L=L+1
IF(L.LE.NISO) GO TO 602
CALL STAT(REMOVE,BETA,RHOSQD,NISO,LOCHN,NHICHN,R,C,O,AT,A,STDEV,RH
1 O,NVAR,NR,RAT,ST,NITER,F,RCORR,P,MAXCON,RHOSUM,XNORM,MVAR,ISO,III)
CALL WRIT(NITER,ITER,O,BETA,STDEV,ST,NISO,F,RCORR)
52 IF(ITER.EQ.0) GO TO 262
IF(NITER.GE.ITER) GO TO 262
CORRG 164
CORRG 165
CORRG 166
CORRG 167
CORRG 168
CORRG 169
CORRG 170
CORRG 171
CORRG 172
CORRG 173
CORRG 174
CORRG 175
CORRG 176
CORRG 177
CORRG 178
CORRG 179
CORRG 180
CORRG 181
CORRG 182
CORRG 183
CORRG 184
CORRG 185
CORRG 186
CORRG 187
CORRG 188
CORRG 189
CORRG 190
CORRG 191
CORRG 192
CORRG 193
CORRG 194
CORRG 195
CORRG 196
CORRG 197
CORRG 198
CORRG 199
CORRG 200
CORRG 201
CORRG 202
CORRG 203
CORRG 204
CORRG 205
CORRG 206
CORRG 207
CORRG 208
CORRG 209
CORRG 210
CORRG 211
CORRG 212
CORRG 213
CORRG 214
CORRG 215
CORRG 216
CORRG 217
CORRG 218
CORRG 219
CORRG 220
CORRG 221
CORRG 222
CORRG 223
CORRG 224
CORRG 225
CORRG 226
CORRG 227
CORRG 228
CORRG 229
CORRG 230
CORRG 231
CORRG 232
CORRG 233
CORRG 234
CORRG 235
CORRG 236
CORRG 237
CORRG 238
CORRG 239
CORRG 240
CORRG 241
CORRG 242
CORRG 243
CORRG 244
CORRG 245

```



```

CALL CHI(CALC,RES,MAXCON,P,0,LOCHN,NHICHN,BETA,AT,C,RHO,NISO,CHIS,CORR) 246
1A) NITER=NITER+1 247
NCTIN=0 248
JPMX=IP*SQRT(2./FLOAT(NHICHN-LOCHN-NISO))+1. 249
IF(CHIS.LE.JPMX)GO TO 252 250
CALL VARNCE(A,OMG,STDEV,AT,LOCHN,NHICHN,NISO,MAXCON,C,P,0,BETA,XNO 251
1RM) 252
NFORM=2 253
GO TO 99 254
262 ISOCHN=NISO 255
NITER=0 256
AMIN=ST(1) 257
DO 950 I=1,NISO 258
IF(AMIN.LT.ST(I))GO TO 950 259
L=I 260
AMIN=ST(I) 261
950 CONTINUE 262
IF(AMIN.GT.T)GO TO 951 263
IF(NISO.GT.1)GO TO 953 264
WRITE(6,35)T 265
GO TO 954 266
953 CALL DELETE(NISO,REMOVE,ISOCHN,BETA,R,C,0,THETA,L) 267
CALL STAT(REMOVE,BETA,RHOSQD,NISO,LOCHN,NHICHN,R,C,0,AT,A,STDEV,RH 268
10,NVAR,NR,RAT,ST,NITER,F,RCORR,P,MAXCON,RHOSUM,XNORM,MVAR,ISO,III) 269
CALL WRIT(NITER,ITER,0,BETA,STDEV,ST,NISO,F,RCORR) 270
IF(NISO.GT.0)GO TO 951 271
WRITE(6,35)T 272
NISO=MAXCON 273
CALL SAVE(NISO,0I,0,RII,R) 274
GO TO 1002 275
51 IE(NISO.LT.ISOCHN)GO TO 952 276
954 JES=NHICHN-LOCHN-NISO 277
IF(NRES.EQ.0)GO TO 959 278
CALL CHI(CALC,RES,MAXCON,P,0,LOCHN,NHICHN,BETA,AT,C,RHO,NISO,CHIS,CORR) 279
1A) REMOVE=0. 280
DO 971 J=LOCHN,NHICHN 281
REMOVE=REMOVE+RES(J)**2 282
REMOVE=SQRT(REMOVE) 283
NOP=NHICHN-LOCHN+1 284
TOTAL=0.0 285
DO 980 J=LOCHN,NHICHN 286
TOTAL=TOTAL+ABS(RES(J)) 287
TOTAL=TOTAL/FLOAT(NOP) 288
WRITE(6,38)TOTAL,NOP 289
IF(NRES.NE.2)GO TO 959 290
WRITE(6,21) 291
WRITE(6,20)(J,RHO(J),CALC(J),RES(J),J=LOCHN,NHICHN) 292
59 CALL WRIT1(0,BETA,STDEV,NISO) 293
NISO=MAXCON 294
WRITE(6,14)T 295
WRITE(6,9) 296
CALL SAVE(NISO,0I,0,RII,R) 297
NHICHN=NH 298
02 CONTINUE 299
GO TO 198 300
END 301
SUBROUTINE REFS(NISO,NHICHN,ISO,XNORM,0,AT,BKG,NW,DECAY,COUNT,III, 302
1NFAC,FRRHL,D,H,MC,NO,RHO,A) 303
DIMENSIONISO(15),H(15),D(15),XNORM(15),O(15),AT(400,15),BKG(400), 304
1RES(400),P(15),CALC(400),NFAC(15),CT(400,5),T(5),C(5),RHO(400),A(4 305
200) 306
10 FORMAT(25I3) 307
11 FORMAT(5E14.8) 308
12 FORMAT(8F10.0) 309
13 FORMAT(15A4) 310
14 FORMAT(I4,(1X,10E12.5)) 311
15 FORMAT(A4,9I3) 312
16 FORMAT(6X,12F5.0,2X) 313
17 FORMAT(/154CU REF AREA =,E14.8/) 314
18 FORMAT(/6X,9(4X,A4,4X)) 315
19 FORMAT(8F10.3) 316
REWIND 10 317
MM=1 318
DO 121 I=1,NISO 319
IF(III.GT.1)GO TO 3194 320
READ(5,24)O(I),ISO(I),NPK,LPP,NORF,MOE,NFAC(I),NO 321
READ(5,12)H(I),D(I),FRRHL,XNORM(I) 322
IF(NO.EQ.0)GO TO 91 323
READ(5,12)(T(J),C(J),J=1,NO) 324
ISOCHN=NO OF CHN IN REF 325

```

NPCK = TAG NAME PKF REF SPECT
 LPPCK = INDEX OF PIVOTAL PKF SPECT FOR SHIFTING
 NORF = 1 NORM OF REF REQ
 MODF = 1 USE 14, 8 INPUT FORMAT
 2 6X, 12F5, J, 2X INPUT FORMAT
 3 8F10.0 INPUT FORMAT
 NFAC = 0 REF IS AT TIME ZERO
 NFAC = 1 REF IS AT TIME OF COMPLEX
 3T 1 REF IS AT TIME OF DECAY CORRECTION
 NO = NO OF SPECT IN THE SET OF REF SPECT
 H = HALF LIFE OF REF SPECT
 D = IRRAD TIME OF REF SPECT (MIN)
 FRPHL = FLUX REF HALF LIFE FOR REF (MIN)
 T = DECAY TIME (MIN) OF REF SPECT
 C = COUNT TIME (MIN) OF REF SPECT

```

91 PO=0(I)
   ISOCHN=ISO(I)
14 DO 1114 J=1,400
   AT(J,I)=0.0
   IF(NORF.NE.1)GO TO 256
   READ(5,10)NJBKG,NCJLL,NCUUL,NCHCU,NCHBK
   READ(5,12)COUN,DECA,BKGTIM,CUDEC,WT,WTCU,CUCT
   DO 114 M=1,NO
   DO 115 J=1,ISOCHN
15 CT(J,M)=0.0
   GO TO(116,117,118),MODE
16 READ(5,11)(CT(J,M),J=1,ISOCHN)
   GO TO 114
17 READ(5,34)(CT(J,M),J=1,ISOCHN)
   GO TO 114
18 READ(5,12)(CT(J,M),J=1,ISOCHN)
14 CONTINUE
   FACT=0.0
   IF(NORF.EQ.0)GO TO 257
   IF(NURKG.EQ.0)GO TO 262
   READ(5,15)(P(J),J=1,15)
   READ(5,34)(BKG(J),J=1,NCHBK)
62 DO 259 K=1,NO
   FACT=C(K)/BKGTIM
   READ(5,10)NPCK,LPP
   DO 258 J=1,ISOCHN
   CALC(J)=CT(J,K)-FACT*BKG(J)
58 IF(CALC(J).LT.0.)CALC(J)=0.0
60 ICHN=ISOCHN
   CALL SHIFT(ICHN,NPK,LPP,0,CALC,PO)
   IF(ICHN.GT.NHICHN)ICHN=NHICHN
   DO 259 J=1,ICHN
59 CT(J,K)=CALC(J)
   ISOCHN=ICHN
   ISO(I)=ISOCHN
   GO TO 257
56 IF(NORF.NE.1)GO TO 93
   READ(5,10)NJBKG,NCJLL,NCUUL,NCHCU,NCHBK
   READ(5,12)COUN,DECA,BKGTIM,CUDEC,WT,WTCU,CUCT
93 GO TO(21,22,23),MODE
21 READ(5,11)(AT(J,I),J=1,ISOCHN)
   GO TO 224
22 READ(5,34)(AT(J,I),J=1,ISOCHN)
   GO TO 224
23 READ(5,12)(AT(J,I),J=1,ISOCHN)
24 IF(NORF.NE.1)GO TO 910
   FACT=0.0
   FACT=COUN/BKGTIM
   IF(NURKG.EQ.0)GO TO 264
   READ(5,15)(P(J),J=1,15)
   READ(5,34)(BKG(J),J=1,NCHBK)
   GO TO 264
64 DO 260 J=1,ISOCHN
   CALC(J)=AT(J,I)-FACT*BKG(J)
60 IF(CALC(J).LT.0.)CALC(J)=0.0
   CALL SHIFT(ISOCHN,NPK,LPP,0,CALC,PO)
   IF(ISOCHN.GT.NHICHN)ISOCHN=NHICHN
   DO 261 J=1,ISOCHN
61 AT(J,I)=CALC(J)
   ISO(I)=ISOCHN
67 IF(NORF.NE.1)GO TO 910
94 CALL BEFIN(AT,III,T,C,CT,NO,I,DECA,COUNT,ISO,RHO,A)
   WRITE(6,39)O(I)
   IC=ISO(I)
   WRITE(6,18)ISO(I),(AT(J,I),J=1,IC)
   IF(III.LE.1)GO TO 910
  
```

2270
 2271
 2272
 2273
 2274
 2275
 2276
 2277
 2278
 2279
 2280
 2281
 2282
 2283
 2284
 2285
 2286
 2287
 2288
 2289
 2290
 2291
 2292
 2293
 2294
 2295
 2296
 2297
 2298
 2299
 2300
 2301
 2302
 2303
 2304
 2305
 2306
 2307
 2308
 2309
 2310
 2311
 2312
 2313
 2314
 2315
 2316
 2317
 2318
 2319
 2320
 2321
 2322
 2323
 2324
 2325
 2326
 2327
 2328
 2329
 2330
 2331
 2332
 2333
 2334
 2335
 2336
 2337
 2338
 2339
 2340
 2341
 2342
 2343
 2344
 2345
 2346
 2347
 2348
 2349
 2350
 2351
 2352
 2353
 2354
 2355
 2356
 2357
 2358
 2359
 2360
 2361
 2362
 2363
 2364
 2365
 2366
 2367
 2368
 2369
 2370
 2371
 2372
 2373
 2374
 2375
 2376
 2377
 2378
 2379
 2380
 2381
 2382
 2383
 2384
 2385
 2386
 2387
 2388
 2389
 2390
 2391
 2392
 2393
 2394
 2395
 2396
 2397
 2398
 2399
 2400
 2401
 2402
 2403
 2404
 2405
 2406
 2407
 2408
 2409
 2410
 2411
 2412
 2413
 2414
 2415
 2416
 2417
 2418
 2419
 2420
 2421
 2422
 2423
 2424
 2425
 2426
 2427
 2428
 2429
 2430
 2431
 2432
 2433
 2434
 2435
 2436
 2437
 2438
 2439
 2440
 2441
 2442
 2443
 2444
 2445
 2446
 2447
 2448
 2449
 2450
 2451
 2452
 2453
 2454
 2455
 2456
 2457
 2458
 2459
 2460
 2461
 2462
 2463
 2464
 2465
 2466
 2467
 2468
 2469
 2470
 2471
 2472
 2473
 2474
 2475
 2476
 2477
 2478
 2479
 2480
 2481
 2482
 2483
 2484
 2485
 2486
 2487
 2488
 2489
 2490
 2491
 2492
 2493
 2494
 2495
 2496
 2497
 2498
 2499
 2500

```

MM=MM+1
IF (MM.LE.MC) GO TO 121
REWIND 10
RETURN
910 IF (NORF.EQ.1) GO TO 9913
READ(5,15) (P(J), J=1,15)
READ(5,34) (RES(J), J=1,NCHCU)
CALL PHOTO(RES,XO,NCOLL,NCUUL,0.511,AREA)
IF (CUDEC.EQ.0.) GO TO 263
AREA=AREA*EXP(CUDEC*ALOG(2.)/FRRHL)
263 WRITE(6,36) AREA
FACT=1.0
IF (NFAC(I).GT.1) GO TO 4022
XL=ALOG(2.)/H(I)
TB=DECA=ALOG((1.-EXP(-XL*COUN))/(COUN*XL))/XL
FACT=EXP(XL*TB)
4022 XN=AREA*WT7(WTCU*CJCT*1.0E+06*FACT)
DO 3019 J=1,ISOCHN
XNORM(I)=XN*COUN
3019 AT(J,I)=AT(J,I)/(XNORM(I))
WRITE(6,39) O(I)
WRITE(6,18) ISOCHN, (AT(J,I), J=1,ISOCHN)
IF (NO.LE.1) GO TO 9913
MODE=1
DO 913 K=1,NO
DO 912 J=1,ISOCHN
912 CT(J,K)=CT(J,K)/XNORM(I)
IF (NW.EQ.0) GO TO 913
NORE=0
MODE=1
WRITE(2,24) O(I), ISO(I), NPK, LPP, NORF, MODE, NFAC(I), NO
WRITE(2,51) H(I), D(I), FRRHL, XNORM(I)
WRITE(2,11) (CT(J,K), J=1,ISOCHN)
9913 CONTINUE
9913 IF (III.GE.1) GO TO 122
IF (NORF.EQ.1) GO TO 911
IF (NO.LE.1) GO TO 121
911 NORF=0
MODE=1
WRITE(2,24) O(I), ISO(I), NPK, LPP, NORF, MODE, NFAC(I), NO
WRITE(2,51) H(I), D(I), FRRHL, XNORM(I)
WRITE(2,11) (AT(J,I), J=1,ISOCHN)
121 CONTINUE
122 WRITE(10) T, C, CT, NO, AT, ISO
REWIND 10
RETURN
END
SUBROUTINE REFIN(AT,III,T,C,CT,NO,M,DECAY,COUNT,ISO,RHO,A)
DIMENSION AT(400,15), CT(400,5), T(5), C(5), ISO(15), IS(15), RHO(400), A(
1400)
USFS THE LAGRANGE METHOD TO INTERPOLATE BETWEEN SPECTRA
ISOCHN=ISO(M)
IF (III.LE.1) GO TO 154
ISOCHN=IS(M)
154 READ(10) T, C, CT, NO, AT, ISO
X=DECAY+COUNT/2.
DO 150 I=1,NO
150 T(I)=T(I)+C(I)/2.
DO 151 J=1,ISOCHN
AT(I,M)=0.
DO 152 J=1,NO
CT(I,J)=CT(I,J)/C(J)
XNUM=1.0
XDEN=1.0
DO 153 L=1,NO
IF (L.EQ.J) GO TO 153
XNUM=XNUM*(X-T(L))
XDEN=XDEN*(T(J)-T(L))
153 CONTINUE
152 AT(I,M)=XNUM*CT(I,J)/XDEN+AT(I,M)
151 CONTINUE
IS(M)=ISOCHN
RETURN
END
SUBROUTINE PHOTO(RES,XO,IPPL,IPPUL,EDP,AREA)
DIMENSION B(3), Q(15,15), RES(400), RE(3), D2(50), C(10)
4 FORMAT(35H LOWER LIMIT OF PHOTOPEAK FROM CHAN. 13.84 TO CHAN. 13 / 3
10 15H UPPER LIMIT OF PHOTOPEAK FROM CHAN. 13.84 TO CHAN. 13)
9 FORMAT( // 30H EXPERIMENTAL PHOTOPEAK POINTS // )
10 FORMAT( // 28H CALCULATED PHOTOPEAK POINTS // )
15 FORMAT(1X,5E14,7)
16 FORMAT( // 2V, 13V, 24X, 0.13X, 2HRO, 15X, 2HYO, 13X, 4HAREA // )

```

| | |
|-----|-----|
| 107 | 107 |
| 108 | 108 |
| 109 | 109 |
| 110 | 110 |
| 111 | 111 |
| 112 | 112 |
| 113 | 113 |
| 114 | 114 |
| 115 | 115 |
| 116 | 116 |
| 117 | 117 |
| 118 | 118 |
| 119 | 119 |
| 120 | 120 |
| 121 | 121 |
| 122 | 122 |
| 123 | 123 |
| 124 | 124 |
| 125 | 125 |
| 126 | 126 |
| 127 | 127 |
| 128 | 128 |
| 129 | 129 |
| 130 | 130 |
| 131 | 131 |
| 132 | 132 |
| 133 | 133 |
| 134 | 134 |
| 135 | 135 |
| 136 | 136 |
| 137 | 137 |
| 138 | 138 |
| 139 | 139 |
| 140 | 140 |
| 141 | 141 |
| 142 | 142 |
| 143 | 143 |
| 144 | 144 |
| 145 | 145 |
| 146 | 146 |
| 147 | 147 |
| 148 | 148 |
| 149 | 149 |
| 150 | 150 |
| 151 | 151 |
| 152 | 152 |
| 153 | 153 |
| 154 | 154 |
| 155 | 155 |
| 156 | 156 |
| 157 | 157 |
| 158 | 158 |
| 159 | 159 |
| 160 | 160 |
| 161 | 161 |
| 162 | 162 |
| 163 | 163 |
| 164 | 164 |
| 165 | 165 |
| 166 | 166 |
| 167 | 167 |
| 168 | 168 |
| 169 | 169 |
| 170 | 170 |
| 171 | 171 |
| 172 | 172 |
| 173 | 173 |
| 174 | 174 |
| 175 | 175 |
| 176 | 176 |
| 177 | 177 |
| 178 | 178 |
| 179 | 179 |
| 180 | 180 |
| 181 | 181 |
| 182 | 182 |
| 183 | 183 |
| 184 | 184 |
| 185 | 185 |
| 186 | 186 |
| 187 | 187 |
| 188 | 188 |
| 189 | 189 |
| 190 | 190 |
| 191 | 191 |
| 192 | 192 |
| 193 | 193 |
| 194 | 194 |
| 195 | 195 |
| 196 | 196 |
| 197 | 197 |
| 198 | 198 |
| 199 | 199 |
| 200 | 200 |

```

17 FORMAT(F6.3,7X,F8.3,7X,F8.3,7X,E12.5,7X,E12.5 )
18 FORMAT(9X,4H+OR=,F8.3,3X,4H+OR=,F8.3,3X,4H+OR=,E12.5,3X,4H+OR=,E12.5)
19 (X)=YO*EXP(-((X-XO)**2)/BO)
1 DO 50 I=1,3
1 DO 50 J=1,3
1 DO 50 K=1,3
212 XN=0
DO 321 I=1,3
DO 321 J=1,3
321 DO 30 I=1,3
DO 30 J=1,3
XN2=XN*XN
XN3=XN2*XN
XN4=XN3*XN
ZN=ALOG(RES(I))
Q(1,1)=Q(1,1)+XN4*RES(I)
Q(1,2)=Q(1,2)+XN3*RES(I)
Q(1,3)=Q(1,3)+XN2*RES(I)
Q(2,3)=Q(2,3)+XN*RES(I)
Q(3,3)=Q(3,3)+RES(I)
B(1)=B(1)+XN2*ZN*RES(I)
B(2)=B(2)+XN*ZN*RES(I)
B(3)=B(3)+ZN*RES(I)
30 XN=XN+1
Q(2,1)=Q(1,2)
Q(3,1)=Q(1,3)
Q(2,2)=Q(1,3)
Q(3,2)=Q(2,3)
CALL SMINV(0,3)
DO 110 J=1,3
RE(J)=0
DO 110 I=1,3
RE(I)=RE(I)+Q(J,I)*R(I)
XO=(-RE(2)/(2.*RE(1)))
YO=1./RE(1)
YO=EXP(XO*YO/BO+RE(3))
IF(BO.LT.0.)BO=ABS(BO)
AREA=SQRT(3.1416*BO)*YO
XO=XO+FLOAT(IPPLL)
DO 80 I=1,3
B(I)=RE(I)
ZN=0
STDE=0
XN=0
DO 40 I=IPPLL,IPPUL
STDE=STDE+(ALOG(RES(I))-B(3)-B(2)*XN-B(1)*XN*XN)**2)
RESI=RES(I)+((RES(I)-F(FLOAT(I)))**2)
40 XN=XN+1
IF(RESI.LT.0.)RESI=ABS(RESI)
STDE=SQRT(RESI/FLOAT(IPPUL-IPPLL-2))
IF(1.5*STDE.GT.ABS(RES(IPPLL)-F(FLOAT(IPPLL))))GO TO 50
IPPLL=IPPLL+1
50 IF(1.5*STDE.GT.ABS(RES(IPPUL)-F(FLOAT(IPPUL))))GO TO 60
IPPUL=IPPUL-1
60 IF(NEQ.0.1)GO TO 212
IF(IPPULO.NE.1)GO TO 70
IF(IPPULO.EQ.1)GO TO 71
WRITE(6,4)IPPLO,IPPLL,IPPUL,IPPUL
70 ZN=1./((4.*3(1)**2)
XN=(B(2)/B(1))**2)
XN2=1./B(1)**4)
XN3=XN/4
THE FOLLOWING STATEMENT HAS BEEN CHANGED BY THE ENEA CPL SO THAT IT
CORRESPONDS TO THE FORMULA WHICH IS IN THE REPORT.
XN4=(XN**2)/FLOAT(NCHN)
C(1)=XN*(Q(2,2)+XN*Q(1,1))
C(2)=XN2*Q(1,1)
C(3)=(YO**2)*(Q(3,3)+XN3*Q(2,2)+XN4*Q(1,1))
C(4)=3.1416*(BO*C(3)+C(2)*YO**2/(2.*30))
DO 100 I=1,4
C(I)=SQRT(C(I))
100 WRITE(6,16)
WRITE(6,17)EPP,XO,BO,YO,AREA
WRITE(6,18)C(1),C(2),C(3),C(4)
RETURN
SUBROUTINE SHIFT(NCHN,NPK,LPP,NWRITE,Y,PO)

```

FOTO 9
FOTO 10
FOTO 11
FOTO 12
FOTO 13
FOTO 14
FOTO 15
FOTO 16
FOTO 17
FOTO 18
FOTO 19
FOTO 20
FOTO 21
FOTO 22
FOTO 23
FOTO 24
FOTO 25
FOTO 26
FOTO 27
FOTO 28
FOTO 29
FOTO 30
FOTO 31
FOTO 32
FOTO 33
FOTO 34
FOTO 35
FOTO 36
FOTO 37
FOTO 38
FOTO 39
FOTO 40
FOTO 41
FOTO 42
FOTO 43
FOTO 44
FOTO 45
FOTO 46
FOTO 47
FOTO 48
FOTO 49
FOTO 50
FOTO 51
FOTO 52
FOTO 53
FOTO 54
FOTO 55
FOTO 56
FOTO 57
FOTO 58
FOTO 59
FOTO 60
FOTO 61
FOTO 62
FOTO 63
FOTO 64
FOTO 65
FOTO 66
FOTO 67
FOTO 68
FOTO 69
FOTO 70
FOTO 71
FOTO 72
FOTO 73
FOTO 74
FOTO 75
FOTO 76
FOTO 77
FOTO 78
FOTO 79
FOTO 80
FOTO 81
FOTO 82
FOTO 83
FOTO 84
FOTO 85
FOTO 86
SHIFT 1


```

22 X=X+1.
   I1=I1+1.
   AD=(X/RATIO)+RLSHFT
   IF(XP.LE.1.)GO TO 24
   IF(XP.GE.2.)GO TO 25
   GO TO 22
24 S(I1)=SPECT(1)/RATIO
   GO TO 22
26 X1=1.
   X2=2.
   X3=3.
   Y1=SPECT(1)
   Y2=SPECT(2)
   Y3=SPECT(3)
   GO TO 27
25 IX2=XP+.5
   IF(IX2.GE.NCHN)GO TO 28
   IX1=IX2-1
   IX3=IX2+1
   Y1=SPECT(IX1)
   Y2=SPECT(IX2)
   Y3=SPECT(IX3)
   X1=IX1
   X2=IX2
   X3=IX3
27 AA=(Y3-2.*Y2+Y1)/2.
   BB=(-Y3*(X2+X1)+2.*Y2*(X1+X3)-Y1*(X3+X2))/2.
   CC=((X1*X2*Y3)-2.*(X1*X3*Y2)+(X2*X3*Y1))/2.
28 S(I1)=(AA*XP+BB*XP+CC)/RATIO
   IF(S(I1).LT.0.)S(I1)=0.0
   IF(I1.LT.NCHN2)GO TO 22
   RETURN
END
SUBROUTINE WRIT1(O,BETA,STDEV,NISO)
DIMENSION O(15),BETA(15),STDEV(15)
40 FORMAT(8X,6(4X,A4,4X))
42 FORMAT(8HSOL(PPM),6(F12.4))
43 FORMAT(8HSTD DEV,6(F12.4))
LO=1
MAX=6
44 IF(NISO.LT.MAX)MAX=NISO
WRITE(2,40)(O(I),I=LO,MAX)
WRITE(2,42)(BETA(I),I=LO,MAX)
WRITE(2,43)(STDEV(I),I=LO,MAX)
IF(NISO.EQ.MAX)GO TO 50
LO=LO+6
MAX=MAX+6
GO TO 44
50 RETURN
END
SUBROUTINE VARNCE(A,OMG,STDEV,AT,LOCHN,NHICHN,NISO,MAXCON,C,P,O,BETA)
1IA,XNORM)
DIMENSION RES(400),OMG(400),STDEV(400),AT(400,15),P(15),O(15),BETA
1(15),C(15),A(400),XNORM(15)
DO 105 J=LOCHN,NHICHN
105 RES(J)=0.0
K=1
DO 101 I=1,MAXCON
IF(P(I).NE.0(K))GO TO 101
DO 102 J=LOCHN,NHICHN
102 RES(J)=RES(J)+((BETA(K)/C(K))**2)*AT(I,I)/XNORM(I)+(AT(J,I)*STDEV(
1(K)/C(K))**2
K=K+1
IF(K.GT.NISO)GO TO 103
101 CONTINUE
103 DO 104 J=LOCHN,NHICHN
104 A(J)=OMG(J)+RES(J)
RETURN
END
SUBROUTINE DELETE(NISO,REMOVE,ISOCHN,BETA,R,C,O,THETA,L)
DIMENSION RES(400),BETA(15),C(15),O(15),THETA(15),R(15,15)
30 FORMAT(/174 ISOTOPE REMOVED , A4 /)
WRITE(6,30)O(L)
NISO=NISO-1
BETAR=BETA(L)/C(L)
DO 60 J=1,ISOCHN
60 RES(J)=R(J,L)
REMOVE=REMOVE-(BETA(L)/C(L))**2/R(L,L)
DO 604 K=1,NISO
BETA(K)=(BETA(K)/C(K)-RES(K)*BETAR/RES(L))*C(K)
604 R(K,J)=R(K,J)-RES(K)*RES(J)/RES(L)

```

```

S HI F 1 15
S HI F 2 16
S HI F 3 17
S HI F 4 18
S HI F 5 19
S HI F 6 20
S HI F 7 21
S HI F 8 22
S HI F 9 23
S HI F 10 24
S HI F 11 25
S HI F 12 26
S HI F 13 27
S HI F 14 28
S HI F 15 29
S HI F 16 30
S HI F 17 31
S HI F 18 32
S HI F 19 33
S HI F 20 34
S HI F 21 35
S HI F 22 36
S HI F 23 37
S HI F 24 38
S HI F 25 39
S HI F 26 40
S HI F 27 41
S HI F 28 42
S HI F 29 43
S HI F 30 44
S HI F 31 45
S HI F 32 46
S HI F 33 47
S HI F 34 48
S HI F 35 49
S HI F 36 50
S HI F 37 51
S HI F 38 52
S HI F 39 53
S HI F 40 54
S HI F 41 55
S HI F 42 56
S HI F 43 57
S HI F 44 58
S HI F 45 59
S HI F 46 60
S HI F 47 61
S HI F 48 62
S HI F 49 63
S HI F 50 64
S HI F 51 65
S HI F 52 66
S HI F 53 67
S HI F 54 68
S HI F 55 69
S HI F 56 70
S HI F 57 71
S HI F 58 72
S HI F 59 73
S HI F 60 74
S HI F 61 75
S HI F 62 76
S HI F 63 77
S HI F 64 78
S HI F 65 79
S HI F 66 80
S HI F 67 81
S HI F 68 82
S HI F 69 83
S HI F 70 84
S HI F 71 85
S HI F 72 86
S HI F 73 87
S HI F 74 88
S HI F 75 89
S HI F 76 90
S HI F 77 91
S HI F 78 92
S HI F 79 93
S HI F 80 94
S HI F 81 95
S HI F 82 96
S HI F 83 97
S HI F 84 98
S HI F 85 99
S HI F 86 100
S HI F 87 101
S HI F 88 102
S HI F 89 103
S HI F 90 104
S HI F 91 105
S HI F 92 106
S HI F 93 107
S HI F 94 108
S HI F 95 109
S HI F 96 110
S HI F 97 111
S HI F 98 112
S HI F 99 113
S HI F 100 114
S HI F 101 115
S HI F 102 116
S HI F 103 117
S HI F 104 118
S HI F 105 119
S HI F 106 120
S HI F 107 121
S HI F 108 122
S HI F 109 123
S HI F 110 124
S HI F 111 125
S HI F 112 126
S HI F 113 127
S HI F 114 128
S HI F 115 129
S HI F 116 130
S HI F 117 131
S HI F 118 132
S HI F 119 133
S HI F 120 134
S HI F 121 135
S HI F 122 136
S HI F 123 137
S HI F 124 138
S HI F 125 139
S HI F 126 140
S HI F 127 141
S HI F 128 142
S HI F 129 143
S HI F 130 144
S HI F 131 145
S HI F 132 146
S HI F 133 147
S HI F 134 148
S HI F 135 149
S HI F 136 150
S HI F 137 151
S HI F 138 152
S HI F 139 153
S HI F 140 154
S HI F 141 155
S HI F 142 156
S HI F 143 157
S HI F 144 158
S HI F 145 159
S HI F 146 160
S HI F 147 161
S HI F 148 162
S HI F 149 163
S HI F 150 164
S HI F 151 165
S HI F 152 166
S HI F 153 167
S HI F 154 168
S HI F 155 169
S HI F 156 170
S HI F 157 171
S HI F 158 172
S HI F 159 173
S HI F 160 174
S HI F 161 175
S HI F 162 176
S HI F 163 177
S HI F 164 178
S HI F 165 179
S HI F 166 180
S HI F 167 181
S HI F 168 182
S HI F 169 183
S HI F 170 184
S HI F 171 185
S HI F 172 186
S HI F 173 187
S HI F 174 188
S HI F 175 189
S HI F 176 190
S HI F 177 191
S HI F 178 192
S HI F 179 193
S HI F 180 194
S HI F 181 195
S HI F 182 196
S HI F 183 197
S HI F 184 198
S HI F 185 199
S HI F 186 200
S HI F 187 201
S HI F 188 202
S HI F 189 203
S HI F 190 204
S HI F 191 205
S HI F 192 206
S HI F 193 207
S HI F 194 208
S HI F 195 209
S HI F 196 210
S HI F 197 211
S HI F 198 212
S HI F 199 213
S HI F 200 214
S HI F 201 215
S HI F 202 216
S HI F 203 217
S HI F 204 218
S HI F 205 219
S HI F 206 220
S HI F 207 221
S HI F 208 222
S HI F 209 223
S HI F 210 224
S HI F 211 225
S HI F 212 226
S HI F 213 227
S HI F 214 228
S HI F 215 229
S HI F 216 230
S HI F 217 231
S HI F 218 232
S HI F 219 233
S HI F 220 234
S HI F 221 235
S HI F 222 236
S HI F 223 237
S HI F 224 238
S HI F 225 239
S HI F 226 240
S HI F 227 241
S HI F 228 242
S HI F 229 243
S HI F 230 244
S HI F 231 245
S HI F 232 246
S HI F 233 247
S HI F 234 248
S HI F 235 249
S HI F 236 250
S HI F 237 251
S HI F 238 252
S HI F 239 253
S HI F 240 254
S HI F 241 255
S HI F 242 256
S HI F 243 257
S HI F 244 258
S HI F 245 259
S HI F 246 260
S HI F 247 261
S HI F 248 262
S HI F 249 263
S HI F 250 264
S HI F 251 265
S HI F 252 266
S HI F 253 267
S HI F 254 268
S HI F 255 269
S HI F 256 270
S HI F 257 271
S HI F 258 272
S HI F 259 273
S HI F 260 274
S HI F 261 275
S HI F 262 276
S HI F 263 277
S HI F 264 278
S HI F 265 279
S HI F 266 280
S HI F 267 281
S HI F 268 282
S HI F 269 283
S HI F 270 284
S HI F 271 285
S HI F 272 286
S HI F 273 287
S HI F 274 288
S HI F 275 289
S HI F 276 290
S HI F 277 291
S HI F 278 292
S HI F 279 293
S HI F 280 294
S HI F 281 295
S HI F 282 296
S HI F 283 297
S HI F 284 298
S HI F 285 299
S HI F 286 300
S HI F 287 301
S HI F 288 302
S HI F 289 303
S HI F 290 304
S HI F 291 305
S HI F 292 306
S HI F 293 307
S HI F 294 308
S HI F 295 309
S HI F 296 310
S HI F 297 311
S HI F 298 312
S HI F 299 313
S HI F 300 314
S HI F 301 315
S HI F 302 316
S HI F 303 317
S HI F 304 318
S HI F 305 319
S HI F 306 320
S HI F 307 321
S HI F 308 322
S HI F 309 323
S HI F 310 324
S HI F 311 325
S HI F 312 326
S HI F 313 327
S HI F 314 328
S HI F 315 329
S HI F 316 330
S HI F 317 331
S HI F 318 332
S HI F 319 333
S HI F 320 334
S HI F 321 335
S HI F 322 336
S HI F 323 337
S HI F 324 338
S HI F 325 339
S HI F 326 340
S HI F 327 341
S HI F 328 342
S HI F 329 343
S HI F 330 344
S HI F 331 345
S HI F 332 346
S HI F 333 347
S HI F 334 348
S HI F 335 349
S HI F 336 350
S HI F 337 351
S HI F 338 352
S HI F 339 353
S HI F 340 354
S HI F 341 355
S HI F 342 356
S HI F 343 357
S HI F 344 358
S HI F 345 359
S HI F 346 360
S HI F 347 361
S HI F 348 362
S HI F 349 363
S HI F 350 364
S HI F 351 365
S HI F 352 366
S HI F 353 367
S HI F 354 368
S HI F 355 369
S HI F 356 370
S HI F 357 371
S HI F 358 372
S HI F 359 373
S HI F 360 374
S HI F 361 375
S HI F 362 376
S HI F 363 377
S HI F 364 378
S HI F 365 379
S HI F 366 380
S HI F 367 381
S HI F 368 382
S HI F 369 383
S HI F 370 384
S HI F 371 385
S HI F 372 386
S HI F 373 387
S HI F 374 388
S HI F 375 389
S HI F 376 390
S HI F 377 391
S HI F 378 392
S HI F 379 393
S HI F 380 394
S HI F 381 395
S HI F 382 396
S HI F 383 397
S HI F 384 398
S HI F 385 399
S HI F 386 400
S HI F 387 401
S HI F 388 402
S HI F 389 403
S HI F 390 404
S HI F 391 405
S HI F 392 406
S HI F 393 407
S HI F 394 408
S HI F 395 409
S HI F 396 410
S HI F 397 411
S HI F 398 412
S HI F 399 413
S HI F 400 414
S HI F 401 415
S HI F 402 416
S HI F 403 417
S HI F 404 418
S HI F 405 419
S HI F 406 420
S HI F 407 421
S HI F 408 422
S HI F 409 423
S HI F 410 424
S HI F 411 425
S HI F 412 426
S HI F 413 427
S HI F 414 428
S HI F 415 429
S HI F 416 430
S HI F 417 431
S HI F 418 432
S HI F 419 433
S HI F 420 434
S HI F 421 435
S HI F 422 436
S HI F 423 437
S HI F 424 438
S HI F 425 439
S HI F 426 440
S HI F 427 441
S HI F 428 442
S HI F 429 443
S HI F 430 444
S HI F 431 445
S HI F 432 446
S HI F 433 447
S HI F 434 448
S HI F 435 449
S HI F 436 450
S HI F 437 451
S HI F 438 452
S HI F 439 453
S HI F 440 454
S HI F 441 455
S HI F 442 456
S HI F 443 457
S HI F 444 458
S HI F 445 459
S HI F 446 460
S HI F 447 461
S HI F 448 462
S HI F 449 463
S HI F 450 464
S HI F 451 465
S HI F 452 466
S HI F 453 467
S HI F 454 468
S HI F 455 469
S HI F 456 470
S HI F 457 471
S HI F 458 472
S HI F 459 473
S HI F 460 474
S HI F 461 475
S HI F 462 476
S HI F 463 477
S HI F 464 478
S HI F 465 479
S HI F 466 480
S HI F 467 481
S HI F 468 482
S HI F 469 483
S HI F 470 484
S HI F 471 485
S HI F 472 486
S HI F 473 487
S HI F 474 488
S HI F 475 489
S HI F 476 490
S HI F 477 491
S HI F 478 492
S HI F 479 493
S HI F 480 494
S HI F 481 495
S HI F 482 496
S HI F 483 497
S HI F 484 498
S HI F 485 499
S HI F 486 500
S HI F 487 501
S HI F 488 502
S HI F 489 503
S HI F 490 504
S HI F 491 505
S HI F 492 506
S HI F 493 507
S HI F 494 508
S HI F 495 509
S HI F 496 510
S HI F 497 511
S HI F 498 512
S HI F 499 513
S HI F 500 514
S HI F 501 515
S HI F 502 516
S HI F 503 517
S HI F 504 518
S HI F 505 519
S HI F 506 520
S HI F 507 521
S HI F 508 522
S HI F 509 523
S HI F 510 524
S HI F 511 525
S HI F 512 526
S HI F 513 527
S HI F 514 528
S HI F 515 529
S HI F 516 530
S HI F 517 531
S HI F 518 532
S HI F 519 533
S HI F 520 534
S HI F 521 535
S HI F 522 536
S HI F 523 537
S HI F 524 538
S HI F 525 539
S HI F 526 540
S HI F 527 541
S HI F 528 542
S HI F 529 543
S HI F 530 544
S HI F 531 545
S HI F 532 546
S HI F 533 547
S HI F 534 548
S HI F 535 549
S HI F 536 550
S HI F 537 551
S HI F 538 552
S HI F 539 553
S HI F 540 554
S HI F 541 555
S HI F 542 556
S HI F 543 557
S HI F 544 558
S HI F 545 559
S HI F 546 560
S HI F 547 561
S HI F 548 562
S HI F 549 563
S HI F 550 564
S HI F 551 565
S HI F 552 566
S HI F 553 567
S HI F 554 568
S HI F 555 569
S HI F 556 570
S HI F 557 571
S HI F 558 572
S HI F 559 573
S HI F 560 574
S HI F 561 575
S HI F 562 576
S HI F 563 577
S HI F 564 578
S HI F 565 579
S HI F 566 580
S HI F 567 581
S HI F 568 582
S HI F 569 583
S HI F 570 584
S HI F 571 585
S HI F 572 586
S HI F 573 587
S HI F 574 588
S HI F 575 589
S HI F 576 590
S HI F 577 591
S HI F 578 592
S HI F 579 593
S HI F 580 594
S HI F 581 595
S HI F 582 596
S HI F 583 597
S HI F 584 598
S HI F 585 599
S HI F 586 600
S HI F 587 601
S HI F 588 602
S HI F 589 603
S HI F 590 604
S HI F 591 605
S HI F 592 606
S HI F 593 607
S HI F 594 608
S HI F 595 609
S HI F 596 610
S HI F 597 611
S HI F 598 612
S HI F 599 613
S HI F 600 614
S HI F 601 615
S HI F 602 616
S HI F 603 617
S HI F 604 618
S HI F 605 619
S HI F 606 620
S HI F 607 621
S HI F 608 622
S HI F 609 623
S HI F 610 624
S HI F 611 625
S HI F 612 626
S HI F 613 627
S HI F 614 628
S HI F 615 629
S HI F 616 630
S HI F 617 631
S HI F 618 632
S HI F 619 633
S HI F 620 634
S HI F 621 635
S HI F 622 636
S HI F 623 637
S HI F 624 638
S HI F 625 639
S HI F 626 640
S HI F 627 641
S HI F 628 642
S HI F 629 643
S HI F 630 644
S HI F 631 645
S HI F 632 646
S HI F 633 647
S HI F 634 648
S HI F 635 649
S HI F 636 650
S HI F 637 651
S HI F 638 652
S HI F 639 653
S HI F 640 654
S HI F 641 655
S HI F 642 656
S HI F 643 657
S HI F 644 658
S HI F 645 659
S HI F 646 660
S HI F 647 661
S HI F 648 662
S HI F 649 663
S HI F 650 664
S HI F 651 665
S HI F 652 666
S HI F 653 667
S HI F 654 668
S HI F 655 669
S HI F 656 670
S HI F 657 671
S HI F 658 672
S HI F 659 673
S HI F 660 674
S HI F 661 675
S HI F 662 676
S HI F 663 677
S HI F 664 678
S HI F 665 679
S HI F 666 680
S HI F 667 681
S HI F 668 682
S HI F 669 683
S HI F 670 684
S HI F 671 685
S HI F 672 686
S HI F 673 687
S HI F 674 688
S HI F 675 689
S HI F 676 690
S HI F 677 691
S HI F 678 692
S HI F 679 693
S HI F 680 694
S HI F 681 695
S HI F 682 696
S HI F 683 697
S HI F 684 698
S HI F 685 699
S HI F 686 700
S HI F 687 701
S HI F 688 702
S HI F 689 703
S HI F 690 704
S HI F 691 705
S HI F 692 706
S HI F 693 707
S HI F 694 708
S HI F 695 709
S HI F 696 710
S HI F 697 711
S HI F 698 712
S HI F 699 713
S HI F 700 714
S HI F 701 715
S HI F 702 716
S HI F 703 717
S HI F 704 718
S HI F 705 719
S HI F 706 720
S HI F 707 721
S HI F 708 722
S HI F 709 723
S HI F 710 724
S HI F 711 725
S HI F 712 726
S HI F 713 727
S HI F 714 728
S HI F 715 729
S HI F 716 730
S HI F 717 731
S HI F 718 732
S HI F 719 733
S HI F 720 734
S HI F 721 735
S HI F 722 736
S HI F 723 737
S HI F 724 738
S HI F 725 739
S HI F 726 740
S HI F 727 741
S HI F 728 742
S HI F 729 743
S HI F 730 744
S HI F 731 745
S HI F 732 746
S HI F 733 747
S HI F 734 748
S HI F 735 749
S HI F 736 750
S HI F 737 751
S HI F 738 752
S HI F 739 753
S HI F 740 754
S HI F 741 755
S HI F 742 756
S HI F 743 757
S HI F 744 758
S HI F 745 759
S HI F 746 760
S HI F 747 761
S HI F 748 762
S HI F 749 763
S HI F 750 764
S HI F 751 765
S HI F 752 766
S HI F 753 767
S HI F 754 768
S HI F 755 769
S HI F 756 770
S HI F 757 771
S HI F 758 772
S HI F 759 773
S HI F 760 774
S HI F 761 775
S HI F 762 776
S HI F 763 777
S HI F 764 778
S HI F 765 779
S HI F 766 780
S HI F 767 781
S HI F 768 782
S HI F 769 783
S HI F 770 784
S HI F 771 785
S HI F 772 786
S HI F 773 787
S HI F 774 788
S HI F 775 789
S HI F 776 790
S HI F 777 791
S HI F 778 792
S HI F 779 793
S HI F 780 794
S HI F 781 795
S HI F 782 796
S HI F 783 797
S HI F 784 798
S HI F 785 799
S HI F 786 800
S HI F 787 801
S HI F 788 802
S HI F 789 803
S HI F 790 804
S HI F 791 805
S HI F 792 806
S HI F 793 807
S HI F 794 808
S HI F 795 809
S HI F 796 810
S HI F 797 811
S HI F 798 812
S HI F 799 813
S HI F 800 814
S HI F 801 815
S HI F 802 816
S HI F 803 817
S HI F 804 818
S HI F 805 819
S HI F 806 820
S HI F 807 821
S HI F 808 822
S HI F 809 823
S HI F 810 824
S HI F 811 825
S HI F 812 826
S HI F 813 827
S HI F 814 828
S HI F 815 829
S HI F 816 830
S HI F 817 831
S HI F 818 832
S HI F 819 833
S HI F 820 834
S HI F 821 835
S HI F 822 836
S HI F 823 837
S HI F 824 838
S HI F 825 839
S HI F 826 840
S HI F 827 841
S HI F 828 842
S HI F 829 843
S HI F 830 844
S HI F 831 845
S HI F 832 846
S HI F 833 847
S HI F 834 848
S HI F 835 849
S HI F 836 850
S HI F 837 851
S HI F 838 852
S HI F 839 853
S HI F 840 854
S HI F 841 855
S HI F 842 856
S HI F 843 857
S HI F 844 858
S HI F 845 859
S HI F 846 860
S HI F 847 861
S HI F 848 862
S HI F 849 863
S HI F 850 864
S HI F 851 865
S HI F 852 866
S HI F 853 867
S HI F 854 868
S HI F 855 869
S HI F 856 870
S HI F 857 871
S HI F 858 872
S HI F 859 873
S HI F 860 874
S HI F 861 875
S HI F 862 876
S HI F 863 877
S HI F 864 878
S HI F 865 879
S HI F 866 880
S HI F 867 881
S HI F 868 882
S HI F 869 883
S HI F 870 884
S HI F 871 885
S HI F 872 886
S HI F 873 887
S HI F 874 888
S HI F 875 889
S HI F 876 890
S HI F 877 891
S HI F 878 892
S HI F 879 893
S HI F 880 894
S HI F 881 895
S HI F 882 896
S HI F 883 897
S HI F 884 898
S HI F 885 899
S HI F 886 900
S HI F 887 901
S HI F 888 902
S HI F 889 903
S HI F 890 904
S HI F 891 905
S HI F 892 906
S HI F 893 907
S HI F 894 908
S HI F 895 909
S HI F 896 910
S HI F 897 911
S HI F 898 912
S HI F 899 913
S HI F 900 914
S HI F 901 915
S HI F 902 916
S HI F 903 917
S HI F 904 918
S HI F 905 919
S HI F 906 920
S HI F 907 921
S HI F 908 922
S HI F 909 923
S HI F 910 924
S HI F 911 925
S HI F 912 926
S HI F 913 927
S HI F 914 928
S HI F 915 929
S HI F 916 930
S HI F 917 931
S HI F 918 932
S HI F 919 933
S HI F 920 934
S HI F 921 935
S HI F 922 936
S HI F 923 937
S HI F 924 938
S HI F 925 939
S HI F 926 940
S HI F 927 941
S HI F 928 942
S HI F 929 943
S HI F 930 944
S HI F 931 945
S HI F 932 946
S HI F 933 947
S HI F 934 948
S HI F 935 949
S HI F 936 950
S HI F 937 951
S HI F 938 952
S HI F 939 953
S HI F 940 954
S HI F 941 955
S HI F 942 956
S HI F 943 957
S HI F 944 958
S HI F 945 959
S HI F 946 960
S HI F 947 961
S HI F 948 962
S HI F 949 963
S HI F 950 964
S HI F 951 965
S HI F 952 966
S HI F 9
```

```

DO 602 K=L,NISO
C(K)=C(K+1)
S(K)=S(K+1)
THETA(K)=THETA(K+1)
BETA(K)=BETA(K+1)
DO 602 J=1,ISOCHN
R(J,K)=R(J,K+1)
DO 603 K=L,NISO
DO 603 J=1,NISO
R(K,J)=R(K+1,J)
RETURN
END
SUBROUTINE DET(R,N,DETERM)
DIMENSION R(15,15)
DETERMINANT EVALUATED USING PIVOTAL CONDENSATION
L = 2
L = 1
DO 20 I=K,N
RATIO = R(I,L)/R(L,L)
DO 20 J=K,N
R(I,J) = R(I,J)-R(L,J)*RATIO
CONTINUE
IF(K=N) 15,21,21
L = K
K = K+1
GO TO 5
DETERM = 1
DO 25 L=1,N
DETERM = DETERM*R(L,L)
RETURN
END
SUBROUTINE WRIT(NITER,ITER,O,BETA,STDEV,ST,NISO,F,RCORR)
DIMENSION O(15),BETA(15),STDEV(400),ST(15)
FORMAT(/10X,10(4X,A4,4X))
FORMAT(/11H ITERATION,I3,4H OF,I3,9H MAX ITER /)
FORMAT(10H SOL (PPM) ,10(E12.4))
FORMAT(10H STD DEV ,10(E12.4))
FORMAT(10H T VALUES ,10(E12.4))
FORMAT(/12HL F VALUE = ,E11.4,14H CORR COEF = ,E11.4/)
WRITE(6,41)NITER,ITER
MAX=NISO
IF(NISO.LE.10)GO TO 4450
MAX=10
WRITE(6,40)(O(I),I=1,MAX)
WRITE(6,42)(BETA(I),I=1,MAX)
WRITE(6,43)(STDEV(I),I=1,MAX)
WRITE(6,44)(ST(I),I=1,MAX)
IF(NISO.LE.10)GO TO 4451
WRITE(6,40)(O(I),I=11,NISO)
WRITE(6,42)(BETA(I),I=11,NISO)
WRITE(6,43)(STDEV(I),I=11,NISO)
WRITE(6,44)(ST(I),I=11,NISO)
WRITE(6,45)F,RCORR
RETURN
END
SUBROUTINE SAVE(NISO,O,OI,RI,RII)
DIMENSION O(15),OI(15),RI(15,15),RII(15,15)
DO 3180 I=1,NISO
OI(I)=O(I)
DO 3180 J=1,NISO
RII(I,J)=RI(I,J)
RETURN
END
SUBROUTINE XFORM(NFORM,NISO,R,RI,MAXCON,P,O,LOCHN,NHICHN,A,AT,NCORR)
1R,NWRIT,NISO,NINTF,RII,BETA,THETA,OI,RHO,III)
DIMENSION R(15,15),RI(15,15),RII(15,15),P(15),O(15),A(400),AT(400)
115),OI(15),O(15,15),S(15,15),BETA(15),THETA(15),RHO(400)
THIS SUB. FORMS OR READS A=TRANS * W * A
FORMAT(1H1)
FORMAT(14,(1X,10E12.5))
FORMAT(1X,A4,(1X,5E12.5))
FORMAT(/18H CORRECTION COEFF /)
FORMAT(/26H INTERFERENCES IN PER CENT /)
FORMAT(/21H CROSS PRODUCT MATRIX /)
FORMAT(3E24.18)
FORMAT(/35H INV CROSS PRODUCT MATRIX /)
FORMAT(/16H IDENTITY MATRIX /)
FORMAT(/13H INV ERROR = ,E14.7 /)
FORMAT(/20H EFF (IN PER CENT) = ,E18.12 /)
FORMAT(/6X,9(4X,A4,4X))
GO TO (98,99,40),NFORM
CALL PEAKS(AT,A,RHO,LOCHN,NHICHN,NISO,III)

```

14
15
16
17
18
19
20
21
22
23
24
25
1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20

```

98  GO TO 99
    READ(5,29)((R(I,J),I=1,NISO),J=1,NISO)
    CALL SAVE(NISO,0,0I,R,RII)
    GO TO 322
99  DO 102 J=1,NISO
    DO 102 I=1,NISO
    R(I,J)=0.
102  L=1
    DO 4052 J=1,MAXCON
    IF(P(J).NE.0(LL))GO TO 4052
    L=1
    DO 4053 K=1,MAXCON
    IF(P(K).NE.0(L))GO TO 4053
    DO 100 I=LOCHN,NHICHN
    IF(A(I).LE.0.001)A(I)=1.0
100  R(LL,L)=R(LL,L)*AT(I,J)*AT(I,K)/A(I)
    L=L+1
    IF(L.GT.NISO)GO TO 4054
4053 CONTINUE
4054 LL=LL+1
    IF(LL.GT.NISO)GO TO 4055
4052 CONTINUE
4055 WRITE(6,23)
    IF(NCORR.NE.0) GO TO 97
    WRITE(6,23)
    DO 311 J=1,NISO
    DO 311 I=1,NISO
311  RI(I,J)=((R(I,J)**2)/(R(I,I)*R(J,J)))
    WRITE(6,18)NISO,((RI(I,J),I=1,NISO),J=1,NISO)
    CALL DET(RI,NISO,DETERM)
    DETERM=100.*DETERM
    WRITE(6,37)DETERM
97  CALL SAVE(NISO,0,0I,R,RII)
    IF(NWRIT.NE.0) GO TO 107
    WRITE(6,9)
    WRITE(6,28)
    WRITE(6,18)NISO,((R(I,J),I=1,NISO),J=1,NISO)
107  CALL SMINV(R,NISO)
    IF(NWRIT.NE.0) GO TO 313
    WRITE(6,31)
    WRITE(6,18)NISO,((R(I,J),I=1,NISO),J=1,NISO)
    CALL SAVE(NISO,0I,0,R,RII)
    CALL SAVE(NISO,0,0I,R,RII)
313  IF(NIDEN.NE.0)GO TO 322
    LL=0
001  ERROR=0.
    DO 324 K=1,NISO
    DO 324 J=1,NISO
    Q(K,J)=0.
    DO 324 I=1,NISO
324  Q(K,J)=Q(K,J)+RI(K,I)*R(I,J)
    DO 323 K=1,NISO
    DO 323 J=1,NISO
323  ERROR=ERROR+ABS(Q(K,J))
    ERROR=(ERROR-FLOAT(NISO))/(FLOAT(NISO)**2)
    IF(LL.GT.3)GO TO 3360
    IF(ERROR.LT.0.000050)GO TO 3360
    ** ITERATE TO IMPROVE ACC
    DO 127 I=1,NISO
    DO 127 J=1,NISO
127  S(I,J)=Q(I,J)
    DO 128 J=1,NISO
128  S(J,J)=2.0+S(J,J)
    DO 129 K=1,NISO
    DO 129 J=1,NISO
    Q(K,J)=0.0
    DO 129 I=1,NISO
129  Q(K,J)=Q(K,J)+R(K,I)*S(I,J)
    DO 130 J=1,NISO
    DO 130 K=1,NISO
130  R(K,J)=Q(K,J)
    LL=LL+1
    GO TO 4001
4001 WRITE(6,33)ERROR
    WRITE(6,32)
    WRITE(6,18)NISO,((Q(J,K),J=1,NISO),K=1,NISO)
322  IF(NINTF.NE.0)RETURN
    WRITE(6,9)
    WRITE(6,25)
    MIN=2
    MAX=NISO=1

```

FORM 2
 FORM 3
 FORM 4
 FORM 5
 FORM 6
 FORM 7
 FORM 8
 FORM 9
 FORM 10
 FORM 11
 FORM 12
 FORM 13
 FORM 14
 FORM 15
 FORM 16
 FORM 17
 FORM 18
 FORM 19
 FORM 20
 FORM 21
 FORM 22
 FORM 23
 FORM 24
 FORM 25
 FORM 26
 FORM 27
 FORM 28
 FORM 29
 FORM 30
 FORM 31
 FORM 32
 FORM 33
 FORM 34
 FORM 35
 FORM 36
 FORM 37
 FORM 38
 FORM 39
 FORM 40
 FORM 41
 FORM 42
 FORM 43
 FORM 44
 FORM 45
 FORM 46
 FORM 47
 FORM 48
 FORM 49
 FORM 50
 FORM 51
 FORM 52
 FORM 53
 FORM 54
 FORM 55
 FORM 56
 FORM 57
 FORM 58
 FORM 59
 FORM 60
 FORM 61
 FORM 62
 FORM 63
 FORM 64
 FORM 65
 FORM 66
 FORM 67
 FORM 68
 FORM 69
 FORM 70
 FORM 71
 FORM 72
 FORM 73
 FORM 74
 FORM 75
 FORM 76
 FORM 77
 FORM 78
 FORM 79
 FORM 80
 FORM 81
 FORM 82
 FORM 83
 FORM 84
 FORM 85
 FORM 86
 FORM 87
 FORM 88
 FORM 89
 FORM 90
 FORM 91
 FORM 92
 FORM 93
 FORM 94
 FORM 95
 FORM 96
 FORM 97
 FORM 98
 FORM 99
 FORM 00


```

DO 315 I=1,MAX
DO 316 J=MIN,NISO
BETA(I,J)=0
316 THETA(J)=((R(I,I)*R(J,J))*2)/(R(I,I)*R(J,J))*100.
WRITE(6,39) (BETA(J),J=MIN,NISO)
WRITE(6,22) O(I), (THETA(J),J=MIN,NISO)
315 MIN=MIN+1
RETURN
END
SUBROUTINE WATE(NVAR,NUMVAR,NCHAN,RHO,XNCOM,III,OMG,A,RES,ISOCHN)
DIMENSION RES(400),RHO(400),OMG(400),A(400)
FORMAT(1X,2A4)
47 FORMAT(/22H AVE WT FACTORS==USES ,I3.4H IN NO /)
48 FORMAT(/19H WT FACTORS = UNITY /)
49 FORMAT(/15H WT FACTORS = Y /)
THIS SUB CALCULATES OR READS THE WT FACTORS
READ OR SET VARIANCES USING (A) AS STORAGE
GO TO (201,203,3010),NVAR
201 DO 2202 J=1,400
2202 RES(J)=0.0
WRITE(6,47) NUMVAR
DO 2201 J=1,NUMVAR
READ(5,11) (A(I),I=1,NCHAN)
DO 2222 K=1,NCHAN
RES(K)=RES(K)+A(K)
2222 CONTINUE
DO 2224 I=1,NCHAN
RES(I)=RES(I)/FLOAT(NUMVAR)
WRITE(2,11) (RES(I),I=1,NCHAN)
DO 2223 I=1,NCHAN
A(I)=RES(I)
203 GO TO 205
213 DO 213 I=1,ISOCHN
A(I)=1.0
WRITE(6,48)
GO TO 205
3010 DO 3014 I=1,ISOCHN
A(I)=RHO(I)/XNCOM
3014 IF (A(I).EQ.0.) A(I)=1.0
WRITE(6,49)
IF (NVAR.EQ.2) GO TO 254
IF (III.GT.1) GO TO 254
256 DO 255 I=1,ISOCHN
255 OMG(I)=A(I)
254 CONTINUE
RETURN
END
SUBROUTINE CHI(CALC,RES,MAXCON,P,O,LOCHN,NHICHN,BETA,AT,C,RHO,NISO)
1,QUAL,A)
DIMENSION CALC(400),RES(400),P(15),O(15),BETA(15),AT(400,15),C(15),
1,RHO(400),A(400)
10 FORMAT(1X,2A4)
14 FORMAT(18H DEG OF FREEDOM = ,F4.0,15H QUAL OF FIT = ,E10.3)
18 FORMAT(/14,(1X,10E12.5))
DO 105 I=LOCHN,NHICHN
CALC(I)=0.
105 RES(I)=0.
K=1
DO 956 I=1,MAXCON
IF (P(I).NE.0.(K)) GO TO 956
DO 970 J=LOCHN,NHICHN
970 CALC(J)=BETA(K)*AT(J,I)/C(K)+CALC(J)
K=K+1
IF (K.GT.NISO) GO TO 957
956 CONTINUE
957 QUAL=0.
DO 104 I=LOCHN,NHICHN
104 RES(I)=RHO(I)-CALC(I)
DO 971 J=LOCHN,NHICHN
971 QUAL=QUAL+(RES(J)**2)/A(J)
QUAL=QUAL/FLOAT(NHICHN-LOCHN-NISO)
DF=NHICHN-LOCHN-NISO
WRITE(6,14) DF,QUAL
RETURN
END
SUBROUTINE PEAKS(AT,A,RHO,LOCHN,NHICHN,NISO,III)
DIMENSION AT(400,15),A(400),RHO(400),NK(40),T(5),C(5),CT(400,5),ISO
41 FORMAT(/14H COMPLEX SPECT /)
42 FORMAT(25I3)
43 FORMAT(1X,5E14.7)

```

```

FORM 10
FORM 103
FORM 104
FORM 105
FORM 106
FORM 107
FORM 108
FORM 109
FORM 110
WATE 111
WATE 112
WATE 113
WATE 114
WATE 115
WATE 116
WATE 117
WATE 118
WATE 119
WATE 120
WATE 121
WATE 122
WATE 123
WATE 124
WATE 125
WATE 126
WATE 127
WATE 128
WATE 129
WATE 130
WATE 131
WATE 132
WATE 133
WATE 134
WATE 135
WATE 136
WATE 137
WATE 138
WATE 139
WATE 140
WATE 141
WATE 142
WATE 143
WATE 144
WATE 145
WATE 146
WATE 147
WATE 148
WATE 149
WATE 150
WATE 151
WATE 152
WATE 153
WATE 154
WATE 155
WATE 156
WATE 157
WATE 158
WATE 159
WATE 160
WATE 161
WATE 162
WATE 163
WATE 164
WATE 165
WATE 166
WATE 167
WATE 168
WATE 169
WATE 170
WATE 171
WATE 172
WATE 173
WATE 174
WATE 175
WATE 176
WATE 177
WATE 178
WATE 179
WATE 180
WATE 181
WATE 182
WATE 183
WATE 184
WATE 185
WATE 186
WATE 187
WATE 188
WATE 189
WATE 190
WATE 191
WATE 192
WATE 193
WATE 194
WATE 195
WATE 196
WATE 197
WATE 198
WATE 199
WATE 200
WATE 201
WATE 202
WATE 203
WATE 204
WATE 205
WATE 206
WATE 207
WATE 208
WATE 209
WATE 210
WATE 211
WATE 212
WATE 213
WATE 214
WATE 215
WATE 216
WATE 217
WATE 218
WATE 219
WATE 220
WATE 221
WATE 222
WATE 223
WATE 224
WATE 225
WATE 226
WATE 227
WATE 228
WATE 229
WATE 230
WATE 231
WATE 232
WATE 233
WATE 234
WATE 235
WATE 236
WATE 237
WATE 238
WATE 239
WATE 240
WATE 241
WATE 242
WATE 243
WATE 244
WATE 245
WATE 246
WATE 247
WATE 248
WATE 249
WATE 250
WATE 251
WATE 252
WATE 253
WATE 254
WATE 255
WATE 256
WATE 257
WATE 258
WATE 259
WATE 260
WATE 261
WATE 262
WATE 263
WATE 264
WATE 265
WATE 266
WATE 267
WATE 268
WATE 269
WATE 270
WATE 271
WATE 272
WATE 273
WATE 274
WATE 275
WATE 276
WATE 277
WATE 278
WATE 279
WATE 280
WATE 281
WATE 282
WATE 283
WATE 284
WATE 285
WATE 286
WATE 287
WATE 288
WATE 289
WATE 290
WATE 291
WATE 292
WATE 293
WATE 294
WATE 295
WATE 296
WATE 297
WATE 298
WATE 299
WATE 300
WATE 301
WATE 302
WATE 303
WATE 304
WATE 305
WATE 306
WATE 307
WATE 308
WATE 309
WATE 310
WATE 311
WATE 312
WATE 313
WATE 314
WATE 315
WATE 316
WATE 317
WATE 318
WATE 319
WATE 320
WATE 321
WATE 322
WATE 323
WATE 324
WATE 325
WATE 326
WATE 327
WATE 328
WATE 329
WATE 330
WATE 331
WATE 332
WATE 333
WATE 334
WATE 335
WATE 336
WATE 337
WATE 338
WATE 339
WATE 340
WATE 341
WATE 342
WATE 343
WATE 344
WATE 345
WATE 346
WATE 347
WATE 348
WATE 349
WATE 350
WATE 351
WATE 352
WATE 353
WATE 354
WATE 355
WATE 356
WATE 357
WATE 358
WATE 359
WATE 360
WATE 361
WATE 362
WATE 363
WATE 364
WATE 365
WATE 366
WATE 367
WATE 368
WATE 369
WATE 370
WATE 371
WATE 372
WATE 373
WATE 374
WATE 375
WATE 376
WATE 377
WATE 378
WATE 379
WATE 380
WATE 381
WATE 382
WATE 383
WATE 384
WATE 385
WATE 386
WATE 387
WATE 388
WATE 389
WATE 390
WATE 391
WATE 392
WATE 393
WATE 394
WATE 395
WATE 396
WATE 397
WATE 398
WATE 399
WATE 400
WATE 401
WATE 402
WATE 403
WATE 404
WATE 405
WATE 406
WATE 407
WATE 408
WATE 409
WATE 410
WATE 411
WATE 412
WATE 413
WATE 414
WATE 415
WATE 416
WATE 417
WATE 418
WATE 419
WATE 420
WATE 421
WATE 422
WATE 423
WATE 424
WATE 425
WATE 426
WATE 427
WATE 428
WATE 429
WATE 430
WATE 431
WATE 432
WATE 433
WATE 434
WATE 435
WATE 436
WATE 437
WATE 438
WATE 439
WATE 440
WATE 441
WATE 442
WATE 443
WATE 444
WATE 445
WATE 446
WATE 447
WATE 448
WATE 449
WATE 450
WATE 451
WATE 452
WATE 453
WATE 454
WATE 455
WATE 456
WATE 457
WATE 458
WATE 459
WATE 460
WATE 461
WATE 462
WATE 463
WATE 464
WATE 465
WATE 466
WATE 467
WATE 468
WATE 469
WATE 470
WATE 471
WATE 472
WATE 473
WATE 474
WATE 475
WATE 476
WATE 477
WATE 478
WATE 479
WATE 480
WATE 481
WATE 482
WATE 483
WATE 484
WATE 485
WATE 486
WATE 487
WATE 488
WATE 489
WATE 490
WATE 491
WATE 492
WATE 493
WATE 494
WATE 495
WATE 496
WATE 497
WATE 498
WATE 499
WATE 500

```



```

IF(I.EQ.0.N) GO TO 60
JJ=I+1
DO 72 N K=JJ,N
IK=INDEX(I,K,N)
IF(I.GT.1) GO TO 62
I2=N*(N+1)/2+K
A(I2)=A(IJ)*A(IK)
GO TO 72
62 I4K=INDEX(I4,K,N)
A(I4K)=A(IJ)*A(IK)
72 CONTINUE
DO 120 N J=JJ,N
DO 120 N K=J,N
JK=INDEX(J,K,N)
IJ=INDEX(I,J,N)
IF(I.GT.1) GO TO 82
I2=N*(N+1)/2+K
A(IJK)=A(IJK)+A(I2)*A(IJ)
GO TO 120
82 I4K=INDEX(I4,K,N)
A(IJK)=A(IJK)+A(I4K)*A(IJ)
120 CONTINUE
CONTINUE
60 DO 240 MM=2,N
I=N-MM+1
I4=I-1
I4I=INDEX(I4,I,N)
II=INDEX(I,I,N)
JJ=I+1
DO 230 N J=JJ,N
IJ=INDEX(I,J,N)
A(IJ)=0.0
DO 220 N J=JJ,N
I4J=INDEX(I4,J,N)
IJ=INDEX(I,J,N)
DO 180 N K=J,N
I4K=INDEX(I4,K,N)
IK=INDEX(I,K,N)
JK=INDEX(J,K,N)
IF(I.EQ.1) GO TO 168
A(IK)=A(IK)+A(I4J)*A(JK)
IF(J.EQ.K) GO TO 180
A(IJ)=A(IJ)+A(I4K)*A(JK)
GO TO 180
168 I2=N*(N+1)/2+K
I3=N*(N+1)/2+J
A(IK)=A(IK)+A(I3)*A(JK)
IF(J.EQ.K) GO TO 180
A(IJ)=A(IJ)+A(I2)*A(JK)
180 CONTINUE
IF(I.EQ.1) GO TO 187
A(II)=A(II)+A(I4J)*A(IJ)
GO TO 220
187 A(II)=A(II)+A(I3)*A(IJ)
220 CONTINUE
IF(I.GT.1) GO TO 131
M=A(I4I)
GO TO 143
131 M=A(I4I)
143 IF(M.EQ.1) GO TO 240
DO 1105 N J=I,N
IJ=INDEX(I,J,N)
JM=INDEX(J,M,N)
BB=A(IJ)
A(IJ)=A(JM)
105 A(JM)=BB
IM=INDEX(I,M,N)
BB=A(IJ)
A(II)=A(IM)
240 A(IM)=BB
CONTINUE
MIN=1
K=0
DO 103 I=1,N
DO 104 J=MIN,N
K=K+1
R(I,J)=A(K)
104 MIN=MIN+1
103 MIN=2
MAX=N-1
DO 115 I=1,MAX
DO 102 J=MIN,N

```

```

SMIN 58
SMIN 59
SMIN 60
SMIN 61
SMIN 62
SMIN 63
SMIN 64
SMIN 65
SMIN 66
SMIN 67
SMIN 68
SMIN 69
SMIN 70
SMIN 71
SMIN 72
SMIN 73
SMIN 74
SMIN 75
SMIN 76
SMIN 77
SMIN 78
SMIN 79
SMIN 80
SMIN 81
SMIN 82
SMIN 83
SMIN 84
SMIN 85
SMIN 86
SMIN 87
SMIN 88
SMIN 89
SMIN 90
SMIN 91
SMIN 92
SMIN 93
SMIN 94
SMIN 95
SMIN 96
SMIN 97
SMIN 98
SMIN 99
SMIN 100
SMIN 101
SMIN 102
SMIN 103
SMIN 104
SMIN 105
SMIN 106
SMIN 107
SMIN 108
SMIN 109
SMIN 110
SMIN 111
SMIN 112
SMIN 113
SMIN 114
SMIN 115
SMIN 116
SMIN 117
SMIN 118
SMIN 119
SMIN 120
SMIN 121
SMIN 122
SMIN 123
SMIN 124
SMIN 125
SMIN 126
SMIN 127
SMIN 128
SMIN 129
SMIN 130
SMIN 131
SMIN 132
SMIN 133
SMIN 134
SMIN 135
SMIN 136
SMIN 137
SMIN 138
SMIN 139
SMIN 140

```

```

105 R(J,I)=R(I,J)
115 MIN=MIN+1
RETURN
300 WRITE (3,6)
STOP
END
FUNCTION INDEX(I8,J8,N)
FOR SYMMETRICAL MATRICES ONLY.
THIS SUBPROGRAM ALLOWS THE PROGRAMMER TO REFER TO A TWO-DIMENSIONAL
MATRIX WITH SUBSCRIPTS (I,J) WHILE THE MATRIX IS STORED AS A VECTOR.
I9= ROW OF MATRIX
JI= COLUMN OF MATRIX
N = ORDER OF MATRIX

I9=I8
JI=J8
IF(I9.LE.JI)GO TO 2
J9=JI
JI=I9
I9=J9
2 MI=0
IF(I9.LT.2)GO TO 4
DO 5 LI=2,I9
5 MI=MI+1-LI
4 INDEX=I9+JI-1+MI
RETURN
END
SUBROUTINE STAT(REMOVE,BETA,RHOSQD,NISO,LOCHN,NHICHN,R,C,O,AT,A,RES)
1 S,RHO,NVAR,NR,RAT,ST,NITER,F,RCORR,P,MAXCON,RHOSUM,XNORM,MVAR,ISO,STAT
2(I,I)
DIMENSION NR=5(400),BETA(15),R(15,15),C(15),O(15),AT(400,15),A(400),
1 RHO(400),S2(15),ST(15),P(15),XNORM(15),ISO(15),ATT(400,15)
28 FORMAT(/124 STD DEV OF ,A4,8H IS NEG /)
34 FORMAT(6X,12F6.0,2X )
IF(IIT.GT.1)GO TO 525
DO 526 J=1,MAXCON
DO 527 I=1,NHICHN
527 ATT(I,J)=0.0
LI=ISO(J)
DO 524 I=1,LI
524 ATT(I,J)=AT(I,J)
526 CONTINUE
VAREST=1.0
RCORR= REMOVE/RHOSQD
IF(RCORR.GT.0.9999)GO TO 502
JF=FLOAT(NHICHN-LOCHN+1)
VAREST=(RHOSQD-RHOSUM/JF)/DF
502 IF(MVAR.EQ.0)GO TO 522
DO 523 J=1,NISO
LI=ISO(J)
523 READ(5,34)(ATT(I,J),I=1,LI)
522 DO 123 I=1,NISO
SN(I)=0.0
LI=1
DO 122 J=1,MAXCON
IF(P(J).NE.0(L))GO TO 122
DO 102 K=LOCHN,NHICHN
102 S2(I)=S2(I)+(R(I,L)*RHO(K)/A(K))*2)*ATT(K,J)/XNORM(J)
LI=LI+1
IF(LI.GT.NISO)GO TO 123
122 CONTINUE
123 CONTINUE
IF(NVAR.EQ.2)GO TO 100
IF(NITER.GT.0)GO TO 107
IF(NVAR.NE.3)GO TO 111
DO 103 I=1,NISO
IF(R(I,I).GE.0.)GO TO 104
WRITE(6,28)O(I)
R(I,I)=ABS(R(I,I))
104 RES(I)=SORT(R(I,I)+S2(I))*C(I)
103 ST(I)=BETA(I)/RES(I)
IF(NVAR.EQ.3)GO TO 105
111 IF(NR.EQ.0)GO TO 107
VAREST=RAT
GO TO 100
107 DO 131 I=1,NISO
RES(I)=0.0
M=1
DO 110 J=1,MAXCON
IF(P(J).NE.0(M))GO TO 119
M=1
DO 120 L=1,MAXCON

```

MIN 40
MIN 41
MIN 42
MIN 43
MIN 44
MIN 45
INDEX 46
INDEX 47
INDEX 48
INDEX 49
INDEX 50
INDEX 51
INDEX 52
INDEX 53
INDEX 54
INDEX 55
INDEX 56
INDEX 57
INDEX 58
INDEX 59
INDEX 60
INDEX 61
INDEX 62
INDEX 63
INDEX 64
INDEX 65
INDEX 66
INDEX 67
INDEX 68
INDEX 69
INDEX 70
INDEX 71
INDEX 72
INDEX 73
INDEX 74
INDEX 75
INDEX 76
INDEX 77
INDEX 78
INDEX 79
INDEX 80
INDEX 81
INDEX 82
INDEX 83
INDEX 84
INDEX 85
INDEX 86
INDEX 87
INDEX 88
INDEX 89
INDEX 90
INDEX 91
INDEX 92
INDEX 93
INDEX 94
INDEX 95
INDEX 96
INDEX 97
INDEX 98
INDEX 99
INDEX 100
INDEX 101
INDEX 102
INDEX 103
INDEX 104
INDEX 105
INDEX 106
INDEX 107
INDEX 108
INDEX 109
INDEX 110
INDEX 111
INDEX 112
INDEX 113
INDEX 114
INDEX 115
INDEX 116
INDEX 117
INDEX 118
INDEX 119
INDEX 120
INDEX 121
INDEX 122
INDEX 123
INDEX 124
INDEX 125
INDEX 126
INDEX 127
INDEX 128
INDEX 129
INDEX 130
INDEX 131
INDEX 132
INDEX 133
INDEX 134
INDEX 135
INDEX 136
INDEX 137
INDEX 138
INDEX 139
INDEX 140
INDEX 141
INDEX 142
INDEX 143
INDEX 144
INDEX 145
INDEX 146
INDEX 147
INDEX 148
INDEX 149
INDEX 150

```

IF (P(L).NE.0(N))GO TO 129
DO 100 K=L,LOCHN,NHICHN
RES(I)=RES(I)+(AT(K,L)*AT(K,J)*RHO(K)*R(I,N)*R(I,M)/(A(K)**2))
NIM+1
IF (N.GT.NISO)GO TO 130
CONTINUE
129 M=M+1
IF (M.GT.NISO)GO TO 131
CONTINUE
131 DO 112 I=1,NISO
RES(I)=RES(I)+S2(I)*C(I)
S2(I)=S2(I)+R(I,I)/RES(I)
GO TO 100
DO 100 I=1,NISO
IF (R(I,I).GE.0.)GO TO 503
WRITE(6,28)O(I)
R(I,I)=ABS(R(I,I))
RES(I)=RES(I)+VAREST+S2(I)*C(I)
S2(I)=S2(I)+R(I,I)/RES(I)
VAREST=(REMOVE/LOAT(NISO))/VAREST
RETURN
END
SUBROUTINE INIT(NITER,RAT,SAV,C,ISOCHN,MODE,NBKGD,NUBKG,NSHFT,NPK,
1 LPB,NORM,NCTIN,NAREA,NCHCU,NCULL,NCUJ,NRHO,NUMVAR,NVAR,NREFS,IT,
2 NR,NW,LOCHN,NHICHN,NCHAN,MC,MVAR,T,AREA,BKGTIM,TIMIRR,FRHL,UP,
3 XNCOM)
DIMENSION C(15)
NITER=10
NITER=100
SAV=1.0
SAV=0.0
DO 15 I=1,15
C(I)=1.0
ISOCHN=256
MODE=2
NBKGD=1
NUBKG=1
NSHFT=1
NPK=3
LPB=1
NORM=1
NCTIN=1
NAREA=1
NCHCU=48
NCULL=33
NCUJ=40
NRHO=1
NUMVAR=1
NVAR=3
NREFS=0
IT=0
NITER=0
NREFS=0
NITER=0
NREFS=0
LOCHN=7
NHICHN=20
NCHAN=20
MC=1
MVAR=0
TIM=0.0
BKGTIM=6.0
TIMIRR=5.0
FRHL=75.0
UP=58.0
XNCOM=0.0
RETURN
END

```

```


STAT 56
STAT 57
STAT 58
STAT 59
STAT 60
STAT 61
STAT 62
STAT 63
STAT 64
STAT 65
STAT 66
STAT 67
STAT 68
STAT 69
STAT 70
STAT 71
STAT 72
STAT 73
STAT 74
STAT 75
STAT 76
STAT 77
STAT 78
INIT 1
INIT 2
INIT 3
INIT 4
INIT 5
INIT 6
INIT 7
INIT 8
INIT 9
INIT 10
INIT 11
INIT 12
INIT 13
INIT 14
INIT 15
INIT 16
INIT 17
INIT 18
INIT 19
INIT 20
INIT 21
INIT 22
INIT 23
INIT 24
INIT 25
INIT 26
INIT 27
INIT 28
INIT 29
INIT 30
INIT 31
INIT 32
INIT 33
INIT 34
INIT 35
INIT 36
INIT 37
INIT 38
INIT 39
INIT 40
INIT 41
INIT 42
INIT 43
INIT 44
INIT 45
INIT 46
INIT 47
INIT 48
INIT 49
INIT 50
INIT 51
INIT 52
INIT 53
INIT 54
INIT 55
INIT 56
INIT 57
INIT 58
INIT 59
INIT 60
INIT 61
INIT 62
INIT 63
INIT 64
INIT 65
INIT 66
INIT 67
INIT 68
INIT 69
INIT 70
INIT 71
INIT 72
INIT 73
INIT 74
INIT 75
INIT 76
INIT 77
INIT 78
INIT 79
INIT 80
INIT 81
INIT 82
INIT 83
INIT 84
INIT 85
INIT 86
INIT 87
INIT 88
INIT 89
INIT 90
INIT 91
INIT 92
INIT 93
INIT 94
INIT 95
INIT 96
INIT 97
INIT 98
INIT 99
INIT 100



```

N

APPENDIX J

-FLOW DIAGRAM OF THE CODE





APPENDIX K

- LISTING OF THE INPUT DATA FOR SAMPLE PROBLEM

The following is the listing of the input data cards used for the analysis of the sample problem data. They are listed in order of appearance in the program. For further use change of some of the control cards is sufficient. The data can be read from punched cards as well as from a storage file or from a magnetic tape. The listing of the input data is also given together with the read formats. The options are given in the comments cards of CORGAM and in the appendix F.

1. Card: (15A4) ID card for the complex spectrum

2. Card: (25I3)

NISO: Number of isotopes in complex spectra used in the analysis

NCASES: Number of complex spectra to be analyzed

NRES: Resonance calculation control parameter (optional)

NFORM: Least-squares fit control parameter (optional)

NWRIT: (X^oWX) matrix writing control parameter (optional)

NCORR: Correlation coefficient control parameter (optional)

NINTF: Interferences calculation control parameter (optional)

NIDEN: Identity matrix control parameter (optional)

3. Card: (25I3)

NCHG1: Change in NAM1 (optional)

NCHG2: Change in NAM2 (optional)

NCHG3: Change in NAM3 (optional)

(NAM1, NAM2, NAM3 are the NAMELIST variables)

4. Card: If change in NAM1 is required, new values are given by this card

5. Card: " NAM2 " " "

6. Card: " NAM3 " " "

7. Card: " " " " "

NAM1 contains ISOCHN, MODE, NBKGD, NUBKG, NSHFT, NPK, LPP, NORM, NCTIN
NHICHN, NCHAN, MC, MVAR

NAM2 contains C(I)

NAM3 contains T, DECAY, COUNT, WT, CUCT, WTCU, AREA, CUDEC, BKGTIM, TIMIRR,
FRHL, UP, XNCOM

The definitions of these variables are given in the comments of

CORGAM. These variables are entered in the program by a DATA statement, if a change is required, this is done in the 3. card and the new values are given in 4., 5., 6. and 7. cards. If a specific isotope wont be used in the analysis -1.0 will be entered in NAM2 for the C(I) value of the Ith isotope. For the complex spectra following the first one (if any), the input parameters are read from the subprogram INIT.

8. Card: (15A4) Complex spectrum ID card

9. Card up to 30. Card (included) (6X, 12(F6.0, 2X)) :

This is the complex spectrum data set. The variable name is RHO(I) where I is from 1 to ISOCHN. ISOCHN is the number of the channel contained in the complex spectrum. Maximum value of I can be 400, and 12 values are punched on each card. If the channel numbers are less than 400 (as in the case of sample problem 256), only the values prespecified are read.

31. Card: (15A4) Background ID card

32. Card up to 53. Card (included) (6X, 12(F6.0, 2X)): Background spectrum deck,
The variable name is BKG(I) where I is from 1 to ISOCHN

54. Card: Flux reference spectrum ID card

55. Card up to 58. Card (included) (6X, 12(F6.0, 2X)): Flux reference spectrum data set. The variable name is RES(I) where I is from 1 to NCHCU.

NCHCU is the channel number contained in the flux reference spectrum.

59. Card: (8F10.0) The gain of the system in MeV/ch. The variable name is GAIN

60. Card: (4(I3, 7X, F10.0)) IXO(I), E(I) I from 1 to NPK, where the variable IXO is the approximate mid-point of photopeak, E is the energy of the gamma photopeak in MeV. Both values are for the complex spectrum. NPK being the number of the photopeaks in the complex spectrum.

61. Card: (A4, 9I3) O(I), ISO(I), NPK, LPP, NORF, MODE, NFAC(I), NO

O(I): Tag name of the Ith isotope

ISO(I): Number of channels in the Ith reference isotope

NPK: is the photo-peaks contained in the complex spectrum

LPP: Index of photo-peak used in shifting

NORF: Specification for normalization of reference spectra

MODE: Reference spectrum input specification

NFAC(I): Decay time specification of reference spectrum

NO: Number of decay spectra of a specific reference isotope

62. Card: (8F10.0) H(I), D(I), FRRHL, XNORM(I)

H(I): Half-life, in minutes, of reference isotope spectrum

D(I): Irradiation time, in minutes, of reference sample spectrum

FRRHL: Reference spectrum, flux reference isotope half-life, in minutes

XNORM(I): Normalization constant of the Ith reference isotope

63.and64.Cards:(8F10.0) T(J),C(J) J:from 1 to NO

These data are for the reference decay spectra.where T(J) is the decay time in minutes of the Jth decay spectra and C(J) is the count time in minutes of the Jth decay spectra.NO being the number of the decay spectra.

65.Card:(25I3) NUBKG,NCULL,NCUUL,NCHCU,NCHBK

NUBKG:Background spectra specification

NCULL:Lower fit limit on copper reference photopeak

NCUUL:Upper " " " " "

NCHCU:Number of channels in copper reference spectrum

NCHBK:Number of channels in background spectrum

66.Card:(8F10.0) COUN,DECA,BKGTIM,CUDEC,WT,WTCU,CUCT

COUN:Counting time,in minutes,for reference spectrum

DECA:Decay time,in minutes,for reference spectrum

BKGTIM:Background counting time,in minutes,for reference spectrum

CUDEC:Decay time,in minutes,for flux reference spectrum

WT:Weight,in grams,of the reference isotope

WTCU:Weight,in grams,of the flux reference isotope

CUCT:Counting time,in minutes,for the reference spectrum,flux referenc spectrum.

67.Card up to 186.Card (5E14.8) CT(J,M) J from 1 to ISOCHN, M from 1 to NO

These are the decay spectra of the reference isotopes,where CT(J,M) is the counts registered at Jth channel of the reference isotope Mth decay time.NO is the number of decay spectra.If an isotope doesn't have a decay spectra M is assumed to be equal to 1.

187.Card:same as card 61

188.Card:same as card 62

189.and 190.Cards:same as cards 63 and 64

191.Card:same as card 65

192.Card:same as card 66

193.Card and the following deck:same as 67. card and the following deck.

The same order repeats itself till all the reference isotopes are read in.In the sample problem data this repetition ends at the 458. card.From this point on the next complex spectrum data is given, using the same formats as in the first spectrum.

459.Card:same as card 3

460.Card:same as card 4

461.Card:same as card 5


462.Card:same as card 6

463. Card: same as card 7

464. Card: same as card 8

465. Card up to 486. Card: same as the 9. card up to 30. card

The same order repeats itself till all the complex spectra are read in.






Table with multiple columns of numbers and text. Includes headers like 'CA 40258', '420120', and various numerical sequences and codes.

| Index | Category | Value | Label |
|-------|----------|-------|-------|
| 1 | AC | 20 | AC |
| 2 | AM | 20 | AM |
| 3 | AM | 20 | AM |
| 4 | AM | 20 | AM |
| 5 | AM | 20 | AM |
| 6 | AM | 20 | AM |
| 7 | AM | 20 | AM |
| 8 | AM | 20 | AM |
| 9 | AM | 20 | AM |
| 10 | AM | 20 | AM |
| 11 | AM | 20 | AM |
| 12 | AM | 20 | AM |
| 13 | AM | 20 | AM |
| 14 | AM | 20 | AM |
| 15 | AM | 20 | AM |
| 16 | AM | 20 | AM |
| 17 | AM | 20 | AM |
| 18 | AM | 20 | AM |
| 19 | AM | 20 | AM |
| 20 | AM | 20 | AM |
| 21 | AM | 20 | AM |
| 22 | AM | 20 | AM |
| 23 | AM | 20 | AM |
| 24 | AM | 20 | AM |
| 25 | AM | 20 | AM |
| 26 | AM | 20 | AM |
| 27 | AM | 20 | AM |
| 28 | AM | 20 | AM |
| 29 | AM | 20 | AM |
| 30 | AM | 20 | AM |
| 31 | AM | 20 | AM |
| 32 | AM | 20 | AM |
| 33 | AM | 20 | AM |
| 34 | AM | 20 | AM |
| 35 | AM | 20 | AM |
| 36 | AM | 20 | AM |
| 37 | AM | 20 | AM |
| 38 | AM | 20 | AM |
| 39 | AM | 20 | AM |
| 40 | AM | 20 | AM |
| 41 | AM | 20 | AM |
| 42 | AM | 20 | AM |
| 43 | AM | 20 | AM |
| 44 | AM | 20 | AM |
| 45 | AM | 20 | AM |
| 46 | AM | 20 | AM |
| 47 | AM | 20 | AM |
| 48 | AM | 20 | AM |
| 49 | AM | 20 | AM |
| 50 | AM | 20 | AM |
| 51 | AM | 20 | AM |
| 52 | AM | 20 | AM |
| 53 | AM | 20 | AM |
| 54 | AM | 20 | AM |
| 55 | AM | 20 | AM |
| 56 | AM | 20 | AM |
| 57 | AM | 20 | AM |
| 58 | AM | 20 | AM |
| 59 | AM | 20 | AM |
| 60 | AM | 20 | AM |
| 61 | AM | 20 | AM |
| 62 | AM | 20 | AM |
| 63 | AM | 20 | AM |
| 64 | AM | 20 | AM |
| 65 | AM | 20 | AM |
| 66 | AM | 20 | AM |
| 67 | AM | 20 | AM |
| 68 | AM | 20 | AM |
| 69 | AM | 20 | AM |
| 70 | AM | 20 | AM |
| 71 | AM | 20 | AM |
| 72 | AM | 20 | AM |
| 73 | AM | 20 | AM |
| 74 | AM | 20 | AM |
| 75 | AM | 20 | AM |
| 76 | AM | 20 | AM |
| 77 | AM | 20 | AM |
| 78 | AM | 20 | AM |
| 79 | AM | 20 | AM |
| 80 | AM | 20 | AM |
| 81 | AM | 20 | AM |
| 82 | AM | 20 | AM |
| 83 | AM | 20 | AM |
| 84 | AM | 20 | AM |
| 85 | AM | 20 | AM |
| 86 | AM | 20 | AM |
| 87 | AM | 20 | AM |
| 88 | AM | 20 | AM |
| 89 | AM | 20 | AM |
| 90 | AM | 20 | AM |
| 91 | AM | 20 | AM |
| 92 | AM | 20 | AM |
| 93 | AM | 20 | AM |
| 94 | AM | 20 | AM |
| 95 | AM | 20 | AM |
| 96 | AM | 20 | AM |
| 97 | AM | 20 | AM |
| 98 | AM | 20 | AM |
| 99 | AM | 20 | AM |
| 100 | AM | 20 | AM |
| 101 | AM | 20 | AM |
| 102 | AM | 20 | AM |
| 103 | AM | 20 | AM |
| 104 | AM | 20 | AM |
| 105 | AM | 20 | AM |
| 106 | AM | 20 | AM |
| 107 | AM | 20 | AM |
| 108 | AM | 20 | AM |
| 109 | AM | 20 | AM |
| 110 | AM | 20 | AM |
| 111 | AM | 20 | AM |
| 112 | AM | 20 | AM |
| 113 | AM | 20 | AM |
| 114 | AM | 20 | AM |
| 115 | AM | 20 | AM |
| 116 | AM | 20 | AM |
| 117 | AM | 20 | AM |
| 118 | AM | 20 | AM |
| 119 | AM | 20 | AM |
| 120 | AM | 20 | AM |
| 121 | AM | 20 | AM |
| 122 | AM | 20 | AM |
| 123 | AM | 20 | AM |
| 124 | AM | 20 | AM |
| 125 | AM | 20 | AM |
| 126 | AM | 20 | AM |
| 127 | AM | 20 | AM |
| 128 | AM | 20 | AM |
| 129 | AM | 20 | AM |
| 130 | AM | 20 | AM |
| 131 | AM | 20 | AM |
| 132 | AM | 20 | AM |
| 133 | AM | 20 | AM |
| 134 | AM | 20 | AM |
| 135 | AM | 20 | AM |
| 136 | AM | 20 | AM |
| 137 | AM | 20 | AM |
| 138 | AM | 20 | AM |
| 139 | AM | 20 | AM |
| 140 | AM | 20 | AM |
| 141 | AM | 20 | AM |
| 142 | AM | 20 | AM |
| 143 | AM | 20 | AM |
| 144 | AM | 20 | AM |
| 145 | AM | 20 | AM |
| 146 | AM | 20 | AM |
| 147 | AM | 20 | AM |
| 148 | AM | 20 | AM |
| 149 | AM | 20 | AM |
| 150 | AM | 20 | AM |
| 151 | AM | 20 | AM |
| 152 | AM | 20 | AM |
| 153 | AM | 20 | AM |
| 154 | AM | 20 | AM |
| 155 | AM | 20 | AM |
| 156 | AM | 20 | AM |
| 157 | AM | 20 | AM |
| 158 | AM | 20 | AM |
| 159 | AM | 20 | AM |
| 160 | AM | 20 | AM |
| 161 | AM | 20 | AM |
| 162 | AM | 20 | AM |
| 163 | AM | 20 | AM |
| 164 | AM | 20 | AM |
| 165 | AM | 20 | AM |
| 166 | AM | 20 | AM |
| 167 | AM | 20 | AM |
| 168 | AM | 20 | AM |
| 169 | AM | 20 | AM |
| 170 | AM | 20 | AM |
| 171 | AM | 20 | AM |
| 172 | AM | 20 | AM |
| 173 | AM | 20 | AM |
| 174 | AM | 20 | AM |
| 175 | AM | 20 | AM |
| 176 | AM | 20 | AM |
| 177 | AM | 20 | AM |
| 178 | AM | 20 | AM |
| 179 | AM | 20 | AM |
| 180 | AM | 20 | AM |
| 181 | AM | 20 | AM |
| 182 | AM | 20 | AM |
| 183 | AM | 20 | AM |
| 184 | AM | 20 | AM |
| 185 | AM | 20 | AM |
| 186 | AM | 20 | AM |
| 187 | AM | 20 | AM |
| 188 | AM | 20 | AM |
| 189 | AM | 20 | AM |
| 190 | AM | 20 | AM |
| 191 | AM | 20 | AM |
| 192 | AM | 20 | AM |
| 193 | AM | 20 | AM |
| 194 | AM | 20 | AM |
| 195 | AM | 20 | AM |
| 196 | AM | 20 | AM |
| 197 | AM | 20 | AM |
| 198 | AM | 20 | AM |
| 199 | AM | 20 | AM |
| 200 | AM | 20 | AM |

01 5
4 5
58 70 72 88 90 94 102 110 51 113 123 125 146 149 189 214 219
BOOK C(3)=0, NAREA=2, NREFS=0, NHICHN=204, NR40=0, REFD
C(3)=1, C(4)=1, C(6)=1, WT=0, C(10)=1, C(8)=1
DECAY=58.5840, COUNT=10, WT=0.0337, CURT=2.0
TCU 0.0338, CUDFC=1058.584, REFD

01 5
4 5
58 70 72 88 90 94 102 110 51 113 123 125 146 149 189 214 219
BOOK C(3)=0, NAREA=2, NREFS=0, NHICHN=204, NR40=0, REFD
C(3)=1, C(4)=1, C(6)=1, WT=0, C(10)=1, C(8)=1
DECAY=58.5840, COUNT=10, WT=0.0337, CURT=2.0
TCU 0.0338, CUDFC=1058.584, REFD

01 5
4 5
58 70 72 88 90 94 102 110 51 113 123 125 146 149 189 214 219
BOOK C(3)=0, NAREA=2, NREFS=0, NHICHN=204, NR40=0, REFD
C(3)=1, C(4)=1, C(6)=1, WT=0, C(10)=1, C(8)=1
DECAY=58.5840, COUNT=10, WT=0.0337, CURT=2.0
TCU 0.0338, CUDFC=1058.584, REFD

E .840 +OR- 55^{X0}.249 +OR- 8^{B0}.337 +OR- .13692+004 +OR- AREA .61093+004
 .136 +OR- .188 +OR- .16361+003 +OR- .74122+003

E 1.810 +OR- 115^{X0}.042 +OR- 19^{B0}.692 +OR- .13283+003 +OR- AREA .10452+004
 .940 +OR- 4.248 +OR- .33023+002 +OR- .50457+003

E 2.120 +OR- 135^{X0}.013 +OR- 19^{B0}.188 +OR- .73560+002 +OR- AREA .57112+003
 1.237 +OR- 5.442 +OR- .24698+002 +OR- .22936+003

OLD GAIN = .1602032-001 INTERCEPT = -.3869112-001
 BASE LINE SHIFT = .2415128+001 NEW GAIN = .1450000-001 RATIO = .1103076+001

E .845 +OR- 58^{X0}.322 +OR- 8^{B0}.198 +OR- .12147+004 +OR- AREA .61645+004
 .127 +OR- .205 +OR- .11347+003 +OR- .58612+003

E .810 +OR- 124^{X0}.247 +OR- 24^{B0}.105 +OR- .12001+003 +OR- AREA .10444+004
 .801 +OR- 3.750 +OR- .24175+002 +OR- .23970+003

E .120 +OR- 146^{X0}.452 +OR- 33^{B0}.923 +OR- .63984+002 +OR- AREA .68053+003
 1.518 +OR- 9.525 +OR- .18211+002 +OR- .22922+003

CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical data, representing a cross product matrix. Values range from 0.0000 to 0.9999.

INV CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical data, representing the inverse of the cross product matrix. Values range from 0.0000 to 0.9999.

INV ERROR = .1109123-004

IDENTITY MATRIX

Table with 10 columns and 10 rows of numerical data, representing an identity matrix. The diagonal elements are 1.0000 and off-diagonal elements are 0.0000.

INTERFERENCES IN PER CENT

| | | | | | | | | | | |
|------------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|------------|
| BR-1 | CU64 .33129+002 | AL28 .34090+000 | S370 .58111+002 | NA24 .42128-001 | CA49 .12584-001 | K420 .30729+000 | CL38 .52735-003 | MN56 .50116+001 | MG27 .16896+001 | |
| CU64 | AL28 .98761-002 | S370 .53347+000 | NA24 .27170+000 | CA49 .61236+000 | K420 .95915-002 | CL38 .14568+000 | MN56 .68528-001 | MG27 .19022+000 | | |
| AL28 | S370 .41628+000 | NA24 .69045+001 | CA49 .10878+001 | K420 .78202+001 | CL38 .94716+001 | MN56 .11703+002 | MG27 .27300+001 | | | |
| S370 | NA24 .16899+002 | CA49 .95725+002 | K420 .21513+000 | CL38 .12249+002 | MN56 .69020+000 | MG27 .32810+000 | | | | |
| NA24 | CA49 .30084+002 | K420 .16606+000 | CL38 .34992+001 | MN56 .13486+001 | MG27 .21339+001 | | | | | |
| CA49 | K420 .26257+000 | CL38 .13208+002 | MN56 .10819+001 | MG27 .66034+000 | | | | | | |
| K420 | CL38 .26989+001 | MN56 .12830+000 | MG27 .25414+001 | | | | | | | |
| CL38 | MN56 .14538+001 | MG27 .11412-001 | | | | | | | | |
| MN56 | MG27 .35122+002 | | | | | | | | | |
| (ATRANS * OMEGA * RHO) | | | | | | | | | | |
| 10 | .31905+007 | .24420+007 | .19313+009 | .28651+005 | .63610+007 | .28482+006 | .51989+005 | .57712+007 | .17032+009 | .15614+007 |

ISOTOPE REMOVED AL28

ISOTOPE REMOVED S370

ISOTOPE REMOVED CA49

ITERATION, 0 OF 0 MAX ITER

| | | | | | | | |
|-----------|-----------|------------|-----------|-----------|-----------|-----------|-----------|
| | BR-1 | CU64 | NA24 | K420 | CL38 | MN56 | MG27 |
| SOL (PPM) | .1256+003 | -.6845+002 | .9456+001 | .2940+005 | .2755+004 | .2571+003 | .6669+005 |
| STD DEV | .8029+002 | .5875+002 | .8508+001 | .2306+004 | .1321+003 | .3332+001 | .2445+006 |
| T VALUES | .1565+001 | -.1131+001 | .1112+001 | .1231+002 | .2066+002 | .7601+002 | .4549+000 |

F VALUE = .1647+002 CORR COEF = .9936+000

ISOTOPE REMOVED CU64

ITERATION, 0 OF 0 MAX ITER

| | | | | | | |
|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| | BR-1 | NA24 | K420 | CL38 | MN56 | MG27 |
| SOL (PPM) | .7461+002 | .9575+001 | .2845+005 | .2739+004 | .2559+003 | .5629+005 |
| STD DEV | .5330+002 | .8508+001 | .2306+004 | .1321+003 | .3332+001 | .1245+006 |
| T VALUES | .1400+001 | .1125+001 | .1234+002 | .2089+002 | .7557+002 | .4554+000 |

F VALUE = .1921+002 CORR COEF = .9934+000

ISOTOPE REMOVED MG27

ITERATION, 0 OF 0 MAX ITER

| | | | | | |
|-----------|-----------|-----------|-----------|-----------|-----------|
| | BR-1 | NA24 | K420 | CL38 | MN56 |
| SOL (PPM) | .7762+002 | .1109+002 | .2864+005 | .2760+004 | .2578+003 |
| STD DEV | .5268+002 | .8176+001 | .2283+004 | .1322+003 | .2582+001 |
| T VALUES | .1474+001 | .1357+001 | .1255+002 | .2088+002 | .9936+002 |

F VALUE = .2305+002 CORR COEF = .9934+000

ISOTOPE REMOVED NA24

ITERATION, 0 OF 0 MAX ITER

| | | | | |
|-----------|-----------|-----------|-----------|-----------|
| | BR-1 | K420 | CL38 | MN56 |
| SOL (PPM) | .7980+002 | .2896+005 | .2788+004 | .2578+003 |
| STD DEV | .5260+002 | .2269+004 | .1305+003 | .2581+001 |
| T VALUES | .1517+001 | .1276+002 | .2134+002 | .9906+002 |

F VALUE = .2681+002 CORR COEF = .9934+000

ISOTOPE REMOVED BR-1

ITERATION, 0 OF 0 MAX ITER

| | | | |
|-----------|-----------|-----------|-----------|
| | K420 | CL38 | MN56 |
| SOL (PPM) | .2932+005 | .2760+004 | .2578+003 |
| STD DEV | .2255+004 | .1305+003 | .2317+001 |
| T VALUES | .1300+002 | .2130+002 | .1113+003 |

F VALUE = .3841+002 CORR COEF = .9932+000

DEG OF FREEDOM = 112. QUAL OF FIT = .204+001

ABS DEV PER PT = .7490734+001 NO OF PTS = 116

8
SAMPLE WT (GM) = .39600-001 IRR TIME (MIN) = .15000+002 DECAY TIME (MIN) = .12655+003 COUNT TIME = .10000+002
FLUX REF WT (GM) = .33900-001 REF DECAY TIME (MIN) = .10586+004 REF COUNT TIME (MIN) = .20000+001 REF HLIFE (MIN) = .77000+001
BKG TIME (MIN) = .60000+002
Y NORM CONSTANT = .25088707+001

E .845 +OR- 55.292 .138 +OR- BO 6.082 .182 +OR- YO .13030+004 +OR- AREA .56959+004
+OR- .16611+003 +OR- .73599+003

E 1.810 +OR- 115.015 .733 +OR- BO 14.820 2.503 +OR- YO .13537+003 +OR- AREA .91004+003
+OR- .34285+002 +OR- .25706+003

E 2.120 +OR- 134.848 1.326 +OR- BO 18.353 5.758 +OR- YO .55091+002 +OR- AREA .42569+003
+OR- .20351+002 +OR- .18110+003

OLD GAIN = 1605694-001 INTERCEPT = -.4161721-001
BASE LINE SHIFT = .2591852+001 NEW GAIN = .1450000-001 RATIO = .1105872+001

E .845 +OR- 58.312 .126 +OR- BO 7.691 .193 +OR- YO .11605+004 +OR- AREA .57041+004
+OR- .11402+003 +OR- .56950+003

E 1.810 +OR- 124.372 .695 +OR- BO 19.695 2.609 +OR- YO .11896+003 +OR- AREA .93571+003
+OR- .25040+002 +OR- .22279+003

E 2.120 +OR- 146.336 1.328 +OR- BO 23.668 6.476 +OR- YO .58180+002 +OR- AREA .45962+003
+OR- .18070+002 +OR- .18319+003

CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical data, representing a cross product matrix. Values range from 0.0000 to 0.0009.

INV CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical data, representing the inverse of the cross product matrix. Values range from 0.0000 to 0.0009.

INV ERROR = .5339334E+05

IDENTITY MATRIX

Table with 10 columns and 10 rows of numerical data, representing an identity matrix. Values are mostly 0.0000, with 1.0000 on the diagonal.

INTERFERENCES IN PER CENT

| | | | | | | | | | | |
|------------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|------------|
| BR71 | CU64 .31393+002 | AL28 .33853+000 | S370 .71415-002 | NA24 .55203-001 | CA49 .15209-001 | K420 .29233+000 | CL38 .77592-002 | MN56 .54066+001 | MG27 .16927+001 | |
| CU64 | AL24 .98517-002 | S370 .50425+000 | NA24 .21223+000 | CA49 .54822+000 | K420 .60307-002 | CL38 .16549+000 | MN56 .74224-001 | MG27 .18237+000 | | |
| AL28 | S370 .47196+000 | NA24 .83994+001 | CA49 .12023+001 | K420 .85460+001 | CL38 .83429+001 | MN56 .11327+002 | MG27 .23395+001 | | | |
| S370 | NA24 .14227+002 | CA49 .96799+002 | K420 .11551+000 | CL38 .16755+002 | MN56 .62462+000 | MG27 .27854+000 | | | | |
| NA24 | CA49 .25855+002 | K420 .76683-001 | CL38 .48933+001 | MN56 .12103+001 | MG27 .22690+001 | | | | | |
| CA49 | K420 .14319+000 | CL38 .18490+002 | MN56 .99180+000 | MG27 .58222+000 | | | | | | |
| K420 | CL38 .20343+001 | MN56 .97560-001 | MG27 .26394+001 | | | | | | | |
| CL38 | MN56 .16975+001 | MG27 .57620-003 | | | | | | | | |
| MN56 | MG27 .33427+002 | | | | | | | | | |
| (ATRANS * OMEGA * RHO) | | | | | | | | | | |
| 10 | .34021+007 | .26039+007 | .20593+009 | .30551+005 | .67827+007 | .30370+006 | .55436+005 | .61539+007 | .18161+009 | .16649+007 |

ISOTOPE REMOVED AL25

ISOTOPE REMOVED S370

ISOTOPE REMOVED CA49

ISOTOPE REMOVED MG27

ITERATION, 0 OF 0 MAX ITER

| | BR-1 | CU64 | NA24 | K420 | CL38 | MN56 |
|-----------|-----------|-------------|-----------|-----------|-----------|-----------|
| SOL (PPM) | .9599+002 | .3198+002 | .2205+002 | .2908+005 | .3064+004 | .2273+003 |
| STD DEV | .7256+002 | .5756+002 | .6845+001 | .2151+004 | .2328+003 | .2929+001 |
| T VALUES | .1321+001 | .5556+000 | .3625+001 | .1052+002 | .1212+002 | .1015+000 |
| F VALUE = | .1421+002 | CORR COEF = | | .9933+000 | | |

ISOTOPE REMOVED CU64

ITERATION, 0 OF 0 MAX ITER

| | BR-1 | NA24 | K420 | CL38 | MN56 |
|-----------------------------------------------|-----------|-------------|-----------|-----------|-----------|
| SOL (PPM) | .1181+003 | .3188+002 | .2904+005 | .3060+004 | .2973+003 |
| STD DEV | .4944+002 | .0835+001 | .2148+004 | .2529+003 | .2929+001 |
| T VALUES | .2388+001 | .3609+001 | .1352+002 | .1210+002 | .1015+000 |
| F VALUE = | .2305+002 | CORR COEF = | | .9933+000 | |
| DEG OF FREEDOM = 110. QUAL OF FIT = .201+001 | | | | | |
| ABS DEV PER PT = .6427585+001 NO OF PTS = 116 | | | | | |

SAMPLE WT (GM) = .35100±001 ERR TIME (MIN) = .15000±002 DECAY TIME (MIN) = .13856±003 COUNT TIME = .10000±002
FLUX REF WT (GR) = .33800±001 REF DECAY TIME (MIN) = .10588±004 REF COUNT TIME (MIN) = .20000±001 REF HLIFE (MIN) =
BKG TIME (MIN) = .60000±002
Y NORM CONSTANT = .21513891±001

E .845 X0 BO YO AREA
 +OR= 55.276 8.149 :99178+003 :53588+004
 .153 :211 :14233+003 :83468+003

 E 1.816 X0 BO YO AREA
 +OR= 115.051 18.869 :99181+002 :72129+003
 .952 :3.700 :29861+002 :26018+003

 E 2.120 X0 BO YO AREA
 +OR= 134.815 21.825 :45182+002 :31786+003
 1.721 :8.959 :17725+002 :18252+003

 OLD GAIN = 16.05597-001 INTERCEPT = -4.144175-001
 BASE LINE SHIFT = .2581081+001 NEW GAIN = .145000-001 RATIO = .11060+001

E .845 X0 BO YO AREA
 +OR= 58.316 7.750 :88458+003 :43649+004
 .144 :223 :99318+002 :49008+003

 E 1.816 X0 BO YO AREA
 +OR= 124.444 25.217 :87978+002 :75136+003
 .939 :4.085 :22449+002 :21336+003

 E 2.120 X0 BO YO AREA
 +OR= 146.316 35.145 :39344+002 :41344+003
 1.964 :13.095 :13595+002 :17957+003

CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical data, representing the cross product matrix.

INV CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical data, representing the inverse cross product matrix.

INV ERROR = 0.1350641-004

IDENTITY MATRIX

Table with 10 columns and 10 rows of numerical data, representing the identity matrix.

INTERFERENCES IN PER CENT

| | | | | | | | | | | |
|------------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|------------|
| BR71 | CU64 .32411+002 | AL28 .37002+000 | S370 .32630+002 | NA24 .40418+001 | CA49 .53588+002 | K420 .35082+000 | CL38 .10183+001 | MN56 .50891+004 | M627 .28298+001 | |
| CU64 | AL28 .14541-001 | S370 .58026+000 | NA24 .17022+000 | CA49 .64588+000 | K420 .71081+002 | CL38 .21149+000 | MN56 .80807+001 | M627 .28992+000 | | |
| AL28 | S370 .69347-001 | NA24 .84343+001 | CA49 .46604+000 | K420 .82571+001 | CL38 .29821+001 | MN56 .11007+002 | M627 .20909+001 | | | |
| S370 | NA24 .53588+001 | CA49 .96512+002 | K420 .10771+000 | CL38 .11029+002 | MN56 .40000+000 | M627 .10788+000 | | | | |
| NA24 | CA49 .13833+002 | K420 .18322+001 | CL38 .34089+001 | MN56 .90915+000 | M627 .20904+001 | | | | | |
| CA49 | K420 .12838+000 | CL38 .20469+002 | MN56 .76411+000 | M627 .32047+000 | | | | | | |
| K420 | CL38 .18593+001 | MN56 .10365+000 | M627 .31512+001 | | | | | | | |
| CL38 | MN56 .14840+001 | M627 .19341-002 | | | | | | | | |
| MN56 | M627 .32414+002 | | | | | | | | | |
| (ATRANS * OMEGA * RHO) | | | | | | | | | | |
| 10 | .29173+007 | .22329+007 | .17659+009 | .26198+005 | .58163+007 | .26043+006 | .47537+005 | .52770+007 | .15578+009 | .14277+007 |

ISOTOPE REMOVED AL28

ISOTOPE REMOVED S370

ISOTOPE REMOVED CA49

ISOTOPE REMOVED MG27

ITERATION, 0 OF 0 MAX ITER

| | | | | | | |
|-----------|-----------|------------|-----------|-----------|-----------|-----------|
| SOL (PPM) | BR-1 | CU64 | NA24 | K420 | CL38 | MN56 |
| STD DEV | .5021+002 | -.4295+002 | .3240+002 | .2822+005 | .2538+004 | .2405+003 |
| T VALUES | .5909+002 | -.4629+002 | .9474+001 | .2121+004 | .2928+003 | .2723+001 |
| | .9512+000 | -.9279+000 | .3420+001 | .1331+002 | .8643+001 | .8831+002 |

F VALUE = .1923+002 CORR COEF = .9944+000

ISOTOPE REMOVED CU64

ITERATION, 0 OF 0 MAX ITER

| | | | | | |
|-----------|-----------|-----------|-----------|-----------|-----------|
| SOL (PPM) | BR-1 | NA24 | K420 | CL38 | MN56 |
| STD DEV | .2559+002 | .3267+002 | .2628+005 | .2536+004 | .2405+003 |
| T VALUES | .4108+002 | .9467+001 | .2118+004 | .2928+003 | .2723+001 |
| | .6229+000 | .3451+001 | .1335+002 | .8663+001 | .8831+002 |

F VALUE = .2307+002 CORR COEF = .9943+000

ISOTOPE REMOVED BR-1

ITERATION, 0 OF 0 MAX ITER

| | | | | |
|-----------|-----------|-----------|-----------|-----------|
| SOL (PPM) | NA24 | K420 | CL38 | MN56 |
| STD DEV | .3287+002 | .2840+005 | .2523+004 | .2405+003 |
| T VALUES | .9462+001 | .2109+004 | .2920+003 | .2544+001 |
| | .3474+001 | .1347+002 | .8641+001 | .9454+002 |

F VALUE = .2884+002 CORR COEF = .9943+000

DEG OF FREEDOM = 111. QUAL OF FIT = .112+001

ABS DEV PER PT = .4230796+001 NO OF PTS = 110

SAMPLE WT (GM) = .34700-001 IRR TIME (MIN) = 15000+002 DECAY TIME (MIN) = .58584+002 COUNT TIME
FLUX RATE WT (GM) = .33800-001 REF DECAY TIME (MIN) = .10586+004 REF COUNT TIME (MIN) = .20000+00
BKG TIME (MIN) = .60000+002
Y NORM CONSTANT = .22553837+001

| | | | | | | | | |
|------------|------|------------------------|------|------------------------|------|--------------------------------|------|----------------------------------|
| E .845 | +OR- | XO 55.270 .145 | +OR- | BO 6.470 .205 | +OR- | YO .12431+004 .15573+003 | +OR- | AREA .56045+004 .71751+003 |
| E 1.810 | +OR- | XO 115.012 .894 | +OR- | BO 18.544 3.817 | +OR- | YO .13245+003 .33133+002 | +OR- | AREA .10109+004 .29259+003 |
| E 2.120 | +OR- | XO 135.616 2.102 | +OR- | BO 31.956 13.915 | +OR- | YO .77483+002 .29524+002 | +OR- | AREA .77635+003 .38033+003 |

OLD GAIN = .1593234-001 INTERCEPT = -.3289121-001
 BAS- LINE SHFT = .2064431+001 NEW GAIN = .1450000-001 RATIO = .1095427+001

| | | | | | | | | |
|------------|------|------------------------|------|-----------------------|------|--------------------------------|------|----------------------------------|
| E .845 | +OR- | XO 58.311 .132 | +OR- | BO 8.073 .211 | +OR- | YO .11186+004 .10986+003 | +OR- | AREA .56235+004 .56299+003 |
| E 1.810 | +OR- | XO 123.744 .821 | +OR- | BO 22.843 3.375 | +OR- | YO .12022+003 .25637+002 | +OR- | AREA .10184+004 .26490+003 |
| E 2.120 | +OR- | XO 146.151 1.351 | +OR- | BO 32.909 8.686 | +OR- | YO .71904+002 .17654+002 | +OR- | AREA .73112+003 .22529+003 |

CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical values, representing a cross product matrix. Values range from approximately 0.1 to 0.8.

INV CROSS PRODUCT MATRIX

Table with 10 columns and 10 rows of numerical values, representing the inverse of the cross product matrix. Values range from approximately -0.7 to 0.9.

INV ERROR = .8257627-005

IDENTITY MATRIX

Table with 10 columns and 10 rows of numerical values, representing an identity matrix. The diagonal elements are 1.0, and off-diagonal elements are 0.0.

INTERFERENCES IN PER CENT

| | | | | | | | | | | |
|------------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|------------|
| BR-1 | CU64 .31370+002 | AL28 .29777+000 | S370 .12574-002 | NA24 .32211-001 | CA49 .78065-002 | K420 .33378+000 | CL38 .15423-001 | MN56 .43482+001 | MG27 .16772+001 | |
| CJ04 | AL28 .75818-002 | S370 .33985+000 | NA24 .23205+000 | CA49 .37415+000 | K420 .25378-002 | CL38 .68428-001 | MN56 .95285-001 | MG27 .17587+000 | | |
| AL28 | S370 .13087+001 | NA24 .69437+001 | CA49 .25843+001 | K420 .70934+001 | CL38 .11727+002 | MN56 .11415+002 | MG27 .30461+001 | | | |
| S370 | NA24 .30098+002 | CA49 .95149+002 | K420 .99334-001 | CL38 .75425+001 | MN56 .33584+000 | MG27 .42914+000 | | | | |
| NA24 | CA49 .49066+002 | K420 .92524-001 | CL38 .29996+001 | MN56 .11807+001 | MG27 .20849+001 | | | | | |
| CA49 | K420 .12595+000 | CL38 .79903+001 | MN56 .74939+000 | MG27 .92140+000 | | | | | | |
| K420 | CL38 .34394+001 | MN56 .20004+000 | MG27 .27401+001 | | | | | | | |
| CL38 | MN56 .15568+001 | MG27 .30060-001 | | | | | | | | |
| MN56 | MG27 .38096+002 | | | | | | | | | |
| (ATRANS * OMEGA * RHO) | | | | | | | | | | |
| 10 | .30583+007 | .23408+007 | .18512+009 | .27464+005 | .60974+007 | .27302+006 | .49835+005 | .55321+007 | .16326+009 | .14967+007 |

ISOTOPE REMOVED AL28

ISOTOPE REMOVED S370

ITERATION, 0 OF 0 MAX ITER

| | | | | | | | | |
|-----------|-----------|-----------|-----------|------------|-----------|-----------|-----------|------------|
| SOL (PPM) | BR-1 | CU64 | NA24 | CA49 | K420 | CL38 | MN56 | MG27 |
| STD DEV | .7947+002 | .5153+002 | .5327+002 | -.2857+005 | .3454+005 | .2770+004 | .2119+003 | -.4250+004 |
| T VALUES | .1646+001 | .9069+000 | .3278+001 | -.4652+001 | .2401+004 | .9354+002 | .3051+001 | -.2359+005 |

F VALUE = .1434+002 CORR COEF = .9891+000

ISOTOPE REMOVED CA49

ITERATION, 0 OF 0 MAX ITER

| | | | | | | | |
|-----------|-----------|-----------|------------|-----------|-----------|-----------|------------|
| SOL (PPM) | BR-1 | CU64 | NA24 | K420 | CL38 | MN56 | MG27 |
| STD DEV | .8202+002 | .4785+002 | -.1536+002 | .3567+005 | .2772+004 | .2112+003 | -.4250+004 |
| T VALUES | .1080+001 | .8425+000 | -.2112+001 | .1462+002 | .9364+002 | .3022+001 | -.2350+005 |

F VALUE = .1638+002 CORR COEF = .9884+000

ISOTOPE REMOVED NA24

ITERATION, 0 OF 0 MAX ITER

| | | | | | | |
|-----------|-----------|-----------|-----------|-----------|-----------|------------|
| SOL (PPM) | BR-1 | CU64 | K420 | CL38 | MN56 | MG27 |
| STD DEV | .8302+002 | .4823+002 | .3450+005 | .2750+004 | .2107+003 | -.4250+004 |
| T VALUES | .1093+001 | .8500+000 | .1456+002 | .2951+002 | .3001+001 | -.2350+005 |

F VALUE = .1911+002 CORR COEF = .9882+000

ISOTOPE REMOVED MG27

ITERATION, 0 OF 0 MAX ITER

| | | | | | |
|-----------|-----------|-----------|-----------|-----------|-----------|
| SOL (PPM) | BR-1 | CU64 | K420 | CL38 | MN56 |
| STD DEV | .8184+002 | .4857+002 | .3482+005 | .2750+004 | .2103+003 |
| T VALUES | .1088+001 | .8577+000 | .1469+002 | .2945+002 | .2181+001 |

F VALUE = .2293+002 CORR COEF = .9882+000

ISOTOPE REMOVED CU64

ITERATION, 0 OF 0 MAX ITER

| | | | | |
|-----------|-----------|-----------|-----------|-----------|
| SOL (PPM) | BR-1 | K420 | CL38 | MN56 |
| STD DEV | .1174+003 | .3477+005 | .2747+004 | .2103+003 |
| T VALUES | .5147+002 | .2368+004 | .9334+002 | .2181+001 |

F VALUE = .2866+002 CORR COEF = .9882+000

DEG OF FREEDOM = 111. QUAL OF FIT = .340+001

ABS DEV PER PT = .7449490+001 NO OF PTS = 116

APPENDIX M-LISTING OF THE UPDATING PROGRAMS

The first program is used to record the content of the magnetic tape obtained from NEA-CPL to the main storage of a UNIVAC-1106 computer used in Bogaziçi University. The main program and the sample problem data were recorded on to the same file. The main program and the sample input data are loaded separately into two different files, CORGAM and DAFCA respectively. The namelist is removed from the 96. line to the 4. line due to the requirements of the UNIVAC computing system. These are done by the help of the 2. and 3. programs. The 4. program is used to update the coding symbols for UNIVAC-1106 Card Reader and Line Printer. In order to execute the program by the sample input data, the 5. program is used to remove the main program into the 2. cycle of the file, because the first cycle contains the program in data form, by removing it to the 2. cycle, it is transformed into an executable form. Finally the 6. program is used to run the program using the sample data as input.

PROGRAM-1 Records the tape in the main storage

@RUN
@ASG,T,J T.,6C9,D80218
@MOVE T.,7
@DELETE,C BOGRUN.
@ASG,UP BOGRUN.
@USE F.,BOGRUN.
@COB,ISE

IDENTIFICATION DIVISION .
PROGRAM-ID. TEYP-OKUMA.
ENVIRONMENT DIVISION.
CONFIGURATION SECTION.
SOURCE-COMPUTER. U-1106
OBJECT-COMPUTER. U-1106
INPUT-OUTPUT SECTION.
FILE-CONTROL.
SELECT T ASSIGN TO UNISERVO T.
DATA DIVISION.
FILE SECTION.
FD T

LABEL RECORDS OMITTED
RECORDING MODE IS F.
01 T-REC PIC X(80).
FD F
RECORDING MODE IS SDF.
01 F-REC PIC X(80).

PROCEDURE DIVISION.
AC. OPEN INPUT T WITH NO REWIND.
OPEN OUTPUT F.
OKU. READ T AT END GO TO SON.
MOVE T-REC TO F-REC.
WRITE F-REC.
GO TO OKU.
SON. CLOSE T WITH NO REWIND.
CLOSE F.
STOP RUN.

@MAP,IN
@XQT
@DATA,L BOGRUN.
@FIN

PROGRAM-2

1.Update of main program

```
@RUN
@ASG,A      BOGRUN.
@DELETE,C   CORGAM.
@ASG,UP     CORGAM.
@DATA,L     BOGRUN.,CORGAM.
-4
```

```
NAMELIST /NAM1/ISOCHN,MODE,NBKGD,NUBKG,NSHFT,LPP,NORM,NCTIN
1NAREA,NCHCU,NCULL,NCUUL,NRHO,NUMVAR,NVAR,NREFS,ITER,NR,NWW,
2LOCHN,NHICHN,NCHAN,MC,MVAR/NAM2/C/NAM3/T,DECAY,COUNT,WT,CUC
3T,WTCU,AREA,CUDEEC,BKGTIM,FHRL,UP,XNCOM
-93,96
v1296,1851
@END
@FIN
```

PROGRAM-3

Load sample problem data

```
@RUN
@ASG,A      BOGRUN.
@DELETE,C   DAFCA.
@ASG,UP     DAFCA.
@DATA,L     BOGRUN.,DAFCA.
-1,1295
@END
@FIN
```

PROGRAM-4

2. Update of main program

```
@RUN
@ASG,A      CORGAM.
@ED,L       CORGAM.,CORGAM.
C /READ(1,/READ(5,/ ALL
C /WRITE(3,/WRITE(6,/ ALL
@END
@FIN
```

PROGRAM - 5

```
@RUN
@ASG,A      CORGAM.
@ASG,UP     CORGAM(+1).
@FTN,US     CORGAM(+1)
@ADD,DP     CORGAM.
@FIN
```

3. Update of main program.

PROGRAM - 6

```
@RUN
@ASG,A      CORGAM.
@ASG,A      DAFCA.
@DELETE,C 2.
@DELETE,C 10.
@ASG,UP     2.
@ASG,UP     10.
@FTN,S      CORGAM.MAIN
@PREP       CORGAM.
@MAP,I      CORGAM.MAIN
IN CORGAM.MAIN
LIB CORGAM.
@XQT        CORGAM.MAIN
@ADD,DP     DAFCA.
@FIN
```

Run of main program with
sample problem inputdata.

APPENDIX N

The following is a list of the computer programs coded for experimental data processing, used in common neutron activation analysis.

A-Computer codes available in ARGONNE NATIONAL LABORATORY (ANL-7411)

- 1- TOAD : "A computer code for processing gamma-ray spectra"
Reactor code abstract: 333, Coded for: UNIVAC-1108
Language: Fortran IV Issued at: June, 1968
- 2- CORGAM : "A correlation algorithm for gamma-ray spectra"
Reactor code abstract: 390, Coded for: IBM 360/50
Language: Fortran IV-G Issued at: May, 1969
- 3- ALPHA-M : "Resolution of gamma-ray spectra" (improved version of ALPHA)
Reactor code abstract: 413, Coded for: IBM 360
Language: Fortran IV, CPL Issued at: December 1970
- 4- GSSLRN1B : "Least squares photopeak spectra code"
Reactor code abstract: 457, Coded for: UNIVAC-1108
Language: Fortran IV, CPL Issued at: June 1972
- 5- GASPAN : "Complex gamma-ray spectra analysis"
Reactor code abstract: 485, Coded for: CDC 6600
Language: Fortran IV Issued at: January 1972

B-Computer codes available in NEA-CPL, Ispra

- 6- GRETEL : "Ge(L1) gamma spectrum unfolding"
Program abstract: E 433 Coded for: IBM 370
Language: Fortran IV, CPL, Assembler
- 7- RADAK : "Neutron and gamma spectrum unfolding"
Program abstract: E 487 Coded for: IBM 370
Language: Fortran IV